

TABLE OF CONTENTS

Fetal Monitors	2
Edan F9 Express Fetal Monitor	3
GE Corometrics 171 Series Fetal Monitor	3
Luxury Trolley for F series Fetal Monitors	4
Philips Agilent Series 50XM Fetal Monitor M1350B	4
Philips Avalon FM30	5
GE/Corometrics 172 Twins Fetal Monitor	5
GE Corometrics 120 Fetal Monitor	6
Edan F6 Twins Fetal Monitor	7
Fetalgard Lite Fetal Monitor	7
GE Corometrics 129 Fetal Monitor w. Twins and maternal monitoring	7
Edan SD5 & SD6 Ultrasonic Tabletop Doppler	8
Edan Sonotrax Series Fetal Dopplers	8



FETAL MONITORS

Auxo Medical prides itself on offering medical professionals the most advanced fetal monitors to date. Dependable during all stages of pregnancy, Auxo Medical's fetal monitoring technology offers measurements and monitoring for the newest addition to your patient's family.


Adaptable for use in all situations, our monitors cover the entire pregnancy continuum, from antepartum to intrapartum. Whether you need a fetal monitor with quick printing, long lasting battery life, high resolution displays, wall-mounting racks, or easy transportation, Auxo Medical offers industry leading machines for every need and budget. Providing accurate data for mother, child, and physician, our fetal monitors offer your facility top of the line technology for your obstetrics and patient care needs.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

REFURBISHED



EDAN F9 EXPRESS FETAL MONITOR

The **Edan F9 Express Fetal Monitor** is a state-of-the-art device designed to provide comprehensive monitoring solutions for fetal and maternal health in various medical settings. Equipped with a vibrant 12.1-inch color touch screen display, this advanced monitor is fully loaded to accommodate twins' fetal heart rate (FHR), maternal parameters, and direct ECG/intrauterine pressure (DECG/IUP) monitoring. The F9 Express is ideal for use in large hospital delivery rooms, private obstetricians' offices, and antepartum clinics, ensuring seamless application across the full spectrum of antepartum, intrapartum, and postpartum care. Boasting an array of features, it includes essential parameters such as twins FHR, TOCO, event marking, and AFM, alongside maternal ECG, non-invasive blood pressure (NIBP), SpO2, and temperature monitoring. Its internal parameters cover IUP and DECG. The monitor's powerful lithium battery ensures four hours of continuous operation, while its 12-crystal pulse wave Doppler waterproof transducer guarantees accurate FHR detection. The device offers signal overlap verification for twins FHR with adjustable offset printing options, and users can benefit from 24-hour waveform playback and quick printing for stored data. With real-time printing speeds of 1, 2, or 3 cm/min, the F9 Express also features an optional built-in wireless module for seamless connectivity to the FTS-3 or Central Nursing Station, and software for efficient data transmission to a PC via Insight Software or CNS Lite.

[Read More](#)

SKU: AM-F9E

GE COROMETRICS 171 SERIES FETAL MONITOR

The Corometrics 170 Series

[Read More](#)

SKU: AM-COROMETRICS-171

Price: \$2,995.00



LUXURY TROLLEY FOR F SERIES FETAL MONITORS

Edan Luxury Trolley provides a great mounting solution for F Series Fetal Monitors. It features 4 wheels with clips and mounting plate with a clamp for steady placement. The storage basket comes in a set.

[Read More](#)

SKU: AM-MT-803



PHILIPS AGILENT SERIES 50XM FETAL MONITOR M1350B



Refurbished **Philips M1350B Fetal monitor**. This Monitor 50XM has five color- and mechanically-coded input sockets: Cardio for fetal heart rate (FHR) and fetal oxygen saturation SPO2 Cardio for fetal heart rate and Toco for uterine activity NIBP for external blood pressure, SpO2 for maternal pulse oximetry. When fetal transducer is connected, the monitor automatically selects the correct operating mode. Philips M1350B can Monitor the followings: SpO2 monitor FHR using both direct ECG (DECG) and ultrasound (US) maternal heart rate using ECG, monitor twins using either DECG and ultrasound or dual ultrasound monitor uterine activity (Toco) externally or internally.

[Read More](#)

SKU: AM-M1350B-50XM

PHILIPS AVALON FM30



The Avalon FM series fetal and maternal monitors are the first and only to offer automated coincidence detection (cross-channel verification) using Smart Pulse. This innovative feature allows automated maternal pulse detection via the Toco MP transducer without the need to monitor maternal SpO2 or ECG separately.

[Read More](#)

SKU: AM-FM30

GE/COROMETRICS 172 TWINS FETAL MONITOR



Whether you perform antepartum fetal monitoring in an office setting or in the hospital, Corometrics* 170 Series Fetal Monitor can help. You can trust our fetal monitors to help you improve the standard of care for your expecting mothers. We designed these compact, lightweight monitors to help you effectively and efficiently care for mom and baby by accurately and noninvasively tracking baby's heart rate for non-stress tests.

[Read More](#)

SKU: AM-COROMETRICS-172T

Price: \$3,683.00

GE COROMETRICS 120 FETAL MONITOR



Corometrics 120 maternal/fetal monitors help you manage the labor, delivery and recovery process for every patient. Every 120 series monitor features the clinical excellence that has set the standard in M/F monitoring for more than 30 years. Patented Smart BP, Heart Beat Coincidence, Spectra Alerts, DINAMAP non-invasive blood pressure algorithm and 9-crystal pulse Doppler ultrasound transducers... You can trust the data as an accurate assessment of mom and baby's progress. Labor can be unpredictable. Be prepared with a monitor that can help you manage the entire process - even in the most difficult cases. Smart BP feature - delays auto determination during contractions Vital Signs History feature saves maternal information in memory for quick review or print summary Recovery mode sets the recorder to standby allowing delivery to print only the vital signs and save paper


Keyboard option for detailed annotations Software Utility upgrade program to install future feature enhancements Song Player serenades both mother and family to celebrate the new arrival Interface with 340 Telemetry for ambulatory monitoring and Hydrotherapy when used with the Nautilus Watertight Transducers Interface capability with the OxiFirst® Fetal Pulse Oximetry system.

[Read More](#)

SKU: AM-COROMETRICS-120

Price: \$3,275.00

EDAN F6 TWINS FETAL MONITOR



The Edan F6 Fetal Monitor contains the power necessary for multi-bed units while maintaining a simplistic design suitable for private practice. This high-end hospital grade monitor offers one of the most advanced integrated monitoring systems for the delivery room.

[Read More](#)

SKU: AM-F6E



FETALGARD LITE FETAL MONITOR

FETALGARD Lite from Analogic can be purchased in a package or configured independently. Unique multifunctional features make it ideal for fetal monitoring in the clinic, office or home. The system is designed and manufactured in the United States with the quality, dependability, performance and sensitivity you expect from Summit Doppler.

[Read More](#)

SKU: AM-AL-FETALGARDLITE

GE COROMETRICS 129 FETAL MONITOR W. TWINS AND MATERNAL MONITORING



The Model 129 Fetal Monitor is a top of the line fully featured fetal monitor. It includes all Model 128 features plus Maternal 3-Lead ECG.

[Read More](#)

SKU: AM-COROMETRICS-129

Price: \$3,342.00

EDAN SD5 & SD6 ULTRASONIC TABLETOP DOPPLER



Ultrasonic tabletop Doppler SD5 and SD6 are the innovative development of Edan's Fetal Doppler category. The application of infrared transmission technology, remote control probe, along with the brand new design, brings out the SD Series to the top of the list.

SD5 tabletop Doppler provides unmatched accuracy in FHR detection, large digital display of FHR, remote control probe and built-in long lasting rechargeable battery, which enables it for the intended antepartum application since the 10th week pregnancy for midwives, private doctors, general practice, OB office and L&D outpatient service.

[Download Brochure](#)

[Read More](#)

SKU: AM-SD6, AM-SD5

EDAN SONOTRAX SERIES FETAL DOPPLERS



SonoTrax Series Ultrasonic Pocket Dopplers provide caregivers with the most cost effective solution to meet the routine obstetric and vascular requirements.

[Download Brochure](#)

[Read More](#)

SKU: AM-E-STV

INDEX

E

Edan F6 Twins Fetal Monitor 6
Edan F9 Express Fetal Monitor 2
Edan SD5 & SD6 Ultrasonic Tabletop Doppler 7
Edan Sonotrax Series Fetal Dopplers 8

F

Fetalgard Lite Fetal Monitor 6

G

GE Corometrics 120 Fetal Monitor 5
GE Corometrics 129 Fetal Monitor w. Twins and maternal monitoring 7
GE Corometrics 171 Series Fetal Monitor 3
GE/Corometrics 172 Twins Fetal Monitor 5

L

Luxury Trolley for F series Fetal Monitors 3

P

Philips Agilent Series 50XM Fetal Monitor M1350B 4
Philips Avalon FM30 4