

TABLE OF CONTENTS

General Biomed 2

 HTP 1500 20" X 57" REUSABLE LARGE WATER MATTRESS 3

 HTP 1500 * 20" X 44"REUSABLE LARGE WATER MATTRESS 3

 HTP 1500 HEAT THERAPY PUMP, NEW 3

 MD5040V UPS system delivers 5000 watts of Single-phase power with a 4.8
 KWH battery bank 4

 InBody 970S Body Composition Analyzer 4



GENERAL BIOMED



HTP 1500 20" X 57" REUSABLE LARGE WATER MATTRESS

The **Adroit Medical Systems HTP 1500** is a reusable large water mattress measuring **20 inches by 57 inches**, designed to deliver effective heat therapy and pressure redistribution in clinical environments. Built for durability and ease of use, this water mattress offers consistent and adjustable warmth to enhance patient comfort during treatment, making it an essential tool for therapeutic and patient management applications.

**Image is for illustration purposes and available sizes only. This listing is for only 1 mattress in the listed size.*

[Read More](#)

SKU: 000.936.0007



HTP 1500 * 20" X 44" REUSABLE LARGE WATER MATTRESS

The **Adroit Medical Systems HTP 1500** is a reusable large water mattress measuring **20 inches by 44 inches**, designed for effective patient warming and pressure redistribution. Crafted for durability and convenience, this water mattress provides consistent, adjustable heat therapy in medical and surgical settings. Its robust construction ensures reliable performance and easy reusability, making it an essential component for patient comfort and safety during treatments.

**Image is for illustration purposes and available sizes only. This listing is for only 1 mattress in the listed size.*

[Read More](#)

SKU: 000.936.0006



HTP 1500 HEAT THERAPY PUMP, NEW

The **HTP 1500 Heat Therapy Pump** by Adroit Medical Systems is a highly precise and reliable device designed for effective heat therapy in medical settings. As a veteran-owned small business, Adroit manufactures a range of patient warming products and is committed to quality and safety, with FDA registration, UL listing, CE marking, and ISO 13485:2003 and ISO 9001:2008 certifications. The HTP 1500 offers set point accuracy within 1.0°F, with a temperature range from 75 to 107°F, and a flow rate of 15 GPM to ensure consistent heat delivery. It features a 50 oz reservoir, hospital-grade 8-foot power cord, uses tap water—mercury-free—and operates at 120 VAC, 2 amps, 180 watts. The device includes a real-time display for easy monitoring, comes with one medium heat pad, and is backed by a 1-year factory warranty, making it an ideal solution for targeted heat therapy applications.

[Read More](#)

SKU: 000.936.0016



MD5040V UPS SYSTEM DELIVERS 5000 WATTS OF SINGLE-PHASE POWER WITH A 4.8 KWH BATTERY BANK

The **MD5040V UPS System** provides a robust 5,000-watt single-phase power solution, equipped with a 4.8 kWh AGM lead-acid battery bank to ensure reliable backup during power outages. Designed for high-capacity performance, it features a pure sine wave output at 120V and a surge capacity of up to 10,000 watts, with a transfer time of just 20 milliseconds to maintain seamless power supply. The system connects via a 10-foot power cord (NEMA L6-30P) and delivers four duplex NEMA 5-20R outlets protected by circuit breakers, supporting an output load of up to 20 amps. Built in a stainless steel enclosure, the UPS has UL listings for both inverter (UL 1741) and battery components (UL 94V-0), and weighs approximately 429 pounds. Additional features include optional power outage alerts, remote battery monitoring, and alarm dry contacts, making it an ideal solution for critical applications that demand uninterrupted power and system management.

[Read More](#)

SKU: AM-MD5040V



INBODY 970S BODY COMPOSITION ANALYZER

The **InBody 970S Body Composition Analyzer** is a high-precision medical and fitness device designed for professionals seeking the most accurate assessments of body composition. It offers advanced measurements—including visceral fat levels, Skeletal Muscle Mass Index (SMI), and detailed segmental analysis—going beyond standard outputs to provide comprehensive insights into body health. The package includes the InBody 970S analyzer, a USB thumb drive for data transfer, a \$300 LookinBody Web credit (expires 90 days after receipt if unused), a box of 300 InBody tissues, three detailed posters, an interpretation poster to help understand results, a laser printer (Brother HL-L2400D) for printing reports, and dedicated result sheets for body composition and water analysis. This setup supports accurate, efficient evaluation for clinical, sports, or wellness applications, empowering professionals with actionable body composition data.

[Read More](#)

SKU: AM-InBody 970S

INDEX

H

HTP 1500 * 20" X 44" REUSABLE LARGE WATER MATTRESS 2

HTP 1500 20" X 57" REUSABLE LARGE WATER MATTRESS 2

HTP 1500 HEAT THERAPY PUMP, NEW 3

I

InBody 970S Body Composition Analyzer 4

M

MD5040V UPS system delivers 5000 watts of Single-phase power with a 4.8 KWH battery bank 3