

TABLE OF CONTENTS

Infusion and Syringe Pumps	3
Medfusion 3500 V 3.0.6 Syringe Pump	4
Alaris Medley 8015 Point of Care Module (Software Version 9.33)	4
Alaris 8120 PCA Syringe Pump	5
Eitan Medical Sapphire Epidural Infusion Pump (New)	7
Eitan Medical Sapphire™ Multi-Therapy Infusion Pump (New)	8
Eitan Medical Sapphire™ Multi-Therapy Infusion Pump with Lockbox/Cradle (Refurbished)	8
Medfusion 4000 Syringe Pump	9
Medfusion 3500 V6 Syringe Pump	10
Baxa Repeater Pump	11
B. Braun Infusomat Space Infusion Pump Wifi Capable, Wifi Batteries (NEW)	12
Foot Pedal for Baxa Repeater Pump	12
Eitan Medical Sapphire™ Multi-Therapy Infusion Pump (Refurbished)	13
B. Braun Infusomat Space Infusion Pump Wifi Capable, Wifi Batteries (Refurbished)	14
Alaris 8300 EtCO2 Module Only	14
Medfusion 3500 Syringe Pump with Version 5 Software	15
Volumat MC Agilia Infusion Pump	16
BodyGuard 323 CME	17
Alaris IVAC MedSystem III 2865B Infusion Pump	17
Eitan Medical Sapphire Epidural Infusion Pump (Refurbished)	18
CMExpress Infusion Pump	19
Plum 360 Infusystem	19
Medfusion 2010i Syringe Pump	20
Alaris Medley 8110 Syringe Pump Module	20
Medfusion 3500 V4 Syringe Infusion Pump	21
Baxter AS50 Syringe Pump	22
Walkmed 350VL Infusion System	23
Hospira Lifecare PCA3	24
Curlin 6000 CMS	24
Curlin 4000 CMS	25
B. Braun Infusomat Space	25
CADD Legacy Plus Pump	26

Baxter Sigma Spectrum Infusion Pump With 6.02 Software	27
Baxter PCA II	27
B. Braun Perfusor Space Syringe Pump	28
Baxter 6301 Infusion Pump	29
Alaris MedSystem 2865	29
B. Braun Outlook 100	30
Alaris 8000 PC Point-of-Care Unit	31
Abbott Triple Channel Plum A+ IV	31
Abbott Plum A+ IV Infusion Pump, Single Channel	32
Alaris Medley 8100 Infusion Pump	33
Alaris Medley 8015 Point of Care Module	33
Baxter Flo-Gard 6201	34
B. Braun Vista Basic Infusion Pump	35
Baxter InfusO.R. Syringe Pump W/ Propofol Smart Label	36



INFUSION AND SYRINGE PUMPS

Auxo Medical understands the importance of choosing the best infusion pumps for your facility. Auxo Medical provides infusion pumps from leading brands such as [Braun](#), [Baxter](#), and [Alaris](#).

These pumps help provide the best in patient care for hospitals, nursing homes, and at-home care situations. In addition, they focus on accurate fluid and drug delivery requirements of today's medical environment. Furthermore, many of our infusion pump options offer both continuous and intermittent infusion options. Built for drug solutions, whole blood, and packed red blood cells, our infusion machines cover a variety of intravenous, intra-arterial, epidural, and subcutaneous options for administration.

Auxo Medical's infusion pumps make it easy for you to deliver top quality patient care and safety. Our selection of innovative infusion units help maintain safety and workflow within any facility by reducing medication errors and streamlining your medication management.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the

equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



MEDFUSION 3500 V 3.0.6 SYRINGE PUMP

The **Medfusion 3500 Syringe Pump with Version 3.0.6 Software** is a versatile, advanced infusion device suitable for use across various hospital departments and physician practices. This refurbished "smart pump" features automatic detection of multiple syringe types from seven different manufacturers, ensuring enhanced safety and precise drug delivery. Designed for ergonomic single-hand loading and easy mounting on IV poles or bassinets, it utilizes FlowSentry technology to promote patient safety and medication accuracy. Capable of infusing at flow rates from 0.01 mL/hr up to 1130 mL/hr, it supports low-flow infusions with larger syringes and is MRI compatible. Sold refurbished with a 6-month warranty, the Medfusion 3500 offers reliable, accurate, and flexible syringe infusion management.

Warranty: 6 Months

[Read More](#)

SKU: AM-MF3500-V3

Price: \$1,795.00



ALARIS MEDLEY 8015 POINT OF CARE MODULE (SOFTWARE VERSION 9.33)

The **Alaris Medley 8015 Point of Care Module (Software Version 9.33)** is a versatile, advanced component of the Alaris® System designed to support adult, pediatric, and neonatal patient care. It features a unified user interface for programming infusion and monitoring functions, enhancing clinical efficiency and safety. The system includes the Guardrails® Suite MX safety software, which provides double-checks of dose, duration, and delivery rate for every infusion, along with wireless connectivity for seamless communication and real-time analytics. Compatible with the Alaris Medley 8100 Infusion Pump, the module offers exceptional versatility, increased IV safety, and robust protection features. Its compact dimensions are 6.9" W x 8.8" H x 9" D with a pole clamp, and it weighs approximately 7.2 pounds. The system's internal memory ensures configuration data retention for at least six months, with module-specific settings stored for 8 hours when powered off. Powering the system requires 100-240V at 50/60 Hz, with a maximum power draw of 150VA.

Warranty: 6 Months

[Compatible with the Alaris Medley 8100 Infusion Pump.](#)

[Read More](#)

SKU: AM-Alaris-8015-9.33-D

Price: \$950.00



ALARIS 8120 PCA SYRINGE PUMP

The **CareFusion Alaris 8120 PCA Syringe Pump** is an advanced infusion device designed for safe and effective patient-controlled analgesia (PCA) management. Integrated within the Alaris® System, it combines programming safety features, barcoding capabilities, and continuous respiratory monitoring into a comprehensive pain control solution. The PCA module offers quick access to total dose history and trending data, enabling healthcare providers to analyze patient responses and optimize pain therapy. Notably, it can automatically pause PCA infusions and deactivate dose requests if the patient's oxygen saturation (SpO₂) or exhaled CO₂ (EtCO₂) levels fall below hospital-defined thresholds when used with the PCA module. The pump allows for different dose ranges for the same drug based on patient diagnosis (e.g., acute vs. chronic pain) and supports barcoding to ensure the accuracy of drug, patient, and clinician identification. Its intuitive interface simplifies programming, reducing training time, while the Guardrails® Suite MX software enables customizable alarm limits and trend monitoring. When combined with the Auto-ID module, it provides reliable barcoding for the right drug, patient, and clinician, enhancing safety and workflow efficiency in pain management.

Warranty: 6 Months

[Read More](#)

SKU: AM-Alaris-8120

Price: \$2,261.00



EITAN MEDICAL SAPPHIRE EPIDURAL INFUSION PUMP (NEW)

This professionally refurbished **Eitan Medical Sapphire Epidural Infusion Pump** is a compact, reliable, and precise device designed specifically for epidural infusion therapy. Weighing approximately 14.7 ounces and measuring about 5.63" high, 3.78" wide, and 1.93" deep, this pump offers seamless portability for use in diverse clinical environments. It features a broad volume capacity of 1 to 9999 mL with high delivery accuracy of $\pm 2.5\%$, supporting flow rates from 0.1 to 99.9 mL/hr (increments of 0.1 mL/hr) and higher rates up to 999 mL/hr. The Sapphire epidural pump is powered by a standard 100-240 V supply, with an impressive battery life of 24 hours at 125 mL/hr, and includes safety alarms for occlusions, cassette issues, flow errors, and insufficient battery to ensure safe operation. Its design incorporates essential accessories such as an AC charger, bolus cord, and mini cradle with pole clamp, making it an efficient, user-friendly solution for delivering precise epidural therapy with confidence.

Warranty: 12 Months

[Read More](#)

SKU: AM-SAPPHIRE-E-LC-N

Price: \$3,695.00



EITAN MEDICAL SAPPHIRE™ MULTI-THERAPY INFUSION PUMP (NEW)

The **Sapphire™ Multi-Therapy Infusion Pump** (OEM: EIT17000-028-0072 / 17000-028-1039) is a versatile, compact, and user-friendly device engineered to deliver a broad range of intravenous therapies with precision across various healthcare settings, including hospitals, ambulances, and home care. Supporting multiple infusion modes such as continuous, intermittent, total parenteral nutrition (TPN), patient-controlled analgesia (PCA), and epidural, it offers unparalleled flexibility to meet diverse patient needs. Renowned for its reliability and accurate dosage delivery, the Sapphire™ pump features an intuitive interface, robust safety alarms—including occlusion detection, cassette misplacement, and flow errors—and drug library access that enhance patient safety and caregiver confidence. With a lightweight design weighing approximately 14.7 ounces and a battery life of approximately 24 hours at 125 mL/hr, it ensures portability and uninterrupted therapy. Its precise flow rate capabilities (0.1 to 99.9 mL/hr with 0.1 mL/hr increments), a broad volume range up to 9999 mL, and a delivery accuracy of $\pm 2.5\%$ make it an essential tool for delivering complex therapies efficiently, improving patient outcomes, and streamlining infusion management.

Warranty: 12 Months

[Read More](#)

SKU: AM-SAPPHIRE-N

Price: \$2,795.00



EITAN MEDICAL SAPPHIRE™ MULTI-THERAPY INFUSION PUMP WITH LOCKBOX/CRADLE (REFURBISHED)

This professionally refurbished **Eitan Medical Sapphire™ Multi-Therapy Infusion Pump with Lockbox/Cradle** is a sophisticated and versatile device designed to provide precise intravenous therapies across diverse clinical settings, including hospitals, ambulances, and home care. Compact and lightweight at approximately 5.63" high, 3.78" wide, and 1.93" deep, and weighing about 14.7 ounces without the battery, this pump offers a broad volume capacity of 1 to 9999 mL with a delivery accuracy of $\pm 2.5\%$. Supporting multiple infusion modes such as continuous, intermittent, total parenteral nutrition (TPN), patient-controlled analgesia (PCA), and epidural, it ensures safe and accurate medication or nutrition administration. Battery life extends up to 24 hours at 125 mL/hr, and safety features include alarms for occlusions, cassette misplacement, flow errors, and low battery, with occlusion alarms capable of withstanding pressures up to 17.4 PSI. The device's intuitive interface, combined with a lockbox and cradle for secure placement and easy transport, makes it an essential tool for delivering complex therapies reliably and safely, helping healthcare professionals improve patient outcomes and streamline infusion management.

Warranty: 6 Months

[Read More](#)

SKU: AM-SAPPHIRE-MS-LC

Price: \$2,724.75



MEDFUSION 4000 SYRINGE PUMP

The **Medfusion 4000 Syringe Pump** is a leading infusion device designed to deliver safe, precise medication administration across a wide range of clinical settings, including critical care, anesthesia, neonatal, and pediatric applications. Building on the trusted platform of the widely used Medfusion 3500, the 4000 offers enhanced pump functionality and an improved drug library, supported by PharmGuard Medication Safety Software and Server Software to ensure safe and straightforward operation. Its advanced features include accurate delivery within $\pm 2\%$ across flow rates from 0.01 mL/hr to 1130 mL/hr, and proprietary smart pump technology that promotes configurable drug libraries, making infusion more reliable than traditional methods. The ergonomic, horizontal design allows clinicians to load syringes easily with one hand, while FlowSentry pressure sensing technology enables early detection of infusion issues for timely clinical intervention. The pump supports both wireless and wired data collection, report generation, and seamless device updates, simplifying quality improvement initiatives. It also boasts extended memory for comprehensive historical data storage, IPX3 fluid resistance, and user-friendly Windows-based software for efficient data reporting. Additionally, PharmGuard Server Software manages drug library configurations, firmware updates, and generates detailed clinical, IT, and biomedical reports, all employing industry-standard protocols and databases. Complete with a pole clamp and AC adapter, the Medfusion 4000 is a versatile, reliable solution for delivering critical medications safely and accurately in any hospital environment.

Warranty: 6 Month

[Read More](#)

SKU: AM-MF4000

Price: \$2,995.00



MEDFUSION 3500 V6 SYRINGE PUMP

The **Medfusion 3500 V6 Syringe Pump**, enhanced by PharmGuard Medication Safety Software, offers an intuitive and safe medication delivery system. Its large screen displays up to 30 characters per drug name, minimizing abbreviation errors, and the bright LCD backlight ensures clear visibility in any lighting. The pump's extensive configurability allows for customizable clinical alerts, promoting best practices. Features include programmable weight, volume, and occlusion limits; support for loading and bolus dose therapies; adjustable alarm settings; and FlowSentry rapid occlusion detection for continuous pressure monitoring. The horizontal syringe plunger simplifies loading and unloading, protecting the barrel. A 10-hour battery life ensures reliable operation. The Medfusion 3500 V6 combines advanced features with ease of use to improve medication safety and efficiency.

Warranty: 6 Months

[Download Brochure](#)

[Read More](#)

SKU: AM-MF3500-V6

Price: \$2,195.00



BAXA REPEATER PUMP

The **Baxa Repeater Compounding Pump** is a versatile Peristaltic pump which has played an important role in routine pharmacy procedures - delivering accuracy and multi-speed settings for fluid transfer processes. Designed to reduce the workload of pharmacy technicians, the Baxa Repeater Pump facilitates the automated fluid transfer.

With compact dimensions, the Repeater Pump offers flexibility in placement options - fit it within a standard hood or mount it on the available pump stand. Accurate to volumes as low as 0.2 mL ($\pm 10\%$), the Repeater Pump was designed to transfer fluids from 0.2 mL to 9,999 mL at flow rates up to 13.5 mL per second with a 16G needle. Designed for quick set up to help streamline pharmacy operations, the versatile Repeater Pump can support any size pharmacy, handling a wide range of fluid transfer and filling applications with efficiency and accuracy.

(Commonly referred to as Baxter Repeater Pump)

[Foot Pedal Sold Separately Here: Foot Pedal for Baxa Repeater Pump](#)



6 Month Warranty Included

\$4995.00

On Backorder: Current Lead Time: 16-18 Weeks

To be added to the waitlist submit a quote request below.

[Read More](#)

SKU: AM-B-REPEATER

Price: Request Quote



B. BRAUN INFUSOMAT SPACE INFUSION PUMP WIFI CAPABLE, WIFI BATTERIES (NEW)

The **B. Braun Infusomat Space Infusion Pump** is designed for acute care adult and pediatric facilities, the compact 9" x 6" Infusomat® Space large volume pump weighs just 3 pounds. The lightweight large volume pump allows for flexibility and easy transport within a number of care settings. Along with industry-leading safety features, the lightweight mobile design allows you to quickly and easily stack up to 24 pumps in a single bedside tower. Equipped with an extensive drug library, easy navigation and intuitive interface, this small infusion pump packs state-of-the-art versatility that promotes efficiencies across numerous care settings. **Model Ref: 8713051U**

Warranty: 12 Months

[Read More](#)

SKU: AM-BRAUN-INFUSOMAT-S-W-N

Price: \$2,495.00



FOOT PEDAL FOR BAXA REPEATER PUMP

The **Foot Pedal for the Baxa Repeater Pump** is a pneumatic foot-operated control device designed to enable hands-free operation of the Baxa Repeater Pump. This pedal allows clinicians, pharmacists, and technicians to activate or control the pump's functions conveniently and hygienically, helping to maintain sterile conditions and improve workflow efficiency during infusion procedures. Its simple pneumatic design ensures durability and reliable performance in clinical and pharmaceutical environments.

[Read More](#)

SKU: AM-IP326NN

Price: \$109.95



EITAN MEDICAL SAPPHIRE™ MULTI-THERAPY INFUSION PUMP (REFURBISHED)

This professionally refurbished **Sapphire™ Multi-Therapy Infusion Pump** (OEM: EIT17000-028-0072 / 17000-028-1039) is a versatile, compact, and user-friendly device engineered to deliver a broad range of intravenous therapies with precision across various healthcare settings, including hospitals, ambulances, and home care. Supporting multiple infusion modes such as continuous, intermittent, total parenteral nutrition (TPN), patient-controlled analgesia (PCA), and epidural, it offers unparalleled flexibility to meet diverse patient needs. Renowned for its reliability and accurate dosage delivery, the Sapphire™ pump features an intuitive interface, robust safety alarms—including occlusion detection, cassette misplacement, and flow errors—and drug library access that enhance patient safety and caregiver confidence. With a lightweight design weighing approximately 14.7 ounces and a battery life of approximately 24 hours at 125 mL/hr, it ensures portability and uninterrupted therapy. Its precise flow rate capabilities (0.1 to 99.9 mL/hr with 0.1 mL/hr increments), a broad volume range up to 9999 mL, and a delivery accuracy of $\pm 2.5\%$ make it an essential tool for delivering complex therapies efficiently, improving patient outcomes, and streamlining infusion management.

Warranty: 6 Months

Lead Time: 8-12 Business Days

[Read More](#)

SKU: AM-SAPPHIRE

Price: \$1,595.00

B. BRAUN INFUSOMAT SPACE INFUSION PUMP WIFI CAPABLE, WIFI BATTERIES (REFURBISHED)

This professionally refurbished **B. Braun Infusomat Space Infusion Pump** is a compact, lightweight, and versatile large-volume infusion system designed for adult and pediatric acute care settings. Measuring just 9" x 6" and weighing only 3 pounds, this pump offers exceptional mobility and ease of transport across various clinical environments. Its user-friendly interface and extensive drug library simplify programming, while audible and visual alarms enhance safety and clarity during critical procedures. The pump's innovative design allows clinicians to stack up to 24 units in a single bedside tower, optimizing space management. Available with or without wireless capabilities, the Infusomat Space features a strong pole clamp capable of securely mounting up to three pumps either horizontally or vertically, providing flexibility and efficiency in busy care areas. Its combination of safety, ease of use, and modularity makes it an ideal infusion solution for modern healthcare facilities.

[Download Brochure](#)**6 Month Warranty**[Read More](#)**SKU:** AM-BRAUN-INFUSOMAT-S-W**Price:** \$1,695.00

ALARIS 8300 ETCO2 MODULE ONLY

The **Alaris 8300 EtCO₂ Module Only** is a noninvasive, continuous respiratory monitoring device designed to measure exhaled carbon dioxide (EtCO₂) and respiratory rate in both intubated and non-intubated patients across all care settings. Incorporating Oridion Microstream technology with Filterline disposables, it provides real-time, accurate respiratory status updates, including trend data and breath alarms, helping to reduce risks associated with opioid administration and respiratory compromise. Its user-friendly design features customizable profiles for versatile application and utilizes sampling lines that eliminate moving parts, external sensors, or routine calibration, ensuring reliable performance with minimal maintenance. Compact at 3.3" x 8.9" x 5.5" and offering up to 5.5 hours of internal power, [**the module requires an Alaris 8015 Point of Care Module to operate**](#) and supports broad patient sampling needs, making it an essential tool for enhancing respiratory safety and clinical decision-making.

Warranty: 6 Months

[Read More](#)

SKU: AM-ALARIS-8300

Price: \$1,586.05





MEDFUSION 3500 SYRINGE PUMP WITH VERSION 5 SOFTWARE

The **Medfusion 3500 Syringe Pump** is a cutting-edge device designed to meet the needs of various departments within hospitals and medical practices. This refurbished pump features intuitive syringe detection capabilities, enhancing safety and accuracy in drug administration. Known as a "smart pump," it accommodates syringes from seven different manufacturers, ensuring versatility for healthcare providers. Its ergonomic design allows for single-hand loading, making it easy to operate, while its mounting options enable secure placement on IV poles or bassinets. Additionally, the Medfusion 3500 employs **FlowSentry technology**, which further safeguards patient safety and ensures precise drug delivery. With an accuracy of $\pm 2\%$ and compatibility with MRI settings, this advanced syringe pump is an essential tool for improving medication administration in any clinical setting. The inclusion of Version 5 Software enhances its functionality, providing healthcare professionals with a reliable and efficient solution for intravenous therapy.

Warranty: 6 Months

[Read More](#)

SKU: AM-MF3500-V5

Price: \$1,795.00



VOLUMAT MC AGILIA INFUSION PUMP

The **Volumat MC Agilia Infusion Pump** is a highly versatile and precise infusion device designed to accommodate a wide range of clinical needs, from routine fluid delivery to complex medication protocols. It features a flow rate range of 1–1500 mL/h in normal mode with 1 mL/h increments, and 0.1–99.9 mL/h in micro mode with 0.1 mL/h increments, with the ability to set soft and hard limits for drug safety via Vigilant®, the IV Medication Safety Solution. Its extensive infusion modes include volumetric, dose-based, ramp-up/down, intermittent, secondary, and drop/min, with advanced programming options such as drug labeling and customized drug libraries. The pump is equipped with the Occlusivity Check System (OCS) for auto-testing infusion line integrity, dynamic pressure warnings through DPS, and dual pressure modes to enhance safety. Its manual or programmable bolus feature provides up to 1500 mL/h, adjustable in 50 mL/h steps, while the rechargeable Li-Ion battery offers at least 8 hours of continuous operation. Compact at approximately 145 x 130 x 135 mm and weighing around 2 kg, the Volumat MC Agilia ensures accurate, safe, and flexible infusion therapy management across diverse healthcare environments.

Warranty: 12 Months

[Read More](#)

SKU: AM-Volumat-MCA-N

Price: \$2,995.00

BODYGUARD 323 CME

The **BodyGuard 323 CME** is a versatile, state-of-the-art infusion pump designed to provide safe, efficient, and flexible medication delivery across a variety of clinical settings. Its six programmable modes—Continuous, 25-Step, TPN, Intermittent Sinusoidal, and PCA—allow clinicians to administer antibiotics, chemotherapy, nutrition, and analgesics with precision. Compact and lightweight, weighing just over 11 ounces without the battery, the pump is ideal for hospital, homecare, transport, and emergency use, offering flexible configurations including pole mounting, carrying pouches, and handheld computer connectivity with customizable drug libraries. Integrated features include a barcode reader, pulse oximeter attachment, patient temperature, and heart rate monitoring, all powered by included AC and standard batteries, making the BodyGuard 323 CME a reliable, multi-therapy infusion solution for diverse patient care environments.

Warranty: 6 Months

[Read More](#)

SKU: AM-CME-BG-323

Price: \$1,295.00





ALARIS IVAC MEDSYSTEM III 2865B INFUSION PUMP

The **Alaris IVAC MedSystem III 2865B** is a compact, multi-channel infusion pump designed to deliver fluids accurately and efficiently across various clinical settings. Featuring three independent fluid delivery systems within a small footprint, it helps reduce bedside clutter and simplifies patient transport. The pump offers six selectable device modes—general purpose, neonatal, controller pressure, O.R., general purpose II, and O.R. II—ensuring versatile performance. Equipped with a drug dose calculator for precise medication management, this lightweight unit (approximately 5.1 lbs.) measures 7.875 inches high, 6 inches wide, and 2.10 inches deep, combining ease of use with advanced functionality to optimize infusion therapy.

Warranty: 6 Months

[Read More](#)

SKU: AM-IVAC-2865B

Price: \$1,695.00



EITAN MEDICAL SAPPHIRE EPIDURAL INFUSION PUMP (REFURBISHED)

This professionally refurbished **Eitan Medical Sapphire Epidural Infusion Pump** is a compact, reliable, and precise device designed specifically for epidural infusion therapy. Weighing approximately 14.7 ounces and measuring about 5.63" high, 3.78" wide, and 1.93" deep, this pump offers seamless portability for use in diverse clinical environments. It features a broad volume capacity of 1 to 9999 mL with high delivery accuracy of $\pm 2.5\%$, supporting flow rates from 0.1 to 99.9 mL/hr (increments of 0.1 mL/hr) and higher rates up to 999 mL/hr. The Sapphire epidural pump is powered by a standard 100-240 V supply, with an impressive battery life of 24 hours at 125 mL/hr, and includes safety alarms for occlusions, cassette issues, flow errors, and insufficient battery to ensure safe operation. Its design incorporates essential accessories such as an AC charger, bolus cord, and mini cradle with pole clamp, making it an efficient, user-friendly solution for delivering precise epidural therapy with confidence.

Warranty: 12 Months

[Read More](#)

SKU: AM-SAPPHIRE-E-LC

Price: \$1,833.15



CMEXPRESS INFUSION PUMP

The **CMExpress Infusion Pump** is a user-friendly device designed for precise delivery of continuous and multi-step infusions, making it ideal for busy oncology infusion suites and outpatient clinics. Its straightforward programming and cost-effective sets streamline setup, while flow rates up to 1200 ml/hour accommodate a wide range of oncology protocols. Supporting both continuous and piggyback infusions, the pump automatically cycles through key therapy information for improved safety and efficiency. Equipped with a pole mount charger that features a large, easily visible display, and backed by an automatic switch to battery power during outages or patient ambulation, the CMExpress provides reliable, flexible infusion management. Included components are the pole mount charger and standard battery for comprehensive use.

Warranty: 6 Months

[Read More](#)

SKU: AM-CME-EX

Price: \$2,705.30



PLUM 360 INFUSYSTEM

The **Plum 360 Infusystem** is a versatile and highly efficient infusion pump designed to meet a wide range of clinical needs. With dimensions of approximately 8 inches in height, width, and 6 inches in depth, and weighing around 10 lbs (4.5 kg) with the battery, it is compact and portable. It operates on a battery life of roughly 3 hours at rates up to 125 mL/hr on a single line, or it can deliver at least 250 mL if the rate exceeds 126 mL/hr. The device requires a power input of 120 V~, 50 to 60 Hz, and 50 VA, with an approximate recharge time of 8 hours while running at 125 mL/hr. The volume range extends from 0.1 to 9999 mL, allowing for precise delivery in 0.1 mL and 1 mL increments. Delivery rate accuracy is maintained within $\pm 5\%$ for rates between 1 to 999 mL/hr and $\pm 10\%$ for rates between 0.1 to 0.9 mL/hr. Its dual-line system can deliver between 0.1 to 999 mL/hr per line, with a concurrent delivery minimum of 0.5 mL/hr per line. The PlumSet can manage a cumulative delivery rate of up to 500 mL/hr for both lines combined. The system includes several alarms for occlusion, air-in-line, low battery, door opened while pumping, lockout violation, VTBI complete, valve/cassette test failure, and more, ensuring safety and reliability. The occlusion alarm pressure can be set between 1-15 psig, with a default limit of 6 psig. This comprehensive set of features

makes the Plum 360 Infusystem a dependable choice for healthcare providers.

Warranty: 6 Months

[Read More](#)

SKU: AM-PLUM-360

Price: \$1,495.00

MEDFUSION 2010I SYRINGE PUMP

The **Medfusion 2010i Syringe Pump** is a highly versatile device designed for precise infusion of blood, lipids, fluids, antibiotics, and other medications, including epidural infusions. It offers 15 different delivery modes, including Body Weight, Mass, Continuous, and Volume over Time, providing flexible options for various clinical needs. Users can program boluses to administer at any point before or during drug delivery, and the device can be customized for specific applications in the field. With the ability to store up to 64 user-defined protocols, the Medfusion 2010i supports continuous, auto or manual intermittent, and volume/time infusions. Supplied with a pole clamp and power cord, it ensures reliable, adaptable infusion management.

Warranty: 6 Months

[Read More](#)

SKU: AM-2010i

Price: \$1,405.00



ALARIS MEDLEY 8110 SYRINGE PUMP MODULE

The **Alaris Medley 8110 Syringe Pump Module** is a precise and reliable infusion device designed to work seamlessly with the Alaris PC unit. It incorporates sophisticated Guardrails® Suite MX software to ensure accurate medication dosing and infusion rates, enhancing patient safety. Ideal for delivering fluids, medications, blood, or blood products via continuous or intermittent routes such as intravenous, intra-arterial, or enteral, this syringe module features a unique pressure sensing disc that quickly detects occlusions, reduces alarm delays, and accurately monitors in-line pressure. Its advanced design minimizes inadvertent bolus administration following occlusion release, providing a safer, more efficient infusion process in clinical settings.

Warranty: 6 Months

[Read More](#)

SKU: AM-Alaris-8110

Price: \$1,895.00





MEDFUSION 3500 V4 SYRINGE INFUSION PUMP

The **Medfusion 3500 V4 Syringe Infusion Pump** is a state-of-the-art device tailored for use in various hospital departments and physician practices. This refurbished syringe pump is engineered to intuitively detect a range of syringes, enhancing both safety and accuracy in drug administration. Recognized as a "smart pump," it accommodates syringes from seven different manufacturers and features an ergonomic design that allows for single-hand loading, making it user-friendly. The Medfusion 3500 can be easily mounted to an IV pole or infant isolette, ensuring versatility in clinical settings. Utilizing **FlowSentry technology**, this pump prioritizes patient safety and medication accuracy, providing healthcare professionals with a reliable infusion solution. Sold refurbished, it includes a **6-month warranty**, giving users peace of mind regarding their investment. With an accuracy of +/- 2% and compatibility with MRI environments, the Medfusion 3500 V4 is an essential addition to any medical facility's equipment lineup.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-MF3500-V4

Price: \$1,595.00



BAXTER AS50 SYRINGE PUMP

The **Baxter AS50 Syringe Pump** is a versatile, compact infusion device designed to accurately deliver fluids, medications, and blood products via various routes—including intravenous, intra-arterial, epidural, or subcutaneous. Equipped with digital controls and flexible drug libraries, it accommodates standard disposable syringes ranging from 1mL to 60mL, ensuring precise continuous or intermittent infusions. Its features include a flow rate accuracy of $\pm 3\%$ and volume accuracy of $\pm 3\%$, with a broad flow rate range from 0.01 mL/hr up to 438 mL/hr depending on syringe size. With approximate dimensions of 3.4" W x 2.6" D x 10" H and weighing about 3.2 lbs, the unit is portable and adaptable. Additional options include an IV pole with caster wheels and hooks for convenient integration into clinical workflows, making it suitable for a variety of infusion needs in today's dynamic healthcare environment.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-AS50

Price: \$1,794.00



WALKMED 350VL INFUSION SYSTEM

The **WalkMed 350VL Infusion System** is a highly reliable ambulatory infusion pump designed for precise delivery of medications and fluids through various routes, including intravenous, subcutaneous, arterial, enteral, and epidural. Renowned for its trouble-free performance and accurate low-rate infusion capabilities, it is particularly well-suited for chemotherapy treatments. Features such as patient lockout, volume tracking, auto-occlusion alarm recovery, and a symbolic keypad ensure safety and ease of operation. Compact and lightweight—measuring 4.4 x 4.0 x 1.5 inches and weighing just 12.8 oz.—the system is intended for use in controlled settings, with contraindications including blood, blood products, insulin, and critical medications where interruption could lead to serious injury or death. Its programmable total volume and delivery tracking, combined with safety features like keypad lockout and alarm silencing, make the WalkMed 350VL a dependable choice for precise ambulatory infusions.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-WM-350VL

Price: \$1,731.00



HOSPIRA LIFECARE PCA3

The **Hospira Lifecare PCA3** is an advanced patient-controlled analgesia (PCA) system designed to improve pain management safety and efficiency, especially when integrated with Hospira MedNet software. It features customizable drug libraries, an integral barcode reader for medication verification, and support for multiple clinical care areas with up to 18 medication rule sets. The system's intuitive, human factors-designed interface includes confirmation screens, stored protocols for quick programming, and pre-filled, color-coded bar-coded medication vials to reduce errors. Wireless connectivity enables real-time infusion monitoring and support for 5-rights verification, ensuring accurate medication delivery. Compact at 8"W x 13"H x 6"D, weighing approximately 10 lbs. with a rechargeable 8V lead-acid battery offering around 4 hours of continuous use, the Lifecare PCA3 combines safety, flexibility, and ease of use for effective analgesia therapy.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-H-PCA3

Price: \$776.00



CURLIN 6000 CMS

The **Curlin 6000 Clinical Management System (CMS™)** is a versatile, multi-therapy ambulatory infusion device designed to deliver a wide range of infusion modes, including PCA, PCEA, subcutaneous, TPN, intermittent, and variable therapies. Its advanced curvilinear peristaltic pump mechanism ensures precise fluid delivery while continuously monitoring line pressure to help detect potential catheter issues such as occlusions, obstructions, or pinched lines. Powered by alkaline batteries or rechargeable external packs, it offers exceptional battery life—up to 10 hours at 400 mL/hr, 30 hours at 125 mL/hr, and 85 hours at 2 mL/hr—making it ideal for complex clinical environments. The system features multiple alarms for critical safety alerts, adjustable alarm volume, and convenient review functions including Immediate On-Demand (IOD) data access for therapy verification. Compact, lightweight at 18.1 oz, and equipped with a clear, backlit display, the Curlin 6000 CMS provides reliable, continuous infusion management with flexible flow rates up to 400 mL/hr and adjustable volume limits for tailored patient care.

[Download Brochure](#)

Warranty: 6 Months[Read More](#)**SKU:** AM-CURLIN-6000CMS**Price:** \$1,895.00

CURLIN 4000 CMS

The **Curlin 4000 CMS** is an advanced, multi-therapy ambulatory infusion pump designed for precise medication delivery across various clinical applications, including continuous, PCA (IV, epidural, subcutaneous), TPN with ramping, intermittent, and variable dosing with up to 24 preset doses. Built with non-DEHP, PVC-free administration sets, it offers superior flow resolution—up to 25 times better at low rates—enabling highly accurate infusion even at minimal flow levels. Its standout feature is real-time line pressure monitoring, aiding in early detection of catheter patency issues such as occlusions or obstructions. The system's configurable alarms, including infusion completion, occlusion alerts, and low battery notifications, enhance safety, while its compact size (5"H x 4"W x 2.5"D) and lightweight design (17.5 oz) ensure portability. Powered by alkaline batteries or rechargeable external packs, the Curlin 4000 CMS delivers an exceptional combination of precision, versatility, and safety, optimizing patient outcomes in diverse clinical settings.

[Download Brochure](#)**Warranty: 6 Months**[Read More](#)**SKU:** AM-CURLIN-4000CMS**Price:** \$1,396.00



B. BRAUN INFUSOMAT SPACE

The **B. Braun Infusomat Space** is a compact, lightweight, large-volume infusion pump designed for both adult and pediatric acute care environments. Measuring just 9" x 6" and weighing approximately 1.4 kg (3 lbs), it offers easy mobility and the ability to stack up to 24 units for streamlined bedside management. The pump features an intuitive user interface with audible and visible differentiated alarms, a strong pole clamp for secure mounting, and an extensive 1,500-entry drug library for quick access and safety. Its single-channel design includes KeyGuard™ for preventing programming errors and supports wireless integration for interoperability across care systems. The backlit graphic display provides a clear view at a ~40° angle, while flow rates range from 0.01 to 1800 mL/h. Powered by an external Space Power Supply or SpaceStation, it offers impressive operational times—up to 16.8 hours at 1 mL/h—making it a versatile, efficient infusion solution in diverse clinical settings.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-BRAUN-INFUSOMAT-S

Price: \$1,195.00



CADD LEGACY PLUS PUMP

The **Legacy Plus Model 6500** is the latest addition to the esteemed CADD family, renowned for its therapy-specific delivery systems. This innovative pump sets a new benchmark for continuous and intermittent ambulatory infusion, particularly in research settings. CADD pumps have enjoyed a stellar reputation for accuracy and durability since their inception in the early 1980s, and the CADD-Legacy PLUS model continues this tradition. Designed for both continuous and intermittent therapies, the pump comes with a single-use 50-/100-ml pump pouch, 2 AA alkaline batteries, and an operator's manual. Among its impressive features are an easy-to-read display for straightforward programming and troubleshooting, a cassette detection sensor, and occlusion sensors to alert users of any interruptions in fluid delivery. Additionally, the air-in-line detector ensures patient safety by detecting air presence in the administration set tubing. The ON/OFF key conserves battery life, while an event memory that records up to 1,000 events can be accessed via the CADD-DIPLOMAT PC Communications System. The pump's reliability means fewer therapy interruptions, and its programming flexibility, with infusion rates up to 125 mL/hr and dose volumes up to 1,000 ml, accommodates high-rate, high-volume therapies. Enhanced delayed start times and dosing parameters provide the adaptability needed for complex therapy regimens, ensuring that the internal clock keeps dosing on

schedule.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-CADD-LP

Price: \$2,374.00



BAXTER SIGMA SPECTRUM INFUSION PUMP WITH 6.02 SOFTWARE

The **Baxter Sigma Spectrum Infusion Pump with 6.02 Software** is a compact, user-friendly, patient-centered infusion device designed to enhance safety and efficiency in clinical settings. Equipped with Baxter's Dose Error Reduction Software (DERS), it significantly helps prevent pump-related medication errors through features like titration error prevention, flow alerts, and automatic medication dose updates. Its versatile infusion profiles support both automated ramped delivery and manual multi-step infusions, making it suitable for a wide range of therapies. The pump's compatibility with standard IV sets streamlines setup and reduces free-flow risks, while its rechargeable lithium-ion battery ensures reliable, portable operation without the need for disposable batteries. With a bright, large-type color display, pump-guided programming, and onboard help screens, the Sigma Spectrum offers quick, intuitive operation. Optional wireless communication and Continuous Quality Improvement (CQI) reporting further support therapy monitoring and compliance. Measuring just 5.8" H x 4.2" W x 2.5" D and weighing only 1.6 lbs, this lightweight infusion pump is ideal for busy healthcare environments requiring precise, safe, and efficient medication delivery.

Battery Options:

- Standard non-wireless battery
- w/ Single Antenna BG Battery
- w/ Dual Antenna ABGN Battery

Warranty: 1 Year

[Read More](#)

SKU: AM-BAXTERSS

Price: \$995.00 - \$1,305.00



BAXTER PCA II

The **Baxter PCA II** is a sophisticated syringe infusion pump designed for safe, precise, and flexible patient-controlled analgesia (PCA) in hospital settings. Capable of delivering both continuous and intermittent infusions, it features internal safety checks to reduce tampering and errors, along with a user-friendly, large back-lit display for easy operation even in low-light conditions. The device offers multiple programmable modes—including PCA only, basal plus PCA, and continuous basal—programmable by milliliters or milligrams, and accommodates both pre-filled and standard syringes for convenience and cost savings. Its built-in history logs track successful injections, attempts, alarms, and operational details, with optional hard-copy reports when connected to a printer. Compact and lightweight at 13 inches high, 6.3 inches wide, and weighing just 4.2 pounds with batteries, the Baxter PCA II provides a reliable and efficient solution for pain management with a flow accuracy of $\pm 3\%$ over a range of 0.1 to 150 mL/hr.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-BPCAII

Price: \$821.14



B. BRAUN PERFUSOR SPACE SYRINGE PUMP

The **B. Braun Perfusor Space Syringe Pump** is a highly advanced, compact medical device designed to deliver precise doses of medication to patients in various healthcare settings. Renowned for its accuracy and reliability, this syringe pump is an integral part of modern infusion therapy protocols. It is engineered to work within the B. Braun Space system, a comprehensive solution that allows for seamless integration with other infusion devices, thereby enabling a flexible, modular approach to patient care. The Perfusor Space Syringe Pump is equipped with sophisticated safety features, including automatic syringe size recognition and drug library compatibility, which significantly minimize the risk of medication errors. Its user-friendly interface ensures ease of use for healthcare professionals, facilitating the efficient administration of critical medications. With its compact size, it offers exceptional versatility, making it suitable for use in intensive care units, operating rooms, and even in patient transport scenarios.

Warranty: 1 Year

[Read More](#)

SKU: AM-BRAUN-PERFUSOR-S

Price: \$1,845.00



BAXTER 6301 INFUSION PUMP

The **Baxter 6301 Infusion Pump** is a versatile, advanced dual-channel volumetric infusion system designed to meet a wide spectrum of clinical needs with high precision and safety. Capable of delivering fluids over a broad range of rates—from as low as 1 mL/hr to a maximum of 9,999 mL/hr in single and multiple channels—the pump supports continuous, secondary (piggyback), and patient-specific protocols with automatic flow rate calculation once volume and time are set. Key features include remote operation and monitoring, an 80% recharge time of 8 hours, and enhanced safety with flow check occlusion alarms and front panel lockouts to prevent tampering. The device automatically restarts after occlusion clears, providing automatic secondary medication piggybacking, and retains infusion data for 5 hours without power. Compact, weighing 8.1 kg (17.9 lbs) with dimensions of 33 cm W x 13 cm D x 29 cm H, and powered by a 12V battery lasting up to 6 hours at high flow rates, the Baxter 6301 offers reliable, precise, and efficient infusion management for diverse patient care environments.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-BAX-6301

Price: \$1,150.00



ALARIS MEDSYSTEM 2865

The **Alaris MedSystem 2865** is a versatile, compact multi-channel infusion pump designed for a range of clinical settings, including ambulatory care, emergency medical services, and air transport. It features three independent fluid delivery channels with customizable settings, including a drug library and a user-adjustable dose calculator, enabling safe and efficient infusion management. Its lightweight design (weighing approximately 5.1 lbs) and small dimensions (8" x 6" x 2") facilitate easy transport and bedside organization. The system supports infusion of IV, intra-arterial, epidural, and subcutaneous medications from various containers, including syringes with its syringe adapter set. Key safety features include set-based, needle-free, free-flow protection. Programming is streamlined with single-page dose rate calculation, user-tunable drug library, and titration options. It operates on 6-watt AC power, with rechargeable NiCad batteries and a non-rechargeable backup, providing a wide delivery range from 0.1 to 999 mL/hr and a KVO rate from 0.1 to 20 mL/hr. The pump comes complete with a charger and pole clamp, making it suitable for flexible, safe infusion delivery in various clinical environments.


[Download Brochure](#)

Warranty: 3 Months

[Read More](#)

SKU: AM-ALARIS-2865

Price: \$3,403.91



B. BRAUN OUTLOOK 100

The **B. Braun Outlook 100** is a compact, user-friendly infusion system designed to enhance medication safety at the point of care. Its large color display and extensive customizable drug library enable clinicians to configure the pump according to specific clinical protocols, helping to reduce medication errors. Featuring a built-in "brain" for each channel, it effectively prevents channel confusion, while the active dose limit display and KeyGuard™ keypad with arrow controls simplify operation and minimize manual input errors. Additional benefits include an air trapping cassette for nuisance alarm prevention, improved battery life, and flexible configurations for loading doses and on-demand boluses for all DoseGuard drugs. Its streamlined single-screen programming limits reliance on memory, making it straightforward to use. The Outlook 100 delivers reliable, accurate infusions with a flow accuracy of $\pm 5\%$, supporting rates from 0.1 to 999.9 mL/h and a volume capacity up to 999.9 mL. Compact, stackable, and weighing just 5 kg, it is ideal for busy clinical environments requiring efficient, safe infusion management.

[Download Brochure](#)

Warranty: 1 Year

[Read More](#)

SKU: AM-BBO100IP-R

Price: \$845.00



ALARIS 8000 PC POINT-OF-CARE UNIT

The **Alaris 8000 PC Point-of-Care Unit** serves as the central hub of the Alaris System, offering a unified user interface that simplifies programming and monitoring of infusion modules, thereby reducing complexity at the point of care. Formerly known as the Medley System, this modular platform provides integrated SpO2 monitoring and ensures precise, automated infusion of a wide variety of intravascular fluids, medications, and blood products. Compact and lightweight at just 7.2 pounds, with dimensions of 6.9" W x 8.8" H x 9" D, it can operate on battery power, with runtime depending on attached modules and activity levels. Equipped with an RS-232 communication port (RJ45 connector) and advanced memory backup, the Alaris 8000 maintains system configuration parameters for up to 6 months and module-specific settings for 8 hours when powered off. Designed for versatility and ease of use, this unit is compatible with the Alaris Medley 8100 Infusion Pump, making it an essential component for efficient, safe patient infusion management in clinical settings.

[Download Brochure](#)

[Compatible with the Alaris Medley 8100 Infusion Pump](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-Alaris-8000

Price: \$660.00



ABBOTT TRIPLE CHANNEL PLUM A+ IV

The **Abbott Plum A+ Triple Channel IV Pump** is an advanced, wireless-enabled infusion system designed to deliver precise, reliable medication management in high-acuity healthcare environments such as ERs, ICUs, and oncology centers. Engineered to support up to six concurrent infusion lines in one compact unit, the Plum A+3 allows clinicians to administer simultaneous or tiered medication infusions with exceptional accuracy and safety. Its distributed system architecture helps safeguard against multiple channel errors, while automated piggyback delivery, loading dose automation, programmable delayed starts, and multistep delivery enhance clinical efficiency and workflow. With a rechargeable battery, hospital-grade 10-foot AC power cord, and the flexibility to deliver secondary infusions via IV container or syringe, this IV infusion pump is built for consistent, proven performance. Measuring 13.75" H x 12.2" W x 7.5" D and weighing 20 lbs with batteries, the Abbott Plum A+ Triple Channel system optimizes bedside space, asset management, and productivity—making it a trusted solution for complex medication dosing regimens.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-ABBOTT-A+T

Price: \$1,921.00



ABBOTT PLUM A+ IV INFUSION PUMP, SINGLE CHANNEL

The **Abbott Plum A+ IV Infusion Pump** is a versatile, reliable device designed to deliver precise, multi-channel infusions across various clinical settings. Featuring three independent channels with individual lines, it supports simultaneous delivery of multiple medications with advanced capabilities like fully automated secondary drug delivery, concurrent flow, and a variety of dosing and delivery options—including 15 preferred rates. Its innovative PlumSet™ technology ensures confident infusions, while air-in-line detection and automated backpriming minimize contamination risks by safely removing trapped air without disconnecting the patient line. The pump can administer secondary infusions via IV container or syringe and offers flexible programming features such as piggyback, multistep delivery, delayed starts, and programmable standby settings. Compact and sturdy, measuring approximately 8" H x 8" W x 6" D and weighing around 9.5 lbs, the Abbott Plum A+ provides steady, proven performance and seamlessly integrates into the entire spectrum of patient care, making it an essential tool for enhancing infusion safety and

efficiency.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-Abbott-A+

Price: \$865.00

ALARIS MEDLEY 8100 INFUSION PUMP

The **Alaris Medley 8100 Infusion Pump Module** is a high-capacity, large-volume infusion device designed to seamlessly integrate into the modular Alaris System, operating in conjunction with the Alaris PC Point of Care Unit (Models 8000 or 8015). Compact and lightweight, measuring just 3.3" W x 8.9" H x 5.5" D and weighing approximately 2.5 lbs, this positive displacement pump provides precise and reliable delivery of fluids across a flow rate range of 0.1 to 999 mL/h, with customizable VTBI up to 9999 mL and adjustable Keep Vein Open rates. Its advanced safety features include Guardrails Suite MX, ultrasonic air-in-line detection, and automatic free-flow protection, ensuring patient safety during complex infusion therapies. The pump supports multiple infusion modes—continuous, intermittent, bolus, and secondary—and features pressure sensing with adjustable occlusion limits from 50 to 525 mmHg, displayed in real time via a dynamic pressure bar graph on the PC Unit. Its compatibility with dedicated Alaris infusion sets, along with environmental resilience to operate between 41°F and 104°F, makes the Alaris 8100 an essential component for hospitals and clinics seeking reliable, safe, and versatile infusion management.

Requires the [Alaris PC Point of Care unit \(Model 8000 or 8015\)](#) to operate

6 Month Warranty

[Read More](#)

SKU: AM-Alaris-8100

Price: \$995.00





ALARIS MEDLEY 8015 POINT OF CARE MODULE

The **Alaris Medley 8015 Point of Care Module** is a versatile, modular infusion and monitoring system designed for adult, pediatric, and neonatal care, providing clinicians with a unified interface for programming and managing multiple channels simultaneously. This compact unit, weighing only 7.2 pounds and measuring 6.9" W x 8.8" H x 9" D (including pole clamp), offers wireless connectivity to enhance IV safety and facilitate seamless clinical workflows. It features powerful analytics, exceptional versatility, and the most comprehensive safety protections available through Guardrails® Suite MX software, which double-checks dose, duration, and delivery rate to prevent errors. Its advanced memory system retains configuration parameters for up to 6 months, with module-specific data stored for 8 hours when powered off. Designed for modern healthcare environments, the Alaris Medley 8015 improves patient safety, clinical efficiency, and adherence to best practices across a broad spectrum of care settings.

[Compatible with the Alaris Medley 8100 Infusion Pump.](#)

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-Alaris-8015

Price: \$725.00



BAXTER FLO-GARD 6201

The **Baxter Flo-Gard 6201 Volumetric Infusion Pump** is a dependable, high-performance IV pump designed to deliver precise fluid administration across a wide range of infusion rates. Ideal for hospitals, surgical centers, and clinical settings, this single-channel infusion pump offers enhanced safety features such as a Flow Check Occlusion Alarm with in-line resistance monitoring, automatic restart after occlusion clearance, and a slide clamp option to prevent inadvertent gravity free flow. The Flo-Gard 6201 automatically calculates flow rates once volume and time are entered, supports programmable delivery profiles with up to 10 steps for customized infusion ramping and tapering, and allows incremental flow rate adjustments from 0.1 mL/hr up to 1,999 mL/hr. With automatic piggyback capability for secondary medications, front panel lock-out protection, and compatibility with Baxter standard IV administration sets to reduce waste and costs, this volumetric infusion pump combines versatility, convenience, and reliability. Compact and portable at 11.4" high, 8" wide, and 11.6 lbs, it operates on a 12-volt sealed lead-acid battery, making the Baxter Flo-Gard 6201 a trusted solution for accurate, controlled IV therapy.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-6201

Price: \$822.00



B. BRAUN VISTA BASIC INFUSION PUMP

The **B. Braun Vista Basic Infusion Pump** is a reliable, easy-to-use volumetric IV pump designed for outpatient, home care, and clinical settings where simplicity and durability matter most. With over 100,000 units sold worldwide, this proven infusion pump delivers continuous, piggyback, and program modes while supporting cost-effective straight-line gravity sets for efficient operation. The Vista Basic features intuitive programming with helpful on-screen prompts, quick tubing installation, and adjustable occlusion settings to enhance patient safety. Built with a robust mechanical design, pump-based and optional set-based free-flow protection, and DEHP-free, needle-free IV sets, it prioritizes both patient and clinician protection. A long-lasting battery with approximately three years of life, oversized display, ergonomic handle, and integrated bag hook ensure portability and convenience without the need for an IV pole. Dependable, durable, and simple to teach and learn, the B. Braun Vista Basic is a trusted infusion pump solution for consistent, safe IV therapy delivery.

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-BBVIP-R

Price: \$1,495.00



BAXTER INFUSO.R. SYRINGE PUMP W/ PROPOFOL SMART LABEL

The **Baxter Infus O.R. Syringe Pump** is a compact, easy-to-use infusion device specifically designed to support the precise delivery of intravenous medications during anesthesia and surgical procedures. Ideal for administering narcotics, muscle relaxants, vasoactive agents, and other anesthesia-related drugs, this lightweight syringe pump weighs just 2 lbs with batteries and operates on four "C" alkaline batteries (5-7 volts) for dependable portability. Featuring magnetically coded, drug-specific **SMART LABELS**, the Infus O.R. helps minimize medication dosing errors, reduce calculation time, and simplify programming by allowing clinicians to quickly select from an available drug library and modify infusion parameters without complex reprogramming steps. Flow rates from 0-600 mL/hr, intuitive front-panel rotary switches for one-handed adjustments, and clear alarms for infusion complete, battery, and malfunction enhance safety and workflow efficiency in the operating room. Measuring 9.2 x 4.5 x 2.0 inches, the Baxter Infus O.R. Syringe Pump delivers reliable, accurate, and convenient syringe-based IV infusion therapy when precision matters most.

[Download Brochure](#)

Warranty: 1 Year

[Read More](#)

SKU: AM-INFUSOR

Price: \$4,395.00

INDEX

A

Abbott Plum A+ IV Infusion Pump, Single Channel 30
Abbott Triple Channel Plum A+ IV 30
Alaris 8000 PC Point-of-Care Unit 29
Alaris 8120 PCA Syringe Pump 4
Alaris 8300 EtCO2 Module Only 13
Alaris IVAC MedSystem III 2865B Infusion Pump 16
Alaris Medley 8015 Point of Care Module 32
Alaris Medley 8015 Point of Care Module (Software Version 9.33) 3
Alaris Medley 8100 Infusion Pump 31
Alaris Medley 8110 Syringe Pump Module 19
Alaris MedSystem 2865 28

B

B. Braun Infusomat Space 24
B. Braun Infusomat Space Infusion Pump Wifi Capable, Wifi Batteries (NEW) 10
B. Braun Infusomat Space Infusion Pump Wifi Capable, Wifi Batteries (Refurbished) 12
B. Braun Outlook 100 28
B. Braun Perfusor Space Syringe Pump 26
B. Braun Vista Basic Infusion Pump 33
Baxa Repeater Pump 9
Baxter 6301 Infusion Pump 27
Baxter AS50 Syringe Pump 20
Baxter Flo-Gard 6201 33
Baxter InfusO.R. Syringe Pump W/ Propofol Smart Label 34
Baxter PCA II 26
Baxter Sigma Spectrum Infusion Pump With 6.02 Software 25
BodyGuard 323 CME 15

C

CADD Legacy Plus Pump 24
CMExpress Infusion Pump 17
Curlin 4000 CMS 23
Curlin 6000 CMS 22

E

Eitan Medical Sapphire Epidural Infusion Pump (New) 5
Eitan Medical Sapphire Epidural Infusion Pump (Refurbished) 16
Eitan Medical Sapphire™ Multi-Therapy Infusion Pump (New) 6
Eitan Medical Sapphire™ Multi-Therapy Infusion Pump (Refurbished) 11
Eitan Medical Sapphire™ Multi-Therapy Infusion Pump with Lockbox/Cradle (Refurbished) 7

F

Foot Pedal for Baxa Repeater Pump 11

H

Hospira Lifecare PCA3 22

M

Medfusion 2010i Syringe Pump 18
Medfusion 3500 Syringe Pump with Version 5 Software 14
Medfusion 3500 V 3.0.6 Syringe Pump 2
Medfusion 3500 V4 Syringe Infusion Pump 20
Medfusion 3500 V6 Syringe Pump 9
Medfusion 4000 Syringe Pump 8

P

Plum 360 Infusystem 17

V

Volumat MC Agilia Infusion Pump 14

W

Walkmed 350VL Infusion System 21