

TABLE OF CONTENTS

Laryngoscopes 2

 Roll Stand for Infinium ClearVue XL Video Laryngoscope 3

 Infinium ClearVue™ VLXL Video Laryngoscope 3

 ClearVue VL3R™ Reusable Video Laryngoscope 4



LARYNGOSCOPES

Video laryngoscopes are advanced medical devices that enhance the intubation process by providing clear visualizations of the airway. Unlike traditional laryngoscopes, these instruments feature a camera and display system, allowing healthcare providers to see the larynx and vocal cords in real-time. This technology is crucial in emergency medicine, anesthesia, and critical care, as it significantly improves the success rate of intubations and reduces the risk of complications, especially in difficult cases.

Key models of video laryngoscopes include the **Verathon GlideScope®, Infinium ClearVue VL3R™ Reusable Video Laryngoscope, and Infinium ClearVue™ VLXL**. The GlideScope® is known for its high-definition video and integrated monitor, making it ideal for both adult and pediatric use. The ClearVue VL3R™ offers a reusable design that is compact and user-friendly, while the Infinium ClearVue™ VLXL features advanced imaging technology and ergonomic design for enhanced performance in various clinical environments. Auxo Medical offers both new and refurbished video laryngoscopes at very competitive low prices.

These devices are widely utilized across several industries, including healthcare, emergency medical services, EMS Transport, and military applications. They are essential for improving airway management, providing training tools for medical education, and ensuring patient safety in high-pressure situations. As technology advances, video laryngoscopes will continue to play a vital role in modern medicine, facilitating better outcomes for patients.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



ROLL STAND FOR INFINIUM CLEARVUE XL VIDEO LARYNGOSCOPE

The **Roll Stand for the Infinium ClearVue XL Video Laryngoscope** is a versatile and essential accessory designed to enhance the functionality and mobility of the medical device. This roll stand features a convenient basket, providing ample storage space for essential tools and accessories, ensuring they are always within reach during critical procedures. Additionally, it is equipped with up to four mounting stations for batons, allowing for organized and secure placement of multiple laryngoscopes or related equipment. The stand is designed with lockable casters, ensuring stability and safety during use, while also allowing for easy transport and repositioning as needed in a busy clinical environment. This thoughtful design ensures that healthcare professionals can maintain optimal workflow efficiency and patient care quality.

[Read More](#)

SKU: 267.300.0011

Price: Request Quote



INFINIUM CLEARVUE™ VLXL VIDEO LARYNGOSCOPE

The **Infinium® ClearVue™ VLXL video laryngoscope** adopts state-of-the-art image technology and possesses the advantage of a 2.0-megapixel full view camera with 8-inch high-resolution monitor. The ClearVue™ features a rechargeable li-ion battery for mobile use. The HD optics offer a unique anti-fog and high-resolution image. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable the clinician's first-attempt success while minimizing any chance of a complication during the intubation process.

[Read More](#)

SKU: 000.300.0009

Price: Request Quote



CLEARVUE VL3R™ REUSABLE VIDEO LARYNGOSCOPE

The **Infinium ClearVue™ video laryngoscope** adopts state-of-the-art image technology and possesses the advantage of a 2.0-megapixel full view camera with 4-inch high-resolution monitor. The ClearVue™ features a rechargeable Li-ion battery for extended use. The HD screen offers a unique anti-fog and anti-shatter capability. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable the clinician's first-attempt success while minimizing any chance of a complication during the intubation process.

[Read More](#)

SKU: 000.300.0006

Price: Request Quote



INDEX

C

ClearVue VL3R™ Reusable Video Laryngoscope 3

I

Infinium ClearVue™ VLXL Video Laryngoscope 3

R

Roll Stand for Infinium ClearVue XL Video Laryngoscope 2