


TABLE OF CONTENTS

Liposuction	2
HK Lipo Valet Cart	2
HK Vein Valet Cart	3
HK Surgical Fat Harvesting Station	4
HK Aspiration Tubing (Box of 10)	5
HK Surgical Single Spike Infiltration Tubing Box of 20	6
HK Klein Infiltration Pump KIP-II+ Plus	6
HK Surgical - Klein TDD Pump for Infiltration (KTDD)	6
HK Klein Infiltration Pump II (Refurbished)	7
MICROAIRE PAL-650 SYSTEM WITH MICROAIRE 1020 CONSOLE	7
Byron Lysonix 3000 Liposuction Unit	8
HK Surgical Aspirator Pump III	8
Aspirator III Liposuction Unit	9



LIPOSUCTION

Auxo Medical offers cosmetic and plastic surgery devices at competitive rates.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



HK LIPO VALET CART

The **HK Lipo Valet Cart** is an ideal, mobile solution tailored for practices performing tumescent liposuction, fat harvesting, and other cosmetic procedures. Designed to securely hold a Klein Pump (models KIP-II, KTP, KTIP, or KTDD), HK Aspiration Pump, foot pedals, and surgical instruments, it features a hospital-grade power strip with four outlets to keep equipment and cords organized. The cart measures 48" x 22" x 22" (with a shelf of 22" x 13.5"), weighs 48 pounds, and includes an IV pole, Mayo tray, one shelf, and holsters for foot pedals, with additional shelves available for customization. Equipped with hospital-grade locking wheels and an innovative back hole for HK Aspirator pole, the Lipo Valet provides a streamlined, organized platform for efficient surgical workflows. Shipping is included.

[Read More](#)

SKU: AM-LV



HK VEIN VALET CART

The **HK Vein Valet Cart** is a specialized, mobile workstation designed for practices, surgery centers, and hospital ORs performing varicose vein procedures using tumescent local anesthesia. Crafted to enhance efficiency and organization, this redesigned cart features a welded storage basket for procedure supplies, a shelf for the Klein Pump, a shelf for the generator and vein light, and a hospital-grade power strip with four outlets to power devices like the Klein infiltration pump and other equipment. Its dimensions are 18" x 23" x 42" (with a shelf size of 22" x 13.5"), equipped with hospital-grade locking wheels for stability and mobility, an IV pole, holsters for foot pedals, a sterile mayo tray, and ample storage for cannulas and tubing. The Vein Valet is an essential accessory for streamlining vein procedures, keeping devices and supplies organized and within reach for maximum efficiency. Shipping is included in the price.

[Read More](#)

SKU: AM-VV



HK SURGICAL FAT HARVESTING STATION

The **HK Surgical Fat Harvesting Station** is a comprehensive, mobile setup designed for cosmetic surgery centers and in-office procedures requiring efficient fat harvesting. The station includes a Klein Infiltration TDD Pump with foot pedal, an AP-III Aspiration Pump with dual foot pedals, the versatile Lipo-Valet Cart, tabletop IV poles, and a Filtron Canister Bracket (canisters sold separately). This integrated system ensures all essential fat harvesting equipment is easily accessible, highly portable, and conveniently organized for use within the operating room or for storage when not in use. Compatible with the Klein KIP-II infiltration pump, HK AP-III aspirator, and tabletop IV poles, this station streamlines workflow, enhances mobility, and optimizes space management during cosmetic and surgical procedures.

[Read More](#)

SKU: AM-HK-FHS




HK ASPIRATION TUBING (BOX OF 10)

The **HK Aspiration Tubing (Box of 10)** is a sterile, high-quality, super-soft aspiration tubing designed for use with all liposuction aspirators, including Wells Johnson Units and HK Surgical Units. Each 10-foot length features a 3/8" inner diameter, 5/8" outer diameter, and 1/8" wall thickness, providing durability and resistance to collapse during high vacuum procedures. The tubing's flexible design ensures reliable and consistent aspiration performance, making it an essential component for safe and effective liposuction and surgical aspiration applications.

[Read More](#)

SKU: AM-AT6-10



HK SURGICAL SINGLE SPIKE INFILTRATION TUBING BOX OF 20

The **HK Surgical Single Spike Infiltration Tubing (Box of 20)** consists of 15.5-foot-long, sterile silicone pressure tubing designed specifically for use with HK Klein infiltration pumps. Each sterile pack includes 20 pieces of standard tubing with a single spike (luer) connector, providing a reliable and convenient solution for delivering tumescent anesthesia during procedures such as liposuction, cosmetic surgery, and vascular treatments. The durable silicone pressure tube ensures consistent flow, compatibility with HK Klein pumps, and ease of use in various surgical settings.

[Read More](#)

SKU: AM-ITS-20



HK KLEIN INFILTRATION PUMP KIP-II+ PLUS

The **HK Klein Infiltration Pump KIP-II+ Plus** is a next-generation, efficient device designed to streamline the delivery of tumescent anesthesia, significantly reducing infiltration and procedure times while enhancing patient and clinician comfort. The pump offers a maximum flow rate of up to 1000 mL/min, with a motor speed adjustable from 0 to 600 RPM, and features two operating modes—Momentary and Toggle—controlled via an easy dial and foot pedal for hands-free operation. Its precise, single-direction flow ensures consistent administration, supported by an audible signal with volume control to monitor operation. Compatible with Klein's infiltration tubing and other accessories, the device is CE marked, adaptable to 115 or 230 volts, and includes all essential components such as power cords, tubing, a foot switch, and instructions. Ideal for procedures like liposuction, cosmetic and plastic surgery, vascular treatments, and skin grafting, the KIP-II+ Plus optimizes infiltration procedures with reliability and ease.

[Read More](#)

SKU: AM-KIP-II+



HK SURGICAL - KLEIN TDD PUMP FOR INFILTRATION (KTDD)

The **HK Klein TDD Infiltration Pump** delivers high-precision tumescent delivery with adjustable pump rates in 1 RPM increments, ranging from 1 to 600 RPM and up to a 1,000 mL/hour flow rate. It offers bi-directional pumping, a membrane touch keypad, and convenient control via foot pedal or finger-touch—plus adjustable sound volume for a quieter, more comfortable workflow. Features include save speed settings, a long life expectancy of 77,000 hours, and a 3-year limited warranty. The package comes complete with a foot pedal, Klein single spike infiltration tubing, power cord, IFU, and an optional additional foot pedal. HK Surgical is ISO 13485:2016 certified, underscoring its commitment to quality and safety.

[Read More](#)

SKU: AM-KTDD



HK KLEIN INFILTRATION PUMP II (REFURBISHED)

The **HK Klein Infiltration Pump II (Refurbished)** is a vital tool designed to efficiently deliver tumescent anesthesia, enhancing patient comfort and reducing infiltration time. With a weight of 26 lbs and compact dimensions of 20" x 16" x 13", this pump is easy to maneuver and integrate into surgical settings. It features a safe, non-conductive roller head for added safety during procedures and operates on 24V DC power at 1 AMP. Capable of delivering up to 1000 ml/min with a maximum head pressure of 25 PSI continuously (and up to 40 PSI intermittently), it offers dual pneumatic control for precise operation. Stackable for space-saving storage and compatible with Klein Aspirator, this pump is a reliable choice from HK Surgical for efficient infiltration.

[Read More](#)

SKU: AM-KIPII-R



MICROAIRE PAL-650 SYSTEM WITH MICROAIRE 1020 CONSOLE

The **Microaire PAL-650 System**, when paired with the Microaire 1020 Console, is a sophisticated setup designed for precision and efficiency in surgical procedures. The system includes the Microaire REF 1020 Instrument Electric Console, which serves as the central control unit, ensuring reliable power and intuitive operation. Complementing the console is the Microaire REF PAL-650 PAL Electric Handpiece, known for its ergonomic design and consistent performance, making it an essential tool for surgeons. Additionally, the package includes the Microaire REF 1006-PALE PAL Electric Instrument Cable, which provides the necessary connectivity and flexibility to ensure seamless integration and functionality of the entire system. Together, these components create a robust and efficient surgical solution, enhancing the quality and outcomes of medical procedures.

[Read More](#)

SKU: AM-PAL650-1020

PSI-TEC III
Aspirator



EZ PUMP
Infiltrator



Programmable system with flow rate of 100-1000 ml/min. Infiltrator or 1000-15000 ml/min. a great method to manage infiltration records.

LYSONIX 3000
ULTRASONIC SYSTEM



Next Generation of precision ultrasound with new dimensions of control. Different energy modes, handpieces, & Probes

BYRON LYSONIX 3000 LIPOSUCTION UNIT

Features:

- PSI-TEC III Aspirator
- EZ Pump Infiltrator
- Lysonix 3000

[Read More](#)

SKU: AM-B-L3000-EZ-PTIII





HK SURGICAL ASPIRATOR PUMP III

The **HK Surgical Aspirator Pump III** is a reliable, portable, and efficient suction device designed for medical applications requiring high torque and sustained vacuum performance. Equipped with a piston/reciprocating pump, it maintains a higher vacuum level with adjustable settings up to 27 inches of mercury (Hg) at sea level, offering maximum flow rates of 3 cubic feet per minute without intake restrictions. Its compact design, weighing approximately 50 pounds and measuring 19" x 13½" x 16", makes it ideal for clinical environments where mobility and reliability are essential. The pump features a minimal control panel, integrated pressure gauge protected from damage, and can be placed atop infiltration pumps for added versatility. Operating quietly with low noise levels, it is supported by a lifetime team support after the sale and offers service and loaner pump programs. Included accessories encompass canisters, foot pedals, brackets, tubing, and filters, making it a practical solution for efficient and effective suction needs in surgical and medical procedures.

[Read More](#)

SKU: AM-APIII



ASPIRATOR III LIPOSUCTION UNIT

The Wells Johnson Aspirator III Liposuction Machine includes new interconnecting tubing, new filter, overflow trap, gauge, foot pedal, and power cable. It is certified to proper operating specifications and in excellent condition. Superior quality inside and out, and state-of-the-art technology in exterior component design with a new space-age polymer casing that is scratch resistant and defies denting. Even with three high-powered dual-cylinder piston pumps, its operation is quiet to the ear. Standard all brass fittings resist corrosion.

[Read More](#)

SKU: AM-WJ-ASPIII

INDEX

A

Aspirator III Liposuction Unit 8

B

Byron Lysonix 3000 Liposuction Unit 7

H

HK Aspiration Tubing (Box of 10) 4

HK Klein Infiltration Pump II (Refurbished) 6

HK Klein Infiltration Pump KIP-II+ Plus 5

HK Lipo Valet Cart 2

HK Surgical - Klein TDD Pump for Infiltration (KTDD) 6

HK Surgical Aspirator Pump III 8

HK Surgical Fat Harvesting Station 4

HK Surgical Single Spike Infiltration Tubing Box of 20 5

HK Vein Valet Cart 3

M

MICROAIRE PAL-650 SYSTEM WITH MICROAIRE 1020 CONSOLE 7