

TABLE OF CONTENTS

Microscopes 2

 Zeiss OPMI Visu 150 Surgical Microscope on S7 Stand 2

 Leica M844 Surgical Microscope on F40 Stand 3

 Zeiss Lumera T S88 Floorstand With Assistant Package 4



MICROSCOPES

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



ZEISS OPMI VISU 150 SURGICAL MICROSCOPE ON S7 STAND

The **Zeiss OPMI Visu 150 Surgical Microscope**, when paired **with the S7 stand**, is a popular choice for cataract surgery due to its excellent Red Reflex and multiple angles of illumination. This microscope offers great optical clarity and features Zeiss apochromatic optics. It includes a Zeiss 0-180 degree inclinable binocular, four 10x widefield eyepieces, and a Zeiss F175 T* objective lens (175mm). The system also comes with a Sony ExwaveHAD camera head and control unit, a Zeiss footswitch, and a new fiber-optic cable. Specifications include a light system with red reflex, field lighting up to 65mm, and a halogen 12v 100w light source with manual focus. The microscope has a powered apochromatic optical zoom ratio of 1:4, XY movement of 100mm x 100mm with automatic centering, and a 50mm stroke motor for approach. It also includes a G 475 filter for eye protection, as well as UV and IR filters.

[Read More](#)

SKU: AM-Z-OPMI-VISU-150-S7

Price: Request Quote



LEICA M844 SURGICAL MICROSCOPE ON F40 STAND

The **Leica M844 F40** is a **premium ophthalmic surgical microscope** designed for high-resolution clarity and optimal contrast with patient-safe low-light levels through direct illumination. It features a quad zoom with four independent beam paths and a double wing for simultaneous observation by three individuals. The microscope includes main and secondary direct illumination using 12V 50W halogen lamps, Ottoflex Red Reflex Enhancement, and APO chromatic corrected optics. It also offers a 0° assistant rotating observer for temporal approach and a built-in rear-facing video port. Designed for various surgical applications such as cataract, vitreo-retinal, and glaucoma procedures, the Leica M844 F40 integrates seamlessly with lasers and includes features like the Leica ToricEyePiece for IOL placement and the Leica DI C800 imaging module for high-resolution digital imaging.

[Read More](#)

SKU: AM-L-M844-F40

Price: Request Quote



ZEISS LUMERA T S88 FLOORSTAND WITH ASSISTANT PACKAGE

The **Zeiss Lumera T S88** is an outstanding microscope designed for a full suite of ophthalmic procedures, excellent for the University Medical Center and specialists in vitreo retinal procedures. If you are looking to achieve the best in red reflex technology, look no further than the Lumera T S88 ophthalmic surgical microscope. With the introduction of Lumera Stereo Coaxial Illumination you see more contrast with this brilliant, stable red reflex – even in decentered or strongly pigmented eyes. The Lumera T S88 floor stand provides superior surgeon control from the ergonomic hand grips to release the magnetic clutch brakes moving the microscope across the surgical field with minimal effort. Program the Lumera T S88 with up to 9 different users with preset motorized speeds, illumination intensity and foot pedal controls.

[Read More](#)

SKU: AM-ZEISS-LUMERA-T-S88

Price: Request Quote

INDEX

L

Leica M844 Surgical Microscope on F40 Stand 2

Z

Zeiss Lumera T S88 Floorstand With Assistant Package 3

Zeiss OPMI Visu 150 Surgical Microscope on S7 Stand 2