

TABLE OF CONTENTS

*Sale	3
Enthermics EC230L Fluid Warmer	3
Steris Dual Compartment Blanket Warmer QDJ02	4
Stretcher, Stryker Advantage, IV pole	4
Cardiac Science Powerheart G5, Semi-Auto, Dual Language	5
Zoll R Series defibrillator, with SP02, NIBP, EtCO2	6
Steris InnoWave™ Unity IW523046NPR Ultrasonic Cleaner	6
GE Datex Ohmeda CardioCap 5 Monitor	6
Philips InvelliVue MP5 Fixed Wall Mount Kit	7
Biodex 058-878 30in Wide Top Imaging Table	7
Morgan MEDesign FLTLT-FW (Float-Top, Free-Wheel Configuration)	8
Mindray Accutor 7, NIBP, Masimo, SPO2, Truetemp, Pulse Rate	8
Steris Harmony LA 585 LED Dual Head OR Lights	9
Valleylab Force FX	9
Verathon GlideScope Monitor with Baton	10
Physio-Control Lifepak 20E Defibrillator	11
Stryker 1068 Eye Surgery Stretcher	11
Welch Allyn SSM Vital Signs Monitor	12
Tuttnuaer 2540EK Kwiklave Electronic Autoclave (Refurbished)	12
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide and Carbon Fiber Extension (Stainless Steel Base)	13
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide (Stainless Steel Base)	14
Midmark Ritter 255 Floor Stand Light	14
Midmark Ritter 355 Floor Stand Light	15
Verathon Glidescope Roll Stand	15
Infinium ClearVue™ VLXL Video Laryngoscope	16
GE Dash 4000 Patient Monitor W/ CO2	17
UMF Power Exam Hi-Lo Table 5060	17
Stryker 1501 Renaissance Series stretcher with IV pole	18
ClearVue VL3R™ Reusable Video Laryngoscope	19
Edan IM70 Patient Monitor (Refurbished)	19
Stryker M Series Model 1015 Stretcher, with Big Wheel	19
Mortara ELI 250 ECG Machine	20
Steris Amsco Model 3085 Surgical Table	20

Mindray Passport V with CO2	21
Midmark 623 Exam Table	21
Steris Caviwave Single Chamber Ultrasonic Cleaner 15 Gal	22
Midmark 641 Procedure Chair	22
Byron Lysonix 3000 Liposuction Unit	23
Tuttnauer EZ10k 220V Autoclave	23
Abbott i-STAT 1 Gen 2 Handheld Blood Analyzer	24
Mindray DPM6 PACU Monitor	24
Steris DJ04 Dual Compartment Blanket Warmer	25
Steris Caviwave Pro Single Chamber 17 GAL Ultrasonic Cleaner	25
Skytron 6600 Surgical Table	26
GE DINAMAP Carescape V100 Vital Signs Monitor	27
Medfusion 3500 Syringe Pump with Version 5 Software	27
Welch Allyn 6000 Series Vital Signs Monitor w/ CO2	28
Aaron Bovie 950 High Frequency Generator	29
Steris Caviwave Rinsing and Drying Ultrasonic Cleaner	29
Philips Intellivue MP50 Patient Monitor	30
GE Dash 4000 Patient Monitor	30
Mindray Accutorr V Vital Signs Monitor	31
Gendron 1190-231 Bariatric Stretcher	31
Steris AMSCO QDJ05 Single Compartment Warming Cabinet	32
Hausted Video Imaging Chair	32
Hill-Rom Transtar 8010 Eye Stretcher	33
Stryker 735 Transport Stretcher with IV Pole	34
Criticare nGenuity CO2	34
Baxter Sigma Spectrum Infusion Pump With 6.02 Software	34
Zoll R Series Defibrillator	35
Midmark Ritter M11 Ultraclave (Autoclave), M11-022	36
Ritter Midmark M9, M9-001	36
Alaris Medley 8100 Infusion Pump	37
Alaris Medley 8015 Point of Care Module	38
Midmark Ritter 75L Exam Chair	38
ConMed Sabre 2400	39
Edan SE-1200 and SE-1200 Express ECG	40
Zoll AED Plus	40
Zimmer ATS 3000	40
Zimmer ATS 2000	41



auxoMEDICAL[®]

SALE ITEMS

**ENJOY SPECIAL
DISCOUNTS ON THE
FOLLOWING ITEMS***

*Restrictions May Apply and while supplies last.

***SALE**

**WHEN AVAILABLE, THIS CATEGORY SHOWS
ALL THE SELECTED EQUIPMENT IN A SALE OR
PROMOTION.**

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the
equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



ENTHERMICS EC230L FLUID WARMER

The **Fully Refurbished Entermics EC230L Fluid Warmer** is a reliable and efficient countertop device designed to keep IV fluids warmed, supporting patient recovery and reducing surgical site infection risks. Capable of holding up to 11 one-liter bottles or 16 one-liter bags, it provides consistent warming in a compact footprint measuring 20.7 inches high, 16.2 inches wide, and 28 inches deep. Constructed with heavy-duty, durable cabinets built for long-term use, the warmer features locked handles for tamper-proof security and an insulated glass door that allows staff to quickly assess fluid levels at a glance. Its stackable design enables multiple units to be easily combined, making it an ideal solution for busy medical environments seeking dependable, straightforward fluid warming.

[Read More](#)

SKU: AM-EC230L



STERIS DUAL COMPARTMENT BLANKET WARMER QDJ02

The **Steris Model QDJ02 Dual Compartment Blanket Warmer** is a completely remanufactured and biomedically inspected unit designed to provide reliable, precise warming for blankets and IVs. It features two independent compartments with digital controls, allowing users to set temperatures within the range of 90 to 160°F for optimal patient comfort and safety. The upper compartment measures 23.5" W x 16" D x 13" H, and the lower compartment is 23.5" W x 16" D x 38" H, all housed within a modern stainless steel chassis. Enhanced with a new circulation fan, digital displays, and overheat indicators, this contemporary warmer offers durable, efficient performance suitable for medical facilities seeking a high-quality, dependable warming solution.

[Read More](#)

SKU: AM-QDJ02



STRETCHER, STRYKER ADVANTAGE, IV POLE

The **Stryker Advantage Stretcher with IV Pole** is a durable, versatile patient transport and treatment platform designed for demanding clinical environments. With an overall length of 83 inches and adjustable height ranging from 22 to 35.5 inches, it accommodates a wide range of patient needs. The stretcher features a spacious 26" x 76" patient surface, siderails that can be raised or lowered for safety and ease of transfer, and sturdy 8-inch casters for smooth mobility. It supports various positioning adjustments, including a knee gatch from 0° to 30°, a backrest from 0° to 90°, and Trendelenburg/reverse tilt of $\pm 18^\circ$, with a maximum patient weight capacity of 500 lbs. Equipped with a fixed IV pole and compatible with models 1001, 1009, 1501, and 1509, this lightweight, reliable stretcher is ideal for patient transport, procedures, and emergency care.

[Read More](#)

SKU: AM-SADV-IV



CARDIAC SCIENCE POWERHEART G5, SEMI-AUTO, DUAL LANGUAGE

The **Cardiac Science Powerheart G5 Semi-Auto AED with Dual Language capability** is a state-of-the-art emergency device designed to maximize survival chances during sudden cardiac arrest. Easily switch between English and Spanish at the push of a button, this AED features RESCUE READY® technology, conducting daily self-checks of all critical components and regular high-voltage system charges to ensure optimal readiness. The device provides clear visual and voice prompts, along with RescueCoach® CPR prompting and a metronome to guide rescuers through effective chest compressions. Its non-polarized electrode pads simplify placement, reducing steps during emergencies. Equipped with an IntelliSense battery that monitors performance and storage conditions, the G5 delivers a shock in less than 10 seconds after CPR, making it a reliable, user-friendly solution for both professional and lay responders seeking to provide immediate, life-saving assistance.

[Read More](#)

SKU: AM-CSG5S-80C-S-R





ZOLL R SERIES DEFIBRILLATOR, WITH SPO2, NIBP, ETCO2

The **Zoll R Series Defibrillator** seamlessly integrates advanced life support features with comprehensive monitoring capabilities to enhance patient care during critical emergencies. Equipped with Real-CPR-Help for ALS, See Thru CPR filtering, and the CPR Compression Index, it guides responders through high-quality compressions with real-time feedback. The innovative OneStep system simplifies electrode and cable application, while Smart Prompts and the Code Readiness testing system support efficient and effective response. Providing familiarization through the Defib Mentor aid, this device features rectilinear biphasic technology, a 3-lead ECG cable, and integrated SPO2, NIBP, and EtCO2 monitoring—all with all accessories included, ensuring the unit is patient-ready. Biocertified for safety and reliability, the Zoll R Series is a comprehensive, user-friendly defibrillator designed to optimize outcomes in emergency settings.

[Read More](#)

SKU: AM-ZR3BIPSBCA



STERIS INNOWAVE™ UNITY IW523046NPR ULTRASONIC CLEANER

The **Steris InnoWave™ Unity IW523046NPR Ultrasonic Cleaner** is a highly efficient, three-phase ultrasonic irrigator designed to meet the demanding needs of modern operating rooms. Standing at 42.5 inches tall, with a width of 34.8 inches and a depth of 25.8 inches, this system can accommodate up to 16 lumened devices per cycle, making it ideal for high-volume surgical environments. Capable of processing approximately 25 pounds of equipment per cycle, it features automated, precise dosing of cleaning chemistry to ensure consistent results and reduce errors. Its customizable cycles allow users to tailor cleaning parameters to meet specific device instructions for use (IFUs), ensuring compliance and optimal cleaning efficacy. Additionally, the Unity IW523046NPR offers digital cycle data download via USB, storing up to 10,000 cycles for easy audit readiness and regulatory compliance.

[Read More](#)

SKU: AM-IW523046NPR



GE DATEX OHMEDA CARDIOCAP 5 MONITOR

The **GE Datex Ohmeda CardioCap 5 Monitor** offers a comprehensive hemodynamic monitoring solution with an intuitive user interface designed for ease of use. It features 5-lead ECG with ST segment analysis, impedance respiration, pulse oximetry (SpO₂), and non-invasive blood pressure (NIBP) measurement. Built for reliable, real-time patient assessment, the CardioCap 5 provides essential vital sign monitoring in a user-friendly platform suitable for various clinical environments.

[Read More](#)

SKU: AM-CC5

Price: ~~\$1,995.00~~ \$1,795.50



PHILIPS INVELLIVUE MP5 FIXED WALL MOUNT KIT

The **Philips IntelliVue MP5 Fixed Wall Mount Kit** features a sturdy 6-inch (15.2 cm) fixed arm designed for secure wall mounting of the patient monitor. It allows for tilt and rotation adjustments to optimize viewing angles and patient access. The kit includes cable hooks positioned beneath the arm for organized hanging of patient leads, enhancing clutter-free workspace and safety. Please note that the 19-inch GCX channel for mounting is sold separately, providing flexibility to customize the mounting setup according to clinical needs.

[Read More](#)

SKU: AM-PH-0061-01

Price: ~~\$195.00~~ \$175.50



BIODESX 058-878 30IN WIDE TOP IMAGING TABLE

The **Biodesx 058-878 30" W Top Imaging Table** is a C-Arm imaging table that is a high-stability, three-axis motorized platform designed for precise, quiet positioning and optimal image access. It features a low-attenuation carbon fiber 30" wide top, a 93" overall length, and a 500 lb weight capacity to accommodate a wide range of patients. The rectangular 90" L x 30" W top sits on high-speed actuators with 30"-40" height adjust, 0-20° Trendelenburg and reverse Trendelenburg, and 0-20° lateral roll for versatile patient positioning. The table is equipped with hand control for all movements, locking 5" swivel casters for mobility, and supports both 115 V and 230 V configurations. A compact base and optional battery operation enhance flexibility in the OR, facilitating rapid setup and stable, vibration-free imaging during procedures.

[Read More](#)

SKU: AM-BD-058-878



MORGAN MEDESIGN FLTLT-FW (FLOAT-TOP, FREE-WHEEL CONFIGURATION)

The **Morgan MEDesign FLTLT-FW Float-Top, Free-Wheel Imaging Table** is a high-precision fluoroscopy-ready C-Arm compatible platform featuring a float-top, radiolucent surface and multidirectional motion. Its free-floating longitudinal and lateral movement enables micro-adjustments during procedures without repositioning the patient, while a powered height-adjust system provides ergonomic, effortless vertical changes. With a free-wheel base and integrated locking, the table moves easily yet stays securely in place when needed. The radiolucent tabletop, coupled with Trendelenburg and reverse Trendelenburg options on select configurations, supports a wide range of interventional and diagnostic applications. Designed to handle up to approximately 500 lbs (227 kg), it measures about 80" (203 cm) in length and 24" (61 cm) in width, with a height range of 28-40" (71-102 cm). The table weighs around 1,250 lb (567 kg), reflecting its robust construction and suitability for demanding clinical environments.

[Read More](#)

SKU: AM-MORGAN-FLTLT-FW



MINDRAY ACCUTOR 7, NIBP, MASIMO, SPO2, TRUETEMP, PULSE RATE

The **MINDRAY ACCUTOR 7** is a non-invasive patient monitoring system that combines NIBP (blood pressure), SpO2, temperature (Truetemp), and pulse rate into one compact unit. It features a color HD LCD display and supports optional Masimo SET or Nellcor OxiMax SpO2, as well as Mindray TrueTemp temperature modules for flexible, accurate vital signs monitoring. Each unit comes with essential accessories including an adult cuff, NIBP hose, power line, SPO2 cable, adult finger sensor (Masimo or Nellcor compatible), oral/axillary temperature probe, probe covers, and an operator's manual, making it ready for quick setup in clinical environments.

[Read More](#)

SKU: AM-ACCUTOR 7-M-T

Price: \$1,340.00 ~~\$1,206.00~~



STERIS HARMONY LA 585 LED DUAL HEAD OR LIGHTS

The **Steris Harmony LED 585** surgical lighting and visualization system is engineered to enhance surgical precision and patient outcomes. Its high-performance optics deliver exceptional surgical site illumination with brighter, whiter light, superior color rendition, and virtually shadow-free visibility thanks to 84 overlapping beams. Beyond illumination, the system's HD-capable suspension provides stunning clarity for video signals, accommodating various HD technologies and allowing for future upgrades and flexibility. By reducing energy consumption, eliminating bulb changes, and minimizing beam heat, the Harmony LED 585 creates a cooler, more comfortable operating room environment while offering years of reliable performance.

[Read More](#)

SKU: AM-LA585-DCM



VALLEYLAB FORCE FX

The **Valleylab Force FX™** generator offers surgeons precise power delivery, even through resistant tissue, thanks to its innovative Instant Response™ technology. This technology optimizes performance at lower power settings, reducing the risk of tissue damage, neuromuscular stimulation, and the need for higher power settings. The generator automatically adjusts to tissue changes, maintaining power delivery and minimizing drag, while also reducing capacitive coupling by 30-50%. Featuring three internal microcontrollers for rapid response and increased processing speed, the Force FX™ offers versatile functionality with three cut modes, four coag modes, and three bipolar modes. Its compatibility with the Computer Motion Hermes™ Voice Command System and lightweight design (under 18 lbs) make it a valuable and adaptable tool for any operating room.

[Download Brochure](#)

[Read More](#)

SKU: AM-VFX



VERATHON GLIDESCOPE MONITOR WITH BATON

The **Verathon GlideScope Monitor** (Model: 0570-0338) **with Baton** (Model: 0570-0307) is a portable video laryngoscopy system designed to provide a consistently clear, real-time view of the airway, enhancing first-pass intubation success and patient safety. The system features a high-resolution 6.4-inch color monitor, image and video recording capabilities for training and documentation, and single-use video batons to eliminate cross-contamination. Weighing just 2.2 lbs (1 kg), the durable monitor is easy to transport and maneuver, making it suitable for a variety of settings, including routine intubations, difficult airway management, emergency situations, operating rooms, ICUs, and emergency departments. Its versatile blade options, user-friendly interface, anti-fogging feature, and battery-powered operation make it an invaluable tool for clinicians managing diverse patient populations, including pediatric cases, and for anesthesia training purposes. The monitor's dimensions are 6.85" x 8.78" x 3.15" (174 mm x 223 mm x 80 mm) and it operates within a temperature range of 50°-104° F (10°-40°C).

Includes Monitor: 0570-0338

Includes Baton: 0570-0307

Stand not included - [order AM-0570-RS](#)

Warranty: 90 Days



[Read More](#)

SKU: AM-0570-0338-0570-0307

Price: ~~\$4,975.00~~ \$4,477.50



PHYSIO-CONTROL LIFEPAK 20E DEFIBRILLATOR

The **Physio-Control LIFEPAK 20e defibrillator/monitor** is designed for in-hospital first responders, offering an intuitive interface with sophisticated features crucial for saving lives during cardiac arrest. It utilizes Physio-Control ADAPTIV™ biphasic technology, delivering energy settings up to 360 joules for difficult-to-defibrillate patients, which has proven effective in terminating shockable heart rhythms in patients requiring multiple shocks. The device adapts to various patient environments and is easy to use, even for those infrequently trained in AED operation. Advanced monitoring parameters, such as ECG, external pacing, and pulse oximetry, are clearly displayed through color-matched waveforms and values. The LIFEPAK 20e performs automated daily self-tests and features a battery gauge to ensure readiness. It includes a hospital-grade power cord, 3-lead ECG cable, Quik Combo cable, SpO2 adult reusable finger sensor, test plug, new battery, and a set of adult pads. **It is biocertified and comes with a 1-year warranty.** The device offers two display options, seven main operating modes, and can print a CODE SUMMARY™ report and transfer data records via IrDA version 1.0.

[Download Brochure](#)

[Read More](#)

SKU: AM-L20EBIPA



STRYKER 1068 EYE SURGERY STRETCHER

The **Stryker 1068 Eye Surgery Stretcher** is expertly refurbished, offering a combination of functionality and comfort with its new pad. Designed for versatility in medical settings, it features a height range of 22" to 35.5" and is equipped with 8" casters for smooth mobility. The stretcher supports a Fowler position from 0° to 90°, and includes a knee gatch adjustment from 0° to 30°, ensuring optimal positioning for various procedures. The headpiece height can be adjusted by +/-3° with a tilt of +/-45°, providing tailored support for patient needs. With a patient surface width of 26" and a weight capacity of 500 lbs, it offers ample space and strength. Additionally, the stretcher's overall dimensions are 31.5" in width and 90" in length, with side rails measuring 13" x 55", emphasizing its robust design and safety features essential for eye surgery settings.

[Read More](#)

SKU: AM-S1068



WELCH ALLYN SSM VITAL SIGNS MONITOR

The **Welch Allyn SSM 300 Series Vital Sign Monitor** is a versatile and reliable device designed to meet the diverse needs of healthcare providers. Equipped with features such as noninvasive blood pressure (NIBP) measurement, pulse rate monitoring, and mean arterial pressure (MAP) calculation, this monitor ensures precise and comprehensive vital sign assessment. It also incorporates SureTemp® Plus Thermometry and Nellcor pulse oximetry to deliver accurate temperature readings and oxygen saturation levels. Catering to adult, pediatric, and neonatal patients, the device is adaptable across various clinical settings. Its user-friendly design is highlighted by an LCD display that simplifies menu navigation, making it easy for healthcare professionals to operate. Additional features include an automatic blood pressure mode, programmable alarms for patient safety, built-in memory for data storage, and intuitive icons that enhance usability. Lightweight and portable, the SSM 300 Series offers a range of mounting options, ensuring it can be conveniently integrated into any healthcare environment.

[Read More](#)

SKU: AM-53NTO

Price: ~~\$895.00~~ \$805.50



TUTTNAUER 2540EK KWIKLAVE ELECTRONIC AUTOCLAVE (REFURBISHED)

The **Tuttnauer 2540EK Kwiklave Electronic Autoclave**, available in a refurbished model, is a versatile and efficient sterilization device ideal for various medical and laboratory settings. It offers four distinct programs—Unwrapped, Wrapped, Liquids, and a Drying cycle—facilitating comprehensive sterilization processes. Its intuitive one-touch operation allows electronics to manage the entire cycle seamlessly from start to finish, ensuring ease of use and reliability. The autoclave features a spacious 22-gallon chamber volume and a closed-door drying cycle, maintaining the sterility of instruments post-sterilization. Designed with user convenience in mind, it includes an automatic shut-off, audible and visual alerts for cycle interruptions, and a user-friendly keypad. For added utility, an optional printer is available for documenting sterilization cycles. All Tuttnauer autoclaves are A.S.M.E. stamped and U.L. approved, guaranteeing high standards of safety and performance. The autoclave has a chamber size of 10 inches in diameter and 19 inches in depth, accommodating four trays with each tray measuring 16.3 inches by 6.7 inches. Operating on 220-230V, this steam autoclave (Class S) is compact yet robust, with overall dimensions of 21.5 inches in depth, 20 inches in width, and 14.4 inches in height, and has a right-side hinged door swing.


[Read More](#)


SKU: AM-2540EK



MINDRAY HYBASE V6 SURGICAL TABLE W/ LONGITUDINAL SLIDE AND CARBON FIBER EXTENSION (STAINLESS STEEL BASE)

The **Mindray HyBase V6 Surgical Table with Longitudinal Slide and Carbon Fiber Extension (Stainless Steel Base)** is designed to meet the complex demands of modern surgical procedures. With an impressive load capacity of up to 460 kg (1014 lbs), it accommodates a wide range of patients, ensuring stability and safety during operations. The table permits all positions up to an overall load of 250 kg (551.1 lbs), offering versatility for various surgical needs. Its advanced design includes a longitudinal shift of up to 350mm (13.8"), facilitating seamless access to C-arm imaging, complemented by head-end and foot-end shifts of 6.3" and 7.5" respectively. The table's height adjusts from a low 600mm (23.6") for neurosurgery to a high 1,050 mm (41.3") for orthopedic procedures. It features 12 electric movements, including table top adjustments and tilts, all controlled via a Bluetooth hand control for ease of use in the operating room. Ensuring maximum safety, the table includes an independent manual override system, a triple-layer





decompression pad for patient comfort, and IPX5 water resistance. Its bellows-free design, coupled with four double-swivel castors, enhances maneuverability. The package includes arm boards, an anesthesia frame, memory foam pads, a stainless steel base cover, a longitudinal slide, a kidney lift, and a carbon fiber extension that optimizes imaging areas at both the head (64.9") and foot (56.7") ends, making it a comprehensive solution for healthcare providers.

[Read More](#)

SKU: AM-121-002775-00



MINDRAY HYBASE V6 SURGICAL TABLE W/ LONGITUDINAL SLIDE (STAINLESS STEEL BASE)

The **Mindray HyBase V6 Surgical Table with Longitudinal Slide (Stainless Steel Base)** is designed to meet the diverse needs of surgical teams. With an impressive load capacity of up to 460 kg (1014 lbs) and the ability to accommodate all positions up to an overall load of 250 kg (551.1 lbs), this table ensures stability and safety during complex procedures. It features a longitudinal shift of up to 350mm (13.8"), providing easy access to C-arm imaging, crucial for precision in surgical interventions. The table offers a head-end shift of 6.3" and a foot-end shift of 7.5", with a versatile height range from 600mm (23.6") for neurosurgical procedures to 1,050mm (41.3") for orthopedic operations. With 12 electric movements, including table top adjustments and tilt functions, the HyBase V6 adapts seamlessly to various surgical requirements. Its Bluetooth hand control offers unobstructed management in the operating room, while an independent manual override system enhances safety. The table is equipped with a triple-layer decompression pad for optimal patient comfort and weight distribution, and boasts IPX5 protection against water jets. Its bellows-free design and four double-swivel castors enhance mobility. Inclusions such as arm boards, an anesthesia frame, memory foam pads, a stainless steel base cover, a longitudinal slide, and a kidney lift add to its functionality, with additional options like a carbon fiber table top available for customization.



[Read More](#)

SKU: AM-121-002774-00



MIDMARK RITTER 255 FLOOR STAND LIGHT

The **Midmark Ritter 255 Floor Stand Light** is designed to deliver superior, white, shadow-free illumination, ensuring safety and coolness to the touch. The innovative yoke design allows for smooth and effortless maneuverability, making it easy to position exactly where needed. With an exceptional color rendering index of greater than 95 and an impressive 7,500 footcandles, the light provides excellent visibility. Its color temperature of 4,400 K ensures clarity and true color representation. The large reflector, featuring 1,040 reflecting facets, guarantees even light distribution and minimal shadows, offering uninterrupted illumination. Additionally, the integrated controls on both sides of the light head allow for convenient adjustments, enhancing user experience and efficiency.

[Read More](#)

SKU: AM-R255-FS



MIDMARK RITTER 355 FLOOR STAND LIGHT

The **Midmark Ritter 355 Single Procedure Light**, a perfect blend of quality and affordability tailored to meet the needs of medical professionals. This light is meticulously counterbalanced, ensuring precise positioning without any drifting, which facilitates accurate procedures. With its 4,200 K color temperature and a large pattern size, the light provides an even distribution that enhances shadow control and delivers superior color rendition, crucial for improved diagnosis. The five-caster base design requires less than two pounds of force to move, allowing for effortless repositioning. Additionally, the handle is fully autoclavable, ensuring high standards of hygiene. The Ritter 355's true color temperature generates less heat, which is beneficial for both patient comfort and prolonged use. Its impressive specifications, including the ability to achieve up to 5,000 foot-candles at a 24-inch focal length, make it an indispensable tool in any medical setting.

[Read More](#)

SKU: AM-R355-FS

VERATHON GLIDESCOPE ROLL STAND

The **Verathon GlideScope Roll Stand** is designed to enhance the usability and accessibility of the GlideScope video laryngoscope in clinical environments. Key features include a sturdy and adjustable height design, allowing healthcare professionals to position the device at an optimal viewing level during procedures. The stand is equipped with smooth-rolling casters, enabling easy mobility between different areas of the facility. Additionally, it includes hanging brackets and a basket for blades and accessories, ensuring that everything needed for airway management is organized and within reach. The durable construction of the roll stand provides stability, making it a reliable companion for emergency situations and routine intubations alike.

Monitor and Baton not included - [Order: AM-0570-0338-0570-0307](#)

[Read More](#)

SKU: AM-0570-RS





INFINIUM CLEARVUE™ VLXL VIDEO LARYNGOSCOPE

The **Infinium® ClearVue™ VLXL video laryngoscope** adopts state-of-the-art image technology and possesses the advantage of a 2.0-megapixel full view camera with 8-inch high-resolution monitor. The ClearVue™ features a rechargeable li-ion battery for mobile use. The HD optics offer a unique anti-fog and high-resolution image. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable the clinician's first-attempt success while minimizing any chance of a complication during the intubation process.

[Read More](#)

SKU: 000.300.0009



GE DASH 4000 PATIENT MONITOR W/ CO2

The **GE Dash 4000 Patient Monitor with CO2 capabilities** is a comprehensive and versatile medical device designed to meet the diverse monitoring needs of healthcare professionals. Featuring a 10.4-inch screen, this full-featured monitor provides clear and detailed displays of vital patient data, including ECG, NIBP, SpO2, temperature, and more. The package includes essential accessories such as ECG cables and leads, a Nellcor finger sensor with an extension, a NiBp cuff and hose, and a roll of paper for the integrated recorder. Additionally, the Dash 4000 comes equipped with the Dash CapnoFlex LF CO2 module, enhancing its functionality by allowing accurate carbon dioxide monitoring. This robust and reliable patient monitor ensures precise and continuous monitoring, making it an invaluable tool in critical care settings.

[Read More](#)

SKU: AM-GE-DASH4000-CO2

Price: \$5,078.00 \$4,570.20





UMF POWER EXAM HI-LO TABLE 5060

The **UMF 5060 procedure table** features a power base, power backrest, two function hand control, ultra comfort seamless top, stirrups, 400 lb. patient capacity and more. One of UMF's most economical power exam tables, the UMF 5060 Power Exam Table lowers to a height of 24 inches (61.0 cm) from the floor and can support a patient weighing up to 400 pounds. It utilizes a two function hand control for both hi-low and power back functions. It can also be used as a chair in its fully upright position. The power exam table also features a stainless steel leg section with a removable upholstered leg cushion. The leg section elevates 10 degrees for added patient comfort.

[Read More](#)

SKU: AM-UMF-5060



STRYKER 1501 RENAISSANCE SERIES STRETCHER WITH IV POLE

The fully refurbished **Stryker 1501 Renaissance Series Stretcher with IV pole** is a versatile and reliable solution for patient transport and care in medical settings. This stretcher boasts a range of features designed to enhance both patient comfort and caregiver convenience. It includes a new 3-inch (8 cm) Enhanced Comfort Mattress, ensuring a restful experience for patients. The stretcher is equipped with 8-inch (20 cm) Omni Surface™ casters, making it easy to maneuver across various floor types. With a 29-inch (74 cm) wide patient surface and a 76-inch long patient surface, it provides ample space for patients of different sizes. The overall length of the stretcher is 83 inches, accommodating the needs of most medical environments. Its dual cylinder pneumatic backrest and dual pedestal hydraulics allow for smooth and precise adjustments, while the unilower pedal simplifies the lowering process. The durable construction includes a Fiberesin backrest and a steel foot-end litter, complemented by fold-down side rails for added safety. The inclusion of an IV pole and adjustable knee gatch further enhances its functionality. With a patient capacity of up to 500 lbs, this stretcher is a robust and dependable choice for healthcare facilities.

[Read More](#)

SKU: AM-S1501



CLEARVUE VL3R™ REUSABLE VIDEO LARYNGOSCOPE

The **Infinium ClearVue™ video laryngoscope** adopts state-of-the-art image technology and possesses the advantage of a 2.0-megapixel full view camera with 4-inch high-resolution monitor. The ClearVue™ features a rechargeable Li-ion battery for extended use. The HD screen offers a unique anti-fog and anti-shatter capability. The overall design is compact and ergonomic. A high-resolution view of the glottis will enable the clinician's first-attempt success while minimizing any chance of a complication during the intubation process.

[Read More](#)

SKU: 000.300.0006



EDAN IM70 PATIENT MONITOR (REFURBISHED)

Embedded with the latest improved designs, iM70 is meant to add splendors to the i Series portfolio and answer for the latest clinical requirements across the applications from general cares to anesthesia monitoring, while keeping connected to hospital systems in all divisions. The high resolution color TFT-LCD together with the up-to-date customizable shortcut menu design brings easy one-touch access to frequent used functions, allowing more attentions to the patients. Without making any noise or causing any dust deposition, the no-fan design of iM70 is especially suitable for night time operations and the divisions which require a dust-free environment.

[Download Brochure](#)

[Read More](#)

SKU: AM-IM70-R

Price: ~~\$1,795.00~~ \$1,615.50



STRYKER M SERIES MODEL 1015 STRETCHER, WITH BIG WHEEL

The **Stryker M-Series SM 204 (1015) with Big Wheel** is a highly adaptable and fully-equipped stretcher that offers a range of mobility solutions to meet individual needs. Designed to prioritize both patient care and clinician comfort, these stretchers are incredibly mobile and user-friendly. With their big wheel feature, they are able to reduce physical strain on clinicians without compromising on the comfort and care provided to the patient.

The stretcher has a variety of additional features, including an overall length of 84 inches and an overall width of 33.5 inches (with siderails up) or 37 inches (with siderails down). The patient surface is 26 inches by 74.25 inches, and the side rails measure 14 inches by 57.75 inches. Additionally, the stretcher has a maximum weight capacity of 700 lbs and includes a 4-inch mattress and an IV pole. Perfect for healthcare settings, the Stryker M-Series SM 204 (1015) with Big Wheel offers a reliable and adaptable solution for safe and comfortable patient transportation.

[Read More](#)

SKU: AM-SM204B



MORTARA ELI 250 ECG MACHINE

Compact and lightweight, the ELI 250c Electrocardiograph provides comprehensive functionality in an easy-to-use, portable device. Multiple connectivity options ensure patient data transfers easily from the ECG to the location of your choice. Experience the benefits of bidirectional communication via internal modem, LAN or wireless LAN with Mortara's E-Scribe™ and Athena products, as well as third-party EMR, HIS, and PACS systems via XML, PDF and DICOM®. Best 10 technology helps you capture the best 10 seconds of ECG data from system memory, reducing the need for repeat ECGs. Wireless data acquisition helps increase workspace flexibility by gathering 12-lead data without the need for a tethered cable.

[Read More](#)

SKU: AM-ELI250



STERIS AMSCO MODEL 3085 SURGICAL TABLE

The Amsco 3085 General Surgical Table is one of the most reliable and popular surgical tables. Recognized as the industry standard due to its sheer versatility and reliability, the Steris Amsco 3085 has been trusted for years by healthcare facilities across the world.

[Download Brochure](#)

[Technical Specs](#)

[Read More](#)

SKU: AM-3085SP



MINDRAY PASSPORT V WITH CO2

Passport V is built on a solid foundation of experience and knowledge. It follows a proven approach to acquisition and presentation of patient data, thereby delivering impressive reliability and performance. Enhanced capability, more optional features and extraordinary value make Passport V a compelling alternative for any healthcare environment requiring continuous physiological monitoring.

[Download Brochure](#)

[Read More](#)

SKU: AM-PASSPORT5-CO2-SS

Price: ~~\$4,321.00~~ \$3,888.90



MIDMARK 623 EXAM TABLE

The **Midmark 623 Barrier-Free Examination Table** allows patients incapable of getting onto an exam table (whether it be the elderly, or young children) to easily get on, increasing the patient's comfort and satisfaction with the patient experience. To meet the needs of America's disabled patients, as well as the physicians who care for them, the Midmark 623 barrier-free exam tables decrease the chance of injury, helping maintain the patient's dignity and reducing their anxiety. This will help increase exam room efficiency, comfort, and safety.

[Read More](#)

SKU: AM-M623



STERIS CAVIWAVE SINGLE CHAMBER ULTRASONIC CLEANER 15 GAL

- Cleaning frequency of 132kHz that more than triples the commonly used 40kHz
- Toe touch cycle start
- Automatic chemistry injection
- Automatically cleans, flushes, lubricates, and dries depending on unit configuration*
- Flush ports allow tank solution reach into internal passages

[Read More](#)

SKU: AM-CAVIWAVE-SC-15

MIDMARK 641 PROCEDURE CHAIR

Midmark 641 Barrier-Free Power Procedures Chair Features:
Narrow patient support surface, slim back section, and tuck away side rails provide maximum patient accessibility and supports proper posture. Powered seat, back, foot, and tilt offer multi-positioning capabilities to improve comfort and reduce likelihood of injury. Premium Comfort System enhances patient comfort for a chair-like feel.

[Read More](#)

SKU: AM-M641



BYRON LYSONIX 3000 LIPOSUCTION UNIT

Features:

- PSI-TEC III Aspirator
- EZ Pump Infiltrator
- Lysonix 3000

[Read More](#)

SKU: AM-B-L3000-EZ-PTIII





TUTTNAUER EZ10K 220V AUTOCLAVE

The **Tuttnauer EZ10k 220V Autoclave** is a cutting-edge sterilization tool renowned for its efficiency and reliability in medical and laboratory settings. Designed with rapid cycle times, it features a fully automatic, user-friendly interface, making it ideal for high-demand environments where time is crucial. Its advanced closed-door drying system ensures instruments are thoroughly dried post-sterilization, preserving their integrity for immediate use. Constructed with durability in mind, the EZ10k includes robust safety features like automatic shut-off and pressure release mechanisms, which provide operator safety and consistent performance. The spacious chamber accommodates a variety of instrument sizes, allowing for diverse sterilization needs while upholding high infection control standards. The Kwiklave model, a quick cycle variant, reduces total sterilization time by 50%, making it perfect for quick instrument turnover or high-volume settings. Operating on a voltage between 220V and 235V, it may require a buck/boost transformer for optimal performance. With a chamber volume of 6 gallons and a power rating of 2200W, the EZ10k is both powerful and versatile, though it does not include a printer.

[Read More](#)

SKU: AM-EZ10K-R



ABBOTT I-STAT 1 GEN 2 HANDHELD BLOOD ANALYZER

Abbott i-STAT 1 is an advanced blood analyzer that provides lab-quality results in minutes for faster patient care. Easy-to-use, this handheld analyzer operates with the advanced technology of single-use i-STAT test cartridges. This point of care (POC) system provides lab-quality results in minutes and allows healthcare professionals to quickly respond to needs and conditions at the patient's side. It delivers lab-quality results significantly quicker than both bench-top systems and the central laboratory.

[Read More](#)

SKU: AM-i-STAT-1-2

Price: ~~\$6,495.00~~ \$5,845.50



MINDRAY DPM6 PACU MONITOR

The Mindray DPM6 patient monitor offers the functionality to meet the demands of the ICU and OR, such as multi-vector Arrhythmia and ST analysis, Cardiac Output (CO), Continuous Cardiac Output (CCO/SVO2), CO2, up to 8 invasive blood pressures, 5-agent multi-gas analysis with auto-ID and data storage.

[Read More](#)

SKU: AM-MR-DPM6

Price: ~~\$1,995.00~~ \$1,795.50



STERIS DJ04 DUAL COMPARTMENT BLANKET WARMER

Steris Model DJ04

[Read More](#)

SKU: AM-DJ04



STERIS CAVIWAVE PRO SINGLE CHAMBER 17 GAL ULTRASONIC CLEANER

The **Steris Caviwave Pro Single Chamber 17 GAL Ultrasonic Cleaner** is a state-of-the-art cleaning solution designed for the meticulous care of delicate instruments. With a remarkable cleaning frequency of 132kHz, it significantly surpasses the standard 40kHz used by most ultrasonic cleaners, ensuring a deeper and more thorough clean. It boasts a robust sonic power of 2000 watts, making it ideal for handling heavy-duty cleaning tasks. The convenience of a toe touch cycle start, coupled with automatic chemistry injection, streamlines the cleaning process, reducing manual intervention. Depending on the unit configuration, the cleaner can automatically clean, flush, lubricate, and dry instruments, making it a versatile addition to any facility. Flush ports are strategically incorporated to allow the tank solution to penetrate internal passages effectively. **The 17-gallon tank can accommodate up to 35 pounds of instruments, including specialized trays like the Da Vinci Xi instrument Tray CRP1723 and Da Vinci® Single Site® and Endowrist® Instruments (CRP1717).** The unit's electrical specifications cater to various setups, supporting both single and three-phase connections. Efficient drainage is ensured through gravity and pump systems, with maximum flow rates of 10 and 20 gpm, respectively, accommodating water temperatures up to 120°F. The cleaner's design, with dimensions of 43-1/2 x 42-1/4 ± 1 x 36 inches, requires a minimum service access clearance of 24 inches on all sides, emphasizing the need for ample space for optimal operation and maintenance.

[Read More](#)

SKU: AM-CAVIWAVE-PRO-SC-17



SKYTRON 6600 SURGICAL TABLE

The **Skytron 6600 general purpose surgical table** provides unparalleled imaging access for today's most demanding cases while lifting 1000 pounds and articulating 600 pounds. The removable Back and Leg sections facilitate the use of a carbon fiber Back and Leg sections to further extend the imaging window. The rotating top of 210° further enhances C-arm and microscope access along with more leg room for a seated surgeon and in many cases, better surgical access.

[Read More](#)

SKU: AM-SKYTRON-6600



GE DINAMAP CARESCAPE V100 VITAL SIGNS MONITOR

The **GE DINAMAP Carescape V100 Vital Signs Monitor** is a robust, portable device designed to help healthcare professionals accurately and efficiently monitor a patient's vital signs, including blood pressure, heart rate, temperature, and SpO2 levels. Engineered for use in a variety of clinical settings, from hospitals and clinics to emergency response environments, this monitor is known for its reliability and ease of use. Its compact design makes it ideal for mobile healthcare providers, allowing for quick assessments of patient health status on the go. The Carescape V100 integrates seamlessly with other GE healthcare systems, facilitating a comprehensive approach to patient care. With its intuitive interface, long battery life, and durable construction, the GE DINAMAP Carescape V100 Vital Signs Monitor is a valuable tool for medical professionals aiming to deliver high-quality care with accuracy and efficiency.

[Read More](#)

SKU: AM-GE-DV100

Price: ~~\$1,136.00~~ \$1,022.40



MEDFUSION 3500 SYRINGE PUMP WITH VERSION 5 SOFTWARE

The **Medfusion 3500 Syringe Pump** is a cutting-edge device designed to meet the needs of various departments within hospitals and medical practices. This refurbished pump features intuitive syringe detection capabilities, enhancing safety and accuracy in drug administration. Known as a "smart pump," it accommodates syringes from seven different manufacturers, ensuring versatility for healthcare providers. Its ergonomic design allows for single-hand loading, making it easy to operate, while its mounting options enable secure placement on IV poles or bassinets. Additionally, the Medfusion 3500 employs **FlowSentry technology**, which further safeguards patient safety and ensures precise drug delivery. With an accuracy of +/- 2% and compatibility with MRI settings, this advanced syringe pump is an essential tool for improving medication administration in any clinical setting. The inclusion of Version 5 Software enhances its functionality, providing healthcare professionals with a reliable and efficient solution for intravenous therapy.

Warranty: 6 Months

[Read More](#)

SKU: AM-MF3500-V5

Price: ~~\$1,795.00~~ \$1,615.50



WELCH ALLYN 6000 SERIES VITAL SIGNS MONITOR W/ CO2

The Connex Vital Signs Monitor (VSM) 67NCXX-B is an advanced, touch screen monitor featuring bright, vivid colors with minimal knobs or buttons for ease of cleaning—perfect for hospital general care floors and for surgery centers and other ambulatory care offices. The Connex VSM helps clinicians dramatically improve patient safety and the time it takes for vitals capture. A wide range of optional configurations help clinicians to reduce/improve common challenges including falls, pressure ulcers, respiratory distress, EMR documentation, early warnings/detection, white coat hypertension and accuracy of vitals information.

[Read More](#)

SKU: AM-WA-67NCXX-B

Price: ~~\$2,776.00~~ \$2,498.40



AARON BOVIE 950 HIGH FREQUENCY GENERATOR

The **Aaron Bowie 950 High Frequency Generator** is a versatile and comprehensive **electrosurgical unit** designed for precision in medical procedures. This state-of-the-art generator comes equipped with a variety of components, including a power cord, a reusable handpiece, and a set of 10 handpiece drapes, catering to the diverse needs of medical professionals. It also includes five disposable active electrodes and a selection of 50 sharp and 50 blunt dermal tips, offering flexibility in surgical applications. For effective patient safety, the device features a reusable dispersive cord and five disposable dispersive electrodes. The unit is user-friendly, with a detailed user's guide and a wall-mount kit for convenient installation. With cutting, blending, and coagulation capacities of up to 60 watts, and a desiccator fulguration capability of 35 watts, the generator ensures precise tissue manipulation. Additionally, its bipolar function offers 30 watts of power, operating within a line voltage range of 120 VAC \pm 10% for the A950™ model and 230-240 VAC for the A950™-220 variant. The device's output frequency ranges from 350-800 kHz, suitable for a wide array of surgical procedures, and it is compatible with a line frequency of 50-60 Hz. Compact in design, the generator measures 26 cm in width, 17.8 cm in height, and 11.4 cm in depth, making it a practical addition to any medical facility.

[Read More](#)

SKU: AM-B950




STERIS CAVIWAVE RINSING AND DRYING ULTRASONIC CLEANER

The refurbished **Steris Caviwave Rinsing and Drying Ultrasonic Cleaner** is a state-of-the-art cleaning system designed to deliver superior performance in medical and laboratory settings. Operating at a cleaning frequency of 132kHz, it significantly surpasses the standard 40kHz frequency, ensuring thorough and efficient cleaning of instruments. This advanced system features a convenient toe touch cycle start and automatic chemistry injection, streamlining the cleaning process. Depending on the unit's configuration, it can automatically clean, flush, lubricate, and dry instruments, making it versatile for various applications. The system is equipped with a **15-gallon (57 L) tank** capable of holding up to 22 pounds (10 kg) of instruments, and its approximate operating weight is 1311 pounds (595 kg). It requires a 208 V, 50/60 Hz, 3-phase (5 wire) electrical connection and is compatible with hot and cold water inputs. The overall dimensions of this fully equipped console model are 39-1/2 x 42-1/4 \pm 1 x 60 inches (1003 x 1073 \pm 25 x 1524 mm), with a minimum service access clearance of 24 inches (610 mm) on all sides,

ensuring easy maintenance and operation.

[Read More](#)

SKU: AM-CAVIWAVE-DC



PHILIPS INTELLIVUE MP50 PATIENT MONITOR

The **Philips IntelliVue MP50 Patient Monitor** offers a portable and adaptable solution for intermediate care environments, combining a bright 12-inch color SVGA display with flexible measurement and data management features. Equipped with a rugged design and comfortable handle, the monitor supports easy in-hospital transfer and extended use with a 5-hour battery life. Its customizable interface includes 10 preset and user-defined screen configurations, enabling clinicians to view data in graphical or numerical formats, trend measurements, and organize on-screen elements to suit specific workflow needs. Featuring touch-enabled operation (on MP50), advanced network integration via portal technology, and compatibility with a range of SpO2 sensors and systems from Philips, Masimo, and Nellcor, the MP50 provides comprehensive clinical support, seamless connectivity, and space-efficient, quiet operation—making it an efficient, user-friendly tool for patient monitoring across a variety of clinical settings.

[Read More](#)

SKU: AM-MP50

Price: ~~\$3,779.85~~ \$3,401.87



GE DASH 4000 PATIENT MONITOR

The **GE Dash 4000** is a comprehensive patient monitoring system designed to deliver accurate and reliable vital measurements in various clinical settings. Featuring a 10.4" screen, this full-featured monitor provides clear and concise data visualization for ease of use. Equipped with ECG, non-invasive blood pressure (NIBP), SpO2, temperature, and a built-in recorder, the Dash 4000 ensures that healthcare professionals have access to critical patient information at a glance. The system includes essential accessories such as ECG cables and leads, a Nellcor finger sensor with extension, a NIBP cuff and hose, and a roll of paper for immediate documentation needs. Its robust design and comprehensive feature set make the GE Dash 4000 an invaluable tool for efficient and effective patient care.

[Read More](#)

SKU: AM-GE-DASH4000

Price: ~~\$1,995.00~~ \$1,795.50



MINDRAY ACCUTORR V VITAL SIGNS MONITOR

The Mindray Accutorr V vital signs monitor is an easy to use, non-invasive blood pressure monitor with more features and capability. Every unit comes equipped with a recorder, lithium ion battery and your choice of SpO2 technology including Masimo SET®, or Nellcor® OxiMax® SpO2. The Accutorr V is capable of storing up to 1200 measurements easily accessible by patient ID through bar code technology. Plus using the Accutorr CS Charting Solution you can conveniently connect the Accutorr V to your facility's EMR. If maximum flexibility, high value, and reliability match your vital signs needs, then look no further than the Accutorr V monitoring solution.

Includes Adult NIBOP Cuff, SpO2 Sensor

[Read More](#)

SKU: AM-ACCUTORR-V

Price: ~~\$922.00~~ \$829.80



GENDRON 1190-231 BARIATRIC STRETCHER

Gendron develops, designs, and manufactures bariatric patient care products for acute care, long-term care, home care, and rehabilitation environments. The Gendron 1190-231 has a patient weight capacity of up to 1000 pounds and a litter top usable patient surface that measures 37" x 80".

[Download Brochure](#)

[Read More](#)

SKU: AM-G1190-132B



STERIS AMSCO QDJ05 SINGLE COMPARTMENT WARMING CABINET

The Steris Amsco Warming Cabinet is the solution for all warming and storage needs within the operating room. This Steris Cabinet is designed as a blanket warmer as well as a warmer and storage container for sterile IV fluids and irrigation supplies. The newly designed system keeps linens and fluids warm, secure and compliant at the proper temperature settings. The Steris QDJ-05 Single Compartment Warming Cabinet is equipped with digital temperature settings to ensure proper the proper warming temperature.

[Download Brochure](#)

[Read More](#)

SKU: AM-QDJ05



HAUSTED VIDEO IMAGING CHAIR

Chair designed specifically for modified barium swallow studies for AP and Lateral imaging as a procedural chair for diagnosis and treatment. The shortened base is easy to maneuver in small rooms with large imaging equipment. The radiolucent backrest adjusts to a full upright position of 90° to clear, unobstructed view for imaging. Inserts are available for pediatric procedures. This unit adjusts easily from sitting to supine and provides a rotating seat section with a stable work surface during patient treatment. Transport, treatment, procedure, and recovery can be performed on the same unit to improve case turn around and stream-line patient handling.

[Download Brochure](#)

[Read More](#)

SKU: AM-HAUSTED-VIC



HILL-ROM TRANSTAR 8010 EYE STRETCHER

Ideal for eye surgery, the Hill-Rom Surgical Stretcher (also referred to as an eye stretcher) offers outstanding mobility, functionality and quality. A surgical rail allows for easy attachment of gas delivery and drape support devices, PACU extenders and arm boards, and surgical hand tables. Articulating head section also provides enhanced patient positioning and versatility.

[Download Brochure](#)

[Read More](#)

SKU: AM-P8010



STRYKER 735 TRANSPORT STRETCHER WITH IV POLE

The **Stryker 735 Transport Stretcher** is a thoughtfully designed piece of medical equipment that provides reliable and efficient patient transport, particularly well-suited for Ambulatory Surgery Centers (ASCs). Its standout features include a fifth-wheel mobility system that enhances traction and cornering, making navigation through tight spaces smoother and more controlled. Safety and stability are prioritized with brakes on all four 8" Omni Surface™ casters, which lock both swivel and rotation, ensuring secure stationary positioning. The stretcher is equipped with a new mattress and a 26" wide patient surface, offering comfort to patients during transit. Practicality is further enhanced with a range of features, such as an integrated storage tray with an oxygen bottle holder, fold-down siderails, and four IV receptacles. The dual-end activated four-wheel, steel-ring brake system, dual cylinder pneumatic backrest, and dual side-mounted foot control hydraulics ensure ease of use for medical staff. With a weight capacity of 500 lbs, the Stryker 735 is robust and versatile, making it an essential asset in any medical transport scenario.

[Read More](#)

SKU: AM-S735-IV



CRITICARE NGENUITY CO2

The Criticare nGenuity 8100E Series patient monitor is ideal for monitoring pulse oximetry (SpO2), non-invasive blood pressure, heart rate, temperature, respiration and ECG during conscious sedation or post-procedure recovery. Information can easily be read on the crisp 10.4 inch color display, while dedicated function keys allow for quick setup and easy use under demanding time constraints. Full compatibility with popular Electronic Medical Records (EMR) systems makes it easy for any facility or office to transfer patient data.

[Download Brochure](#)

[Read More](#)

SKU: AM-8100EP1



BAXTER SIGMA SPECTRUM INFUSION PUMP WITH 6.02 SOFTWARE

The **Baxter Sigma Spectrum Infusion Pump with 6.02 Software** is a compact, user-friendly, patient-centered infusion device designed to enhance safety and efficiency in clinical settings. Equipped with Baxter's Dose Error Reduction Software (DERS), it significantly helps prevent pump-related medication errors through features like titration error prevention, flow alerts, and automatic medication dose updates. Its versatile infusion profiles support both automated ramped delivery and manual multi-step infusions, making it suitable for a wide range of therapies. The pump's compatibility with standard IV sets streamlines setup and reduces free-flow risks, while its rechargeable lithium-ion battery ensures reliable, portable operation without the need for disposable batteries. With a bright, large-type color display, pump-guided programming, and onboard help screens, the Sigma Spectrum offers quick, intuitive operation. Optional wireless communication and Continuous Quality Improvement (CQI) reporting further support therapy monitoring and compliance. Measuring just 5.8" H x 4.2" W x 2.5" D and weighing only 1.6 lbs, this lightweight infusion pump is ideal for busy healthcare environments requiring precise, safe, and efficient medication delivery.

Warranty: 1 Year

[Read More](#)

SKU: AM-BAXTERSS

Price: ~~\$995.00~~ \$895.50

ZOLL R SERIES DEFIBRILLATOR



The R Series ALS continues to meet the needs of ALS providers the world over, and additional features and parameters are being added on a regular basis. Wi-Fi capability is now available. With the addition of a wireless card, R Series can download and transmit code data wirelessly with just the push of a button. Also available is Defibrillator Dashboard, an optional software program that can further reduce your total cost of ownership. Defibrillator Dashboard will automatically notify clinical engineering if a defibrillator is not ready, and even allows for total remote troubleshooting and monitoring of the defibrillator fleet.

[Read More](#)

SKU: AM-ZR3BIPA



MIDMARK RITTER M11 ULTRACLAVE (AUTOCLAVE), M11-022

The refurbished **Midmark Ritter M11 Ultraclave** (Autoclave) offers a compact yet spacious sterilization solution with its **11" x 18" chamber**, touted as the largest among standard countertop sterilizers. This increased capacity makes it ideal for processing larger instrument packs and equipment, improving workflow efficiency. The M11 features an automatic-opening door, facilitating ease of use and promoting rapid drying times, which is crucial for preventing contamination. Its Steam-Flush Pressure Pulse Air Removal System ensures thorough steam penetration for effective sterilization. The self-program controls allow customization of cycle parameters to meet specific instrument processing needs, offering flexibility for various materials and instruments. Pre-programmed controls provide convenient, standardized cycles for routine sterilization tasks. The green LCD display enhances visibility, allowing for easy monitoring of cycle indicators and parameters. The M11 Ultraclave combines capacity, convenience, and customizable control in a compact countertop design.

[Download Brochure](#)

[Read More](#)

SKU: AM-M11-022



RITTER MIDMARK M9, M9-001

The M9 autoclave couldn't be easier and is time tested for reliability. You simply load it - set it - and forget it. Select your cycle and the time is set for you. Press "Start" and the M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam, giving you fast drying time for your instruments.

[Download Brochure](#)
[Read More](#)

SKU: AM-M9-001



ALARIS MEDLEY 8100 INFUSION PUMP

The **Alaris Medley 8100 Infusion Pump Module** is a high-capacity, large-volume infusion device designed to seamlessly integrate into the modular Alaris System, operating in conjunction with the Alaris PC Point of Care Unit (Models 8000 or 8015). Compact and lightweight, measuring just 3.3" W x 8.9" H x 5.5" D and weighing approximately 2.5 lbs, this positive displacement pump provides precise and reliable delivery of fluids across a flow rate range of 0.1 to 999 mL/h, with customizable VTBI up to 9999 mL and adjustable Keep Vein Open rates. Its advanced safety features include Guardrails Suite MX, ultrasonic air-in-line detection, and automatic free-flow protection, ensuring patient safety during complex infusion therapies. The pump supports multiple infusion modes—continuous, intermittent, bolus, and secondary—and features pressure sensing with adjustable occlusion limits from 50 to 525 mmHg, displayed in real time via a dynamic pressure bar graph on the PC Unit. Its compatibility with dedicated Alaris infusion sets, along with environmental resilience to operate between 41°F and 104°F, makes the Alaris 8100 an essential component for hospitals and clinics seeking reliable, safe, and versatile infusion management.

Requires the [Alaris PC Point of Care unit \(Model 8000 or 8015\)](#) to operate

6 Month Warranty

[Read More](#)

SKU: AM-Alaris-8100

Price: \$995.00 ~~\$895.50~~



ALARIS MEDLEY 8015 POINT OF CARE MODULE

The **Alaris Medley 8015 Point of Care Module** is a versatile, modular infusion and monitoring system designed for adult, pediatric, and neonatal care, providing clinicians with a unified interface for programming and managing multiple channels simultaneously. This compact unit, weighing only 7.2 pounds and measuring 6.9" W x 8.8" H x 9" D (including pole clamp), offers wireless connectivity to enhance IV safety and facilitate seamless clinical workflows. It features powerful analytics, exceptional versatility, and the most comprehensive safety protections available through Guardrails® Suite MX software, which double-checks dose, duration, and delivery rate to prevent errors. Its advanced memory system retains configuration parameters for up to 6 months, with module-specific data stored for 8 hours when powered off. Designed for modern healthcare environments, the Alaris Medley 8015 improves patient safety, clinical efficiency, and adherence to best practices across a broad spectrum of care settings.

[Compatible with the Alaris Medley 8100 Infusion Pump.](#)

[Download Brochure](#)

Warranty: 6 Months

[Read More](#)

SKU: AM-Alaris-8015

Price: ~~\$725.00~~ \$652.50



MIDMARK RITTER 75L EXAM CHAIR

The **Midmark Ritter 75L Exam Chair** is a versatile and highly adjustable piece of medical equipment, designed to enhance patient comfort and improve efficiency in clinical settings. This exam chair features a range of movements including table up/down, tilt, foot up/down, and back up/down adjustments, allowing healthcare providers to easily position patients for various examinations and procedures. The chair also includes an adjustable headrest and hideaway adjustable stirrups, catering to a wide array of medical needs. For added functionality, it is equipped with a pull-out cysto pan, table paper holder, and table paper straps, ensuring cleanliness and convenience during use. The chair can be controlled via a foot switch or hand control, providing flexibility for the operator. Accessory side rails offer additional support for attaching medical instruments or accessories. With a width of 28 inches at the seat and up to 26 inches at the back, the chair can extend up to 83 inches in length with both the foot section and headrest extended. It offers a height range from a minimum of 19 inches to a maximum of 40 inches, accommodating a diverse range of patient heights and medical procedures. The chair is robust, with a patient load rating of 450 pounds, and includes a paper roll holder that accommodates an 18-inch by 3.5-inch paper roll, ensuring it meets the practical needs of busy medical environments. Colors may vary.

[Download Brochure](#)

[Read More](#)

SKU: AM-M75L



CONMED SABRE 2400

The ConMed™ Sabre 2400™ is a full-power electrocautery generator to accommodate the needs of all surgical procedures. The Sabre 2400™ easily meets the high demanding performance needs of electrocautery. Features independent Monopolar outputs & enhanced bipolar performance. It is safe, convenient, cost-effective, and can be programmed for standardization.

[Download Brochure](#)

[Read More](#)

SKU: AM-CM-SABRE2400



EDAN SE-1200 AND SE-1200 EXPRESS ECG

Excellent cost performance with built-in high-resolution thermal printer. You can record 12 ECG traces on A4 size and US Letter-size paper.

[Edan SE-1200 Brochure](#)

[Edan SE-1200 Express Brochure](#)

[Read More](#)

SKU: AM-SE-1200EX



ZOLL AED PLUS

The **Zoll AED Plus** is an advanced automated external defibrillator that stands out with its Real CPR Help feature, providing real-time feedback to guide users through effective CPR. Equipped with Biphasic Technology, the device ensures precise and effective defibrillation during life-saving interventions. The AED Plus comes complete with a durable carry case, making it easy to transport and protect. It includes one set of long-lasting lithium batteries and one set of Stat-Padz, ensuring the device is always ready for use in emergencies. To maintain optimal functionality, the AED Plus also features an AED Check Tag for regular maintenance checks. The package is rounded out with an operator's manual, offering comprehensive instructions for both novice and experienced users.

[Download Brochure](#)

[Read More](#)

SKU: AM-ZAEDP

ZIMMER ATS 3000

The **Zimmer A.T.S.® 3000 Automatic Tourniquet System** is a state-of-the-art dual-port, dual-cuff system designed with microprocessor controls and dedicated ports that supply and measure pressure independently. This advanced system features the innovative Limb Occlusion Pressure (LOP) technology, which accurately detects, calculates, and reports the precise cuff pressure required to achieve complete blood occlusion in the operative limb. The A.T.S.® 3000 Tourniquet exemplifies Zimmer's long-standing commitment to safety, reliability, and convenience, boasting a four-hour battery backup to maintain pressure during power loss, along with visual and audible alarms for enhanced safety. Additional features include self-check calibration, a lightweight and portable design, color-coded cuffs ranging from 8 to 42 inches, and a convenient carrying handle. The system also ensures rapid inflation and deflation, independent pressure monitoring, and a cuff lockout feature that prompts users to confirm deflation commands, reducing the risk of complications from sudden cuff deflation. Built-in safety features and positive locking connectors for leak-free inflation make the A.T.S.® 3000 a versatile and dependable choice for use on a table or IV pole.

[Download Brochure](#)

[Read More](#)

SKU: AM-ATS3000



ZIMMER ATS 2000

The built-in safety, reliability, and convenience of the A.T.S.® 2000 Tourniquet frees you to focus on surgery. Dual displays that monitor actual cuff pressure, precise microprocessor controls, and one-hour battery backup help assure patient safety. User-set alarms warn of unacceptable or hazardous conditions.

[Download Brochure](#)

[Read More](#)

SKU: AM-ATS2000



INDEX

A

Aaron Bovie 950 High Frequency Generator 27
Abbott i-STAT 1 Gen 2 Handheld Blood Analyzer 22
Alaris Medley 8015 Point of Care Module 36
Alaris Medley 8100 Infusion Pump 35

B

Baxter Sigma Spectrum Infusion Pump With 6.02 Software 33
Biodex 058-878 30in Wide Top Imaging Table 6
Byron Lysonix 3000 Liposuction Unit 21

C

Cardiac Science Powerheart G5, Semi-Auto, Dual Language 3
ClearVue VL3R™ Reusable Video Laryngoscope 17
ConMed Sabre 2400 37
Criticare nGenuity CO2 32

E

Edan IM70 Patient Monitor (Refurbished) 17
Edan SE-1200 and SE-1200 Express ECG 38
Enthermics EC230L Fluid Warmer 2

G

GE Dash 4000 Patient Monitor 29
GE Dash 4000 Patient Monitor W/ CO2 15
GE Datex Ohmeda CardioCap 5 Monitor 5
GE DINAMAP Carescape V100 Vital Signs Monitor 25
Gendron 1190-231 Bariatric Stretcher 30

H

Hausted Video Imaging Chair 31
Hill-Rom Transtar 8010 Eye Stretcher 31

I

Infinium ClearVue™ VLXL Video Laryngoscope 15

M

Medfusion 3500 Syringe Pump with Version 5 Software 26
Midmark 623 Exam Table 20
Midmark 641 Procedure Chair 21
Midmark Ritter 255 Floor Stand Light 13
Midmark Ritter 355 Floor Stand Light 13
Midmark Ritter 75L Exam Chair 37
Midmark Ritter M11 Ultraclave (Autoclave), M11-022 34
Mindray Accutor 7, NIBP, Masimo, SPO2, Truetemp, Pulse Rate 7
Mindray Accutorr V Vital Signs Monitor 29
Mindray DPM6 PACU Monitor 23
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide (Stainless Steel Base) 12
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide and Carbon Fiber Extension (Stainless Steel

Base) 11

Mindray Passport V with CO2 19

Morgan MEDesign FLTLT-FW (Float-Top, Free-Wheel Configuration) 6

Mortara ELI 250 ECG Machine 18

P

Philips Intellivue MP50 Patient Monitor 28

Philips InvelliVue MP5 Fixed Wall Mount Kit 5

Physio-Control Lifepak 20E Defibrillator 9

R

Ritter Midmark M9, M9-001 35

S

Skytron 6600 Surgical Table 25

Steris Amsco Model 3085 Surgical Table 19

Steris AMSCO QDJ05 Single Compartment Warming Cabinet 30

Steris Caviwave Pro Single Chamber 17 GAL Ultrasonic Cleaner 24

Steris Caviwave Rinsing and Drying Ultrasonic Cleaner 27

Steris Caviwave Single Chamber Ultrasonic Cleaner 15 Gal 20

Steris DJ04 Dual Compartment Blanket Warmer 23

Steris Dual Compartment Blanket Warmer QDJ02 2

Steris Harmony LA 585 LED Dual Head OR Lights 7

Steris InnoWave™ Unity IW523046NPR Ultrasonic Cleaner 4

Stretcher, Stryker Advantage, IV pole 3

Stryker 1501 Renaissance Series stretcher with IV pole 16

Stryker 735 Transport Stretcher with IV Pole 32

Stryker M Series Model 1015 Stretcher, with Big Wheel 18

Stryker 1068 Eye Surgery Stretcher 10

T

Tuttnauer EZ10k 220V Autoclave 22

Tuttnauer 2540EK Kwiklave Electronic Autoclave (Refurbished) 11

U

UMF Power Exam Hi-Lo Table 5060 16

V

Valleylab Force FX 8

Verathon GlideScope Monitor with Baton 8

Verathon Glidescope Roll Stand 14

W

Welch Allyn 6000 Series Vital Signs Monitor w/ CO2 26

Welch Allyn SSM Vital Signs Monitor 10

Z

Zimmer ATS 2000 39

Zimmer ATS 3000 39

Zoll AED Plus 38

Zoll R Series Defibrillator 34

Zoll R Series defibrillator, with SP02, NIBP, EtCO2 4