

TABLE OF CONTENTS

- Scrub Sinks** 2
 - Steris Manual Dual Bay Flexmatic Scrub Sink 3
 - Dual Bay Flexematic Scrubsink, IR 3
 - Steris Amsco Flexmatic Single Bay Scrub Station 3



SCRUB SINKS

Built with safety in mind, these state-of-the-art surgical scrub sinks offer ease, convenience, and utility for all users. Durable and hands-free, Auxo Medical offers the highest quality scrub sinks from leading brands such as [Steris](#) and Mac Medical.

Offering one-, two-, and three-sink options, Auxo Medical has the right scrub sink for any operating room. Whether your facility needs a scrub sink with deep basins, adjustable spouts, easy to use mixing valves, or optimum water temperature and pressure, our surgical scrub sinks can fit your need. Our refurbished, stainless steel, easy to clean sinks offer user-friendly water controls with convenient knee and foot pedals. Auxo Medical provides only the most popular brands that focus on craftsmanship and industry-leading components.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

REFURBISHED



STERIS MANUAL DUAL BAY FLEXMATIC SCRUB SINK

Dual bay Flexmatic scrub station without timer or water control sensor. Each unit is constructed of welded 16 gauge stainless steel with a matte finish for easy cleaning and low maintenance. The built-in thermostatically controlled mixing valve provides optimum water temperatures up to 115° F (46° C)

[Download Brochure](#)

[Read More](#)

SKU: AM-FLEX-DB-M

Price: Request Quote



DUAL BAY FLEXEMATIC SCRUBSINK, IR

Dual bay Flexmatic scrub station without timer or water control sensor. Each unit is constructed of welded 16 gauge stainless steel with a matter finish for easy cleaning and low maintenance. The built-in thermostatically controlled mixing valve provides optimum water temperatures up to 115° F (46° C). Chair carrier Included

Foot pedal or knee activated soap dispensers must be ordered separately (not included)

[Download Brochure](#)

[Read More](#)

SKU: AM-FLEX-DB-IR

Price: Request Quote



STERIS AMSCO FLEXMATIC SINGLE BAY SCRUB STATION

Designed in accordance with Underwriters Laboratories standards to meet all of your scrub sink needs, Amsco®Flexmatic® Scrub Sinks have stood the test of time. These durable sinks feature continuous welded compartments and coved corners constructed entirely of stainless steel for easy cleaning. Durable, thermoformed Streamline scrub stations are made with built-in Microban antimicrobial protection with a seamless easy- to-clean basin. Easy-to-use knee panels make water control quick and convenient. The built-in thermostatically controlled mixing valve provides water temperatures up to 115° F (46° C). The standard “rose” spray head produces a shower-head-like stream of water without any aeration. A unique eyewash station attachment provides convenient operating room access to an emergency eye wash station. When water is turned on with knee panel, one pull of the center knob releases both orange safety caps to provide instant eye flushing. When knob is not activated, faucet functions as a scrub sink.

[Download Brochure](#)

[Read More](#)

SKU: AM-FLEX-SB-M

Price: Request Quote

INDEX

D

Dual Bay Flexematic Scrubsink, IR 2

S

Steris Amsco Flexmatic Single Bay Scrub Station 3

Steris Manual Dual Bay Flexmatic Scrub Sink 2