

TABLE OF CONTENTS

| | |
|---|----|
| Stainless OR Suite Furniture | 3 |
| Soiled Linen Hamper with Foot Pedal, Single-Bag, Metal Lid with "Soiled Linen" Print, Dimensions: 19" W x 19 1/2" D x 37 3/4" H | 4 |
| Stainless Steel Clark Socket | 4 |
| MAC Medical Pass Through Window 37.5" W x 34.5" H (shelf included) | 5 |
| MAC Medical Pass Through Window 37.5" W x 34.5" H (NO shelf) | 5 |
| Detecto Step-On Can, 12 Gallon, 48 Quart, Red | 5 |
| Detecto Step-On Can, 12 Gallon, 48 Quart, White | 6 |
| Prep and Pack Workstations, 30" D x 72" W (table top dim), Fixed Height | 6 |
| Mac Medical Over Operating Tables, Adjustable, 60"W X 24"D | 7 |
| Mac Medical Over Operating Tables, Adjustable, 34"W X 42"D | 7 |
| Mac Medical Over Operating Tables, Adjustable, 48"W X 24"D | 8 |
| Soiled Linen Hamper with Foot Pedal, Single-Bag, Metal Lid with "Soiled Linen" Print, 20" W x 18" L 38" H | 8 |
| Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glass with Full Shelf Double Column Medical Storage | 9 |
| Stainless Steel Instrument Table with Shelf 32" W x 60" L x 40" H, 3 locking casters, 3 drawers under top | 10 |
| MAC Medical 24" W x 34.5" H (no shelf) Dual door Pass-through Window | 10 |
| Set of 6 Stainless Steel Nesting Tables | 11 |
| Prep and Pack Workstations, 30" D x 72" W (table top dim), electronic adjustable height | 11 |
| Mac Medical OR Supply Cabinet | 12 |
| Harloff Quad Column Mobile Cabinet with Key Locking doors with tempered glass with 14 total Cath storage slides | 12 |
| Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glass with Full Shelf Double Column Medical Storage | 13 |
| Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glass, Right Side Exchange Trays, Left side 2 Cath Storage Slides | 13 |
| Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glasswith 7 Cath Storage Slides | 14 |

| | |
|--|----|
| SS Instrument Table with Shelf 24" W x 48-72" L x 34" H, 12" Drop Leaf on either end | 14 |
| SS Instrument Table with Shelf 24" W x 36-60" L x 34" H, 12" Drop-Leaf on either end | 15 |
| Economy Central Sterile Instrument Wrapping Prep and Pack Station 32 x 60 x 35 | 15 |
| Prep and Pack Workstations | 16 |
| Wall Shelves | 16 |
| Carts | 17 |
| Stainless Steel Instrument / Back Tables with Shelf | 17 |
| Step Stools | 18 |
| Solution Stands | 19 |
| Physician Seating | 19 |
| Patient Transfer | 20 |
| OR Table Accessories | 20 |
| Stainless Steel Foot Control Mayo Stand | 21 |
| Kick Buckets | 21 |
| IV Pole Accessories | 22 |
| Stainless Steel 5-leg Space Saving IV Pole | 22 |
| Hampers | 23 |
| Armboards | 24 |
| Arm/Hand Surgery Tables | 24 |
| MAC Supply Storage Cabinets | 25 |



STAINLESS OR SUITE FURNITURE

High quality medical surgical suite furniture is vital to maintain a sterile and organized environment. Auxo Medical offers a variety of Stainless OR Suite Furniture.

Today's operating room requires an enhanced focus on the needs of your patients and staff, as well as safety. The more sterile the environment, the safer the patient during surgical procedures. Proper equipment and furnishings will help reduce recovery times from illness or injury. Overall, this contributes to the overall efficiency of the healthcare facility. Auxo Medical provides wall shelves, kick buckets, hampers, arm-board and solution stands, for your OR suite needs.

It tackles many demands of your operating room. It also allows you and your team to optimize patient care and outcomes.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[**BRAND NEW**](#)

[**RECONDITIONED**](#)

REFURBISHED



SOILED LINEN HAMPER WITH FOOT PEDAL, SINGLE-BAG, METAL LID WITH "SOILED LINEN" PRINT, DIMENSIONS: 19" W X 19 1/2" D X 37 3/4" H

The **Soiled Linen Hamper with Foot Pedal** is an affordable, durable chrome waste solution designed for easy, hands-free operation. Featuring a flip-top lid and a convenient foot pedal, it allows for hygienic, hassle-free disposal of soiled linens. The spacious design accommodates an 18-inch hamper bag and includes a easily cleaned plastic top. Measuring **19 inches wide, 19.5 inches deep, and 37.75 inches high**, this stylish hamper is perfect for healthcare facilities, laundry rooms, or any environment where efficient linen management is essential.

[Read More](#)

SKU: AM-MCM2010



STAINLESS STEEL CLARK SOCKET

Mid Central Medical's Stainless Steel Clark Socket is a robust, all-stainless component designed to adapt to a variety of OR setups. It accepts posts up to 11/16" in diameter and is compatible with all US-made operating room tables. The socket can be attached at rail ends or left unattached as needed, providing flexible integration with existing table configurations.

[Read More](#)

SKU: AM-MCM160



MAC MEDICAL PASS THROUGH WINDOW 37.5" W X 34.5" H (SHELF INCLUDED)

The **MAC Medical Dual Door Pass-Through Window** measures **37.5 inches in width and 34.5 inches in height and includes a shelf**. It is constructed with 3/16" tempered glass and is designed to fit walls with adjustable widths ranging from 4.5" to 6.125". The window features a 4-sided gasket and a dual swing hinge with an adjustable spring. A decontamination side shelf can be added as an optional feature, but it is not included.

[Read More](#)

SKU: AM-PTW0002




MAC MEDICAL PASS THROUGH WINDOW 37.5" W X 34.5" H (NO SHELF)

The **MAC Medical Dual Door Pass-Through Window** is **37.5 inches wide and 34.5 inches high and does not include a shelf**. It features 3/16" tempered glass and is designed for installation in walls with adjustable widths between 4.5" and 6.125". It includes a 4-sided gasket and a dual swing hinge with adjustable spring. A decontamination side shelf is available as an optional feature.

[Read More](#)

SKU: AM-PTW0001



DETECTO STEP-ON CAN, 12 GALLON, 48 QUART, RED

The **Detecto Step-On Can, with a 12-gallon (48-quart) capacity**, is designed for demanding environments such as hospitals, doctors' offices, patient rooms, and laboratories, and is UL listed. Known for their durability, these cans feature a heavy-duty foot pedal and linkage system, along with a mechanism to securely hold bags in place. Available in either baked epoxy finish (**red** or white) or stainless steel, the cans are flame retardant and feature a self-closing lid to promote sanitation. The unitized construction ensures longevity, and each can includes free biohazard and infectious waste adhesive labels for proper waste management. The overall dimensions of the 12-gallon can are 23.5 x 13.188 x 14 inches (60 x 33 x 36 cm).

[Read More](#)

SKU: AM-P-48-R



DETECTO STEP-ON CAN, 12 GALLON, 48 QUART, WHITE

The **Detecto Step-On Can, with a 12-gallon (48-quart) capacity**, is designed for demanding environments such as hospitals, doctors' offices, patient rooms, and laboratories, and is UL listed. Known for their durability, these cans feature a heavy-duty foot pedal and linkage system, along with a mechanism to securely hold bags in place. Available in either baked epoxy finish (red or **white**) or stainless steel, the cans are flame retardant and feature a self-closing lid to promote sanitation. The unitized construction ensures longevity, and each can includes free biohazard and infectious waste adhesive labels for proper waste management. The overall dimensions of the 12-gallon can are 23.5 x 13.188 x 14 inches (60 x 33 x 36 cm).

[Read More](#)

SKU: AM-P-48



PREP AND PACK WORKSTATIONS, 30" D X 72" W (TABLE TOP DIM), FIXED HEIGHT

The **Mac Medical Prep & Pack Workstation** is designed with a focus on ergonomics and features a modular design for flexible configuration. It is constructed entirely from stainless steel and includes an integrated footrest for added comfort. The workstation measures **30" deep, 72" wide, and 81" high** (to the top of the back frame), with a **table top height of 36"**. It is equipped with 5" casters, two of which swivel for easy maneuverability and two that lock to ensure stability.

[Read More](#)

SKU: AM-TWS-3072-B



MAC MEDICAL OVER OPERATING TABLES, ADJUSTABLE, 60"W X 24"D

The **Mac Medical Over Operating Table** is an adjustable table designed for use in operating rooms, focusing on ergonomic design and functionality. Measuring **60" wide and 24" deep**, this table features a manually adjustable height to accommodate different procedures and user preferences. It is equipped with 3" casters, including two swivel casters for easy movement and two locking casters to ensure stability during use.

[Read More](#)

SKU: T0082-6024



MAC MEDICAL OVER OPERATING TABLES, ADJUSTABLE, 34"W X 42"D

The **Mac Medical Over Operating Table** is an adjustable table designed for operating room use, prioritizing ergonomic functionality. With dimensions of **34" wide and 42" deep**, this table offers a manually adjustable height to suit various procedural needs. It is equipped with 3" casters, two of which swivel for easy maneuvering and two that lock to ensure stability during use.

[Read More](#)

SKU: T0082-3442



MAC MEDICAL OVER OPERATING TABLES, ADJUSTABLE, 48"W X 24"D

The **Mac Medical Over Operating Table** is an adjustable, ergonomically designed table intended for use in operating rooms. Measuring **48" wide and 24" deep**, this table features a manually adjustable height ranging from 36.5" to 48.25". It is equipped with 3" casters, including two swivel casters and two locking casters, for easy maneuverability and secure placement.

[Read More](#)

SKU: T0082-4824



SOILED LINEN HAMPER WITH FOOT PEDAL, SINGLE-BAG, METAL LID WITH "SOILED LINEN" PRINT, 20" W X 18" L 38" H

Durable, soiled linen hamper with open-frame designed to help reduce heavy upward lifting
Hands-free lid opens easily with a foot pedal and shuts automatically
Rolled bags fit in convenient wallaby pouch
Intended to hold 29" x 43" liners

[Read More](#)

SKU: AM-MDS80529



HARLOFF DOUBLE COLUMN MOBILE CABINET WITH KEY LOCKING DOORS (2) WITH TEMPERED GLASS WITH FULL SHELF DOUBLE COLUMN MEDICAL STORAGE

The **Harloff Double Column Mobile Cabinet** signifies the pinnacle of medical storage solutions, offering unmatched durability and a professional grade construction tailored for high-demand medical settings. Its diverse locking mechanisms, including cutting-edge options like badge entry and WiFi capabilities, set a new standard for security in the industry. This cabinet combines functional elegance with practical design, featuring sleek door handles, hinged doors that can swing up to 270 degrees for easy access, and a depth of 24 inches to accommodate substantial storage needs. With the flexibility to customize through an extensive selection of trays, baskets, shelves, or accessories, the cabinet's interior design supports a highly organized system with 47 insertion points per column. Innovations such as modular trays with safety features, optional dimming LED lights for enhanced visibility in surgical environments, and ultra-smooth 5" Tente® casters underscore its advanced functionality. Delivered fully assembled for immediate use, this cabinet is a testament to Harloff's commitment to quality, offering a robust 12-Year Limited Manufacturer's Warranty. Its dimensions and weight reflect a substantial yet efficiently designed unit capable of seamless integration into any medical facility, ensuring that every tool and supply is stored with the utmost care and accessibility.

[Read More](#)

SKU: MSPM82-20GK-MS-SH2ND-24-MS-SH-DIVPK-24



STAINLESS STEEL INSTRUMENT TABLE WITH SHELF 32" W X 60" L X 40" H, 3 LOCKING CASTERS, 3 DRAWERS UNDER TOP

This **Stainless Steel Instrument Table with Shelf**, measuring an ample 32 inches in width, 60 inches in length, and standing at a height of 40 inches, is meticulously designed to enhance the efficiency and organization within medical, laboratory, or industrial settings. It features robust 3-inch locking casters, ensuring both mobility and stability when required. This table is not just about its significant surface area; it also boasts three drawers located directly under the top surface, providing secure and easily accessible storage for tools, instruments, and supplies. The inclusion of a shelf further augments its storage capabilities, allowing for an organized separation of equipment and essentials. Crafted from stainless steel, this table combines durability with the ease of maintenance, making it an invaluable asset in any high-demand environment.

[Read More](#)

SKU: MCM326040U-3DR



MAC MEDICAL 24" W X 34.5" H (NO SHELF) DUAL DOOR PASS-THROUGH WINDOW

MAC Medical 24" W x 34.5" H (no shelf) Dual door Pass-through Window

[Read More](#)

SKU: AM-PTW0003



SET OF 6 STAINLESS STEEL NESTING TABLES

Set of 6 Stainless Steel Nesting Tables

[Read More](#)

SKU: AM-MCM560



PREP AND PACK WORKSTATIONS, 30" D X 72" W (TABLE TOP DIM), ELECTRONIC ADJUSTABLE HEIGHT

Designed with ergonomics in mind, the Prep & Pack Workstation features a modular design, allowing for multiple configurations and adjustability for all body types.

[Read More](#)

SKU: AM-TWS-3072-EAH



MAC MEDICAL OR SUPPLY CABINET

18"D X 30"W X 64"H
STANDARD FEATURES

[Download Brochure](#) □

[Read More](#)

SKU: AM-IC83064



HARLOFF QUAD COLUMN MOBILE CABINET WITH KEY LOCKING DOORS WITH TEMPERED GLASS WITH 14 TOTAL CATH STORAGE SLIDES

Qty x2 MSCATH7-8-SHLF ND - Catheter Slide Shelf with 7 Cath Slides, Double Column width, 8 Hooks per Slide for a Total of 56 Hooks, for MedStorMax Cabinets,

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

[Read More](#)

SKU: MSPM84-20GK-MSCATH7-8-SHLF ND



HARLOFF DOUBLE COLUMN MOBILE CABINET WITH KEY LOCKING DOORS (2) WITH TEMPERED GLASS WITH FULL SHELF DOUBLE COLUMN MEDICAL STORAGE

Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications

MS-SH2ND-24

Shelf for Full Double Column Medical Storage Cart, Flat with Divider Slots, with no Adjustable Dividers included, full 24" depth (variant of MH-SH2)

MS-SH-DIVPK-24

MS-SH SHELF DIVIDER, for MS-SH1-24 and MS-SH2-24 Shelves, Pack of 5

Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

[Read More](#)

SKU: MSPM82-00GK-MS-SH2ND-24-MS-SH-DIVPK-24



HARLOFF DOUBLE COLUMN MOBILE CABINET WITH KEY LOCKING DOORS (2) WITH TEMPERED GLASS, RIGHT SIDE EXCHANGE TRAYS, LEFT SIDE 2 CATH STORAGE SLIDES

Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications
Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

[Read More](#)

SKU: MSPM82-00GK-MSCATH2-8-DIVS HLF-81030S-81032S-81031S



HARLOFF DOUBLE COLUMN MOBILE CABINET WITH KEY LOCKING DOORS (2) WITH TEMPERED GLASS WITH 7 CATH STORAGE SLIDES

Unmatched quality of construction — built to last using professional grade design and manufacturing for medical-grade applications
MS-SH-DIVPK - MS-Shelf Divider for MSCATH, MS-SH1, MS-SH2, Shelves, Pack of 5
Warranty: 12-Year Limited Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

[Read More](#)

SKU: MSPM82-00GK-MSCATH7-8-SHLF ND MS-SH-DIVPK

SS INSTRUMENT TABLE WITH SHELF 24" W X 48-72" L X 34" H, 12" DROP LEAF ON EITHER END



Tired of large tables in your OR? Not to worry anymore, with 2 flip up extensions, these tables are essentially 3 in 1! Get a quote on one to stop storing and start saving in your OR today!
Less space to store. Fit to the size case you are performing. No expensive custom drapes needed.

[Read More](#)

SKU: AM-MCM507-DPL



SS INSTRUMENT TABLE WITH SHELF 24" W X 36-60" L X 34" H, 12" DROP-LEAF ON EITHER END

Tired of large tables in your OR? Not to worry anymore, with 2 flip up extensions, these tables are essentially 3 in 1! Get a quote on one to stop storing and start saving in your OR today!

Less space to store. Fit to the size case you are performing. No expensive custom drapes needed.

[Read More](#)

SKU: AM-MCM505-DPL



ECONOMY CENTRAL STERILE INSTRUMENT WRAPPING PREP AND PACK STATION 32 X 60 X 35

Central Sterile Instrument Wrapping Station

This stainless steel instrument wrapping station has extra large plastic bins. With the stainless steel top and rack for wrapping instruments in sterile wrap, this table is everything you need. Ergonomic 35" work height and heavy duty 4" caster that lock in place ensure this table will last for many years. Standard size is 32" x 60" x 35" high. Can be easily adjusted in size to meet your needs.

[Read More](#)

SKU: AM-MCM596

PREP AND PACK WORKSTATIONS



Designed with ergonomics in mind, the Prep & Pack Workstation features a modular design, allowing for multiple configurations and adjustability for all body types. Choose from the many optional features including adjustable height (manual or electronic), overhead light, magnification light, CPU holder, electrical strip, and more.

[Read More](#)

SKU: AM-TWS-3060-EAH



WALL SHELVES

We offer a wide variety of Wall Shelves for your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



CARTS

We offer a wide variety of Carts for your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



STAINLESS STEEL INSTRUMENT / BACK TABLES WITH SHELF

Mid Central Medical's Stainless Steel Instrument / Back Tables with Shelf are built to endure in the demanding OR environment. Made from heavy-duty 16-gauge #304 stainless steel, these fully welded tables include a lower shelf for extra storage and come with 3" non-marking neoprene casters as standard. Key design features include an 1 1/2" revealed edge for added rigidity, thicker 16-gauge top and shelf (0.0625") for reduced deflection, welded reinforcing supports on longer models with sound-dampening cushions, and square legs for stronger welds and durability. Ships fully assembled and designed with no open cavities to collect debris, ensuring easy cleaning and glove-friendly edges.

Available in a variety of sizes:

- 16"W x 20"L x 34"H
- 16"W x 30"L x 34"H
- 18"W x 33"L x 34"H
- 20"W x 20"L x 34"H
- 20"W x 36"L x 34"H
- 24"W x 36"L x 34"H
- 20"W x 48"L x 34"H
- 24"W x 48"L x 34"H
- 24"W x 60"L x 34"H
- 24"W x 72"L x 34"H
- 24"W x 30"L x 34"H

[Read More](#)

SKU: AM-MCM507

STEP STOOLS

We offer a wide variety of Step Stools for your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



SOLUTION STANDS

We offer a wide variety of Solution Stands for your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



PHYSICIAN SEATING



We offer a wide variety of Physician Seating your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: AM-MCM862

PATIENT TRANSFER



We offer a wide variety of Patient Transfer your Stainless OR Suite needs. Please contact us or request a quote for more information

[Read More](#)

SKU: AM-BF213

OR TABLE ACCESSORIES

We offer a wide variety of OR Table Accessories your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



STAINLESS STEEL FOOT CONTROL MAYO STAND

Mid Central Medical's Stainless Steel Foot Control Mayo Stand provides a versatile, cleanable tray support that extends over the bed or table for convenient access to instruments and supplies near the surgical site. Built in solid stainless steel, it features removable trays for easy cleaning, brass inserts for smoother height adjustment, and a low-profile base that slides easily under an operating room table. The foot-operated lock adjusts the stand's height from 36" to 62", and it rides on 2" 'sayno2rust' composite casters to minimize rust in the OR. **Three tray sizes are available** to suit different needs.

[Read More](#)

SKU: AM-MCM751





KICK BUCKETS

We offer a wide variety of Kick Buckets for your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



IV POLE ACCESSORIES

We offer a wide variety of IV Pole Accessories- both new and refurbished- for your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



STAINLESS STEEL 5-LEG SPACE SAVING IV POLE

Mid Central Medical's Stainless Steel 5-leg Space Saving IV Pole is designed to maximize OR floor space without sacrificing stability. It features a compact 16" diameter base weighing 21 pounds, supported by 3" scuff-resistant casters for smooth mobility. The stainless steel pole measures 1" in diameter and adjusts from 52" to 94" with a no-lose thumb knob for easy height customization. The base comes in five color options to help departmentalize equipment and includes rubber bumpers on the base legs to protect walls and surroundings. It is **available with 2-Hook, 4-Hook, or 6-Hook Tops** to fit varying IV needs.

[Read More](#)

SKU: AM-MCM222





HAMPERS

We offer a wide variety of Hampers your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: AM-8510



ARMBOARDS

We offer a wide variety of Armboards your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A

ARM/HAND SURGERY TABLES

We offer a wide variety of Arm/Hand Surgery Tables your Stainless OR Suite needs. Please contact us or request a quote for more information.

[Read More](#)

SKU: N/A



MAC SUPPLY STORAGE CABINETS

A variety of stations are available with plenty of features and options.

[Download Brochure](#)

[Read More](#)

SKU: N/A



INDEX

A

Arm/Hand Surgery Tables 23

Armboards 22

C

Carts 15

D

Detecto Step-On Can, 12 Gallon, 48 Quart, Red 4

Detecto Step-On Can, 12 Gallon, 48 Quart, White 4

E

Economy Central Sterile Instrument Wrapping Prep and Pack Station 32 x 60 x 35 14

H

Hampers 22

Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glass with Full Shelf
Double Column Medical Storage 7, 11

Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glass, Right Side
Exchange Trays, Left side 2 Cath Storage Slides 12

Harloff Double Column Mobile Cabinet with Key Locking doors (2) with tempered glass with 7 Cath
Storage Slides 12

Harloff Quad Column Mobile Cabinet with Key Locking doors with tempered glass with 14 total Cath
storage slides 11

I

IV Pole Accessories 20

K

Kick Buckets 20

M

MAC Medical 24" W x 34.5" H (no shelf) Dual door Pass-through Window 9

Mac Medical OR Supply Cabinet 10

Mac Medical Over Operating Tables, Adjustable, 34"W X 42"D 6

Mac Medical Over Operating Tables, Adjustable, 48"W X 24"D 6

Mac Medical Over Operating Tables, Adjustable, 60"W X 24"D 5

MAC Medical Pass Through Window 37.5" W x 34.5" H (NO shelf) 3

MAC Medical Pass Through Window 37.5" W x 34.5" H (shelf included) 3

MAC Supply Storage Cabinets 23

O

OR Table Accessories 19

P

Patient Transfer 18

Physician Seating 18

Prep and Pack Workstations 14

Prep and Pack Workstations, 30" D x 72" W (table top dim), electronic adjustable height 10

Prep and Pack Workstations, 30" D x 72" W (table top dim), Fixed Height 5

S

Set of 6 Stainless Steel Nesting Tables 9

Soiled Linen Hamper with Foot Pedal, Single-Bag, Metal Lid with "Soiled Linen" Print, 20" W x 18" L 38"
H 7

Soiled Linen Hamper with Foot Pedal, Single-Bag, Metal Lid with "Soiled Linen" Print, Dimensions: 19" W
x 19 1/2" D x 37 3/4" H 2

Solution Stands 17

SS Instrument Table with Shelf 24" W x 36-60" L x 34" H, 12" Drop-Leaf on either end 13

SS Instrument Table with Shelf 24" W x 48-72" L x 34" H, 12" Drop Leaf on either end 13

Stainless Steel 5-leg Space Saving IV Pole 21

Stainless Steel Clark Socket 2

Stainless Steel Foot Control Mayo Stand 19

Stainless Steel Instrument / Back Tables with Shelf 16

Stainless Steel Instrument Table with Shelf 32" W x 60" L x 40" H, 3 locking casters, 3 drawers under
top 8

Step Stools 17

W

Wall Shelves 15