

TABLE OF CONTENTS

Surgical Lighting	4
Amico ICE Series Surgical Light Disposable Handle Cover for Non-Sterilizable Handle (1 box - 2pcs/pkg, 80pkg/box)	5
Amico Mira Disposable Handle Covers Box of 50	5
Nuvo/DRE/Medical Illumination Surgical and Exam Light Disposable Sterile Handle Covers Case of 25	5
Amico Clarity 50 Exam Light, Roll Stand	6
Steris Harmony LA 585 LED Single Head OR Light	6
Luxtec 90600 Integra DUO LED Headlamp package	7
Amico Leo Series MG 6in Clinical Recessed Downlight	7
Amico Nuvo Vu Dual Head LED Surgical Light	8
Steris Dual Center Mount vLED Surgical Light	8
Steris Harmony LA 585 LED Dual Head OR Lights	9
Amico iCE LED Surgical Light - Dual Ceiling, 30M/30M Low Ceiling Height Model	9
Midmark Ritter 255 Floor Stand Light	10
Midmark Ritter 355 Floor Stand Light	11
Skytron ST9 Single Floor Mount	11
Amico Mira LED 65, Ceiling Mount, Double, Standard Arm	12
Amico Mira LED 65, Ceiling Mount, Single, Standard Arm	13
Amico Vega 80, Ceiling Mount Exam Procedure Light	13
Steris Harmony Air M Series Dual Head LED Surgical Light	13
Spectra G8 Hospital Kit; G8 Headlight, 6 Batt, 4-Bay Charger, Wall Plug/USB 45W, Padding Conv Kit, Carry Case	14
Spectra G8 Standard Kit; G8 Headlight, 2 Batt., 45W Wall Plug/USB 45W, Padding Conv Kit, Carry Case	14
Spectra G6 Hospital Kit; G6 Headlight, 6 Batt, 4-Bay Charger, 45W Wall Plug/USB, Padding Conv Kit, Carry Case	15
Spectra G6 Standard Kit; G6 Headlight, 2 Batt., 45W Wall Plug/USB 45W, Padding Conv Kit, Carry Case	15
Spectra G3 Hospital Kit; G3 Headlight, 6 Batt,4-Bay Charger, 45W Wall Plug/USB, Padding Conv Kit, Carry Case	16
Spectra G3 Standard Kit; G3 Headlight, 2 Batt., 45W Wall Plug/USB, Padding Conversion Kit, Carrying Case	16
Integra Luxtec MLX 300W Light Source	17

Luxtec Integra LED Headlamp Package	17
Skytron Aurora LED 2 Head Surgical Light System	18
Skytron Aurora LED Single Surgical Light System	18
Steris LED385 Exam Light on Floor Stand	19
Amico Vega 80, Floor Stand Exam Procedure Light	19
Amico Nova Motorized Exam Pair 4000K, Single Keypad	20
Amico Mira LED 90, Floor Stand, Single	20
Amico iCE LED Surgical Light -Triple Ceiling, 30M/30M/30M	21
Amico Single Wall Control Panel	21
Amico iCE LED Surgical Light - Dual Ceiling, 30M/Single Monitor Holder	22
Amico iCE LED Surgical Light - Single Ceiling, 30M	23
Amico Mira LED 90, Ceiling Mount, Double, Standard Arm	23
Amico Mira LED 90, Ceiling Mount, Single, Standard Arm	24
Amico Mira LED 90, Ceiling Mount, Single, Extended Arm	24
Amico Gamma LED 60, Floor Stand, Single	25
Amico Gamma LED 50, Floor Stand, Single	25
Amico Mira LED 50, Floor Stand, Single	26
Amico Mira LED 50, Ceiling Mount, Single, Extended Arm	27
Amico Mira LED 65, Ceiling Mount, Single, Extended Arm	27
Nuvo/Medical Illumination Shuttle Track Dual Trolley	27
Nuvo/Medical Illumination Shuttle Track Single Trolley	28
Steris Harmony Air G Series Dual Head LED Surgical Light	28
Amico iCE LED Surgical Light - Dual Ceiling, 30M/30M	29
Amico iCE LED Surgical Light	29
Nuvo VU Triple Head Surgical w/ Dual Light, One Camera, and One Monitor Arm	30
NUVO Vu LED Surgical Light	31
Nuvo Vu OR Lighting Touch Screen Tablet Control (Optional)	31
NUVO Vistor PRO LED Single Head Ceiling Procedure Light	31
Nuvo Vistor EX LED Examination Light Wall Mount	32
Nuvo Vistor EX Portable LED Examination Light Floor Mount	32
Steris VLED Single Head Ceiling Mounted Surgical Light	33
Steris Harmony LC 500 OR Lights	33
Ritter Midmark 354 Light	34
Skytron ST9 Single Ceiling Mount	34
Ritter Midmark 355 Procedure Light	35
Dual Vistor PRO LED Ceiling Mount	35
NUVO Vistor MSLED Dual Arm Lighting System	36
Steris SQ240	36

Steris Harmony Series 37
Skytron Stellar Series 37
Skytron Aurora Series 38



SURGICAL LIGHTING

Exam lights are essential when operating on and examining patients. Auxo Medical offers energy-efficient wall mounted operating room and medical exam lights to ensure you have the illumination you need.

Stocking the highest quality operating and surgical lights from [Steris](#), [Skytron](#), [NUVO](#), and other popular brands. Auxo Medical understands the importance of clarity during procedures. In fact we offer brighter, whiter, and cooler technology, many of our operating room lights omit very little heat. These exam lights offer sterile controls on handles. With long life expectancy and great depth of illumination, our surgical lights ensure your medical practitioners are able to work effectively.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

REFURBISHED



**Amico ICE Series Surgical
Light Disposable Handle
Cover for Non-Sterilizable
Handle (1 box – 2pcs/pkg,
80pkg/box)**

**AMICO ICE SERIES SURGICAL LIGHT
DISPOSABLE HANDLE COVER FOR NON-
STERILIZABLE HANDLE (1 BOX - 2PCS/PKG,
80PKG/BOX)**

The **Amico ICE Series Surgical Light Disposable Handle Covers** are designed to provide a hygienic barrier for non-sterilizable surgical light handles, ensuring infection control during procedures. Each box contains 80 packages, with 2 disposable covers per package, totaling 160 individual pieces. These sterile, single-use covers are easy to apply and remove, helping healthcare providers maintain a sterile environment, prevent cross-contamination, and comply with infection prevention protocols during surgical lighting care.

[Read More](#)

SKU: LT-CN-P-DHC-BOX

Price: \$595.00



**Amico Mira Series
Disposable Handle
Covers Box of 50**

**AMICO MIRA DISPOSABLE HANDLE COVERS
BOX OF 50**

The **Amico Mira Disposable Handle Covers (Box of 50)** are sterile, single-use covers designed to provide a hygienic barrier over medical device handles during procedures. These covers help prevent cross-contamination and maintain infection control standards, ensuring safety for both patients and healthcare providers. Compatible with the Amico Mira system, they offer a convenient solution for maintaining cleanliness and efficiency in clinical environments, supporting compliance with sterilization protocols.

[Read More](#)

SKU: L-X-MLED-HNDL-CVR-50

Price: \$199.95



**Nuvo/DRE/Medical
Illumination Surgical
and Exam Light
Disposable Sterile
Handle Covers
Case of 25**

**NUVO/DRE/MEDICAL ILLUMINATION
SURGICAL AND EXAM LIGHT DISPOSABLE
STERILE HANDLE COVERS CASE OF 25**

The **Nuvo/DRE/Medical Illumination Surgical and Exam Light Disposable Sterile Handle Covers (Case of 25)** are single-use, sterile covers designed to maintain hygiene and prevent cross-contamination during surgical and examination procedures. Compatible with a broad range of Nuvo, DRE, Medical Illumination, and Bovie surgical and exam lights—including Nuvo/Medical Illumination VU, System 2, MSLED, Maxx Luxx, Vision, VistOR PRO, PR, EX, and various MI models—these covers provide a secure barrier over light handles. They are essential for infection control, ensuring safe, efficient, and compliant workflows in medical environments.

[Read More](#)

SKU: AM-NV0001810

Price: \$185.00



**AMICO CLARITY 50 EXAM LIGHT, ROLL
STAND**

The **Amico Clarity 50 Exam Light with Roll Stand** is a lightweight, versatile lighting solution for any hospital room. Equipped with seven powerful LEDs, it delivers up to 50,000 lux of bright, color-accurate illumination at a cool 4400K and a CRI of up to 93, ensuring safe and precise diagnoses. Clinicians can tailor the light with three additional color levels (3300K–4400K) to suit different exam types, and five-stage dimming lets you quickly adjust brightness from a convenient keypad on the back of the head. The Clarity 50 offers flexible arm configurations and mounting options, plus standard and spotlight configurations with a 7-inch head, making it an optimal choice for any examination setting.

[Read More](#)

SKU: AM-L-CL50-AR-FS-00-00



STERIS HARMONY LA 585 LED SINGLE HEAD OR LIGHT

Steris Harmony LA 585 LED Single Head OR Light is a high-performance surgical illumination solution designed for exceptional visualization in the operating room. It delivers brighter, cooler white light with up to 41% more illuminance than competing LEDs, a 4400K color temperature ($\pm 300\text{K}$), and a large bright spot up to 160,000 lux, all while maintaining superior color rendering with a CRI up to 96 and robust R9 values for saturated reds. The HD-capable suspension supports virtually any video signal and fiber-optic transmission to Steris HD monitors, enabling razor-sharp images and vivid colors. Designed for flexibility, the system accommodates a wide range of HD signals (DVI, HD-SDI, RGBHV, etc.), supports future upgrades with modular lightheads and monitors, and benefits from efficient cable management for reliable long-term operation.

[Read More](#)

SKU: AM-LA585-SCM



LUXTEC 90600 INTEGRA DUO LED HEADLAMP PACKAGE

The **Luxtec 90600 Integra DUO LED Headlamp package** delivers brighter, whiter visualization with enhanced color consistency for precise tissue differentiation across the surgical field. Featuring a large beam for comprehensive illumination and a design emphasis on comfort, the ergonomic headband, Cranial Cushion dual-layer padding, and Conforming Occipital Architecture provide balanced weight distribution and reduced cranial pressure. The system offers stable brightness with reduced light fall-off and a 48,000-hour LED life, plus audible low-battery notifications and a 4-hour maximum power delivery from the included batteries. The package includes the 90600 DUO Headlight, two 90622 Li-Ion batteries, a 90621 battery holster, and a nylon storage bag for portability and storage.

[Read More](#)

SKU: AM-LINTEGRA-LEDH-DUO-90600



AMICO LEO SERIES MG 6IN CLINICAL RECESSED DOWNLIGHT

The **Amico Leo Series MG 6" Clinical Recessed Downlight** is a hospital-grade recessed lighting fixture engineered for the demanding needs of operating rooms and clinical spaces. Delivering 1000 to 4000 lumens with a wide beam and glare-free performance, it supports 120/277 VAC, 15-50 W and 0-10 V dimming to 1% for precise light control. With standard 90 CRI (80 CRI available on request) and selectable CCTs from 3500K to 5000K, the Leo MG ensures accurate color rendering across procedures. It features a finsned heat sink for efficient thermal management, is serviceable from beneath the ceiling, and is available in sealed and non-sealed configurations for versatile installation. The fixture is UL-listed for dry and damp locations (ETL), CSA-compliant, and backed by a five-year limited warranty, making it a reliable choice for maintaining optimal illumination in operating rooms and other clinical

[Read More](#)

SKU: L-SDL6-UNV-40-L50-A-S-0-00-0-00-00



AMICO NUVO VU DUAL HEAD LED SURGICAL LIGHT

The **Amico Nuvo Vu Dual Head LED Surgical Light** utilizes innovative Cross Focus Technology (CFT), allowing independent movement of the LEDs to adjust the pattern size from 9.5" to 14" while maintaining consistent light output across all settings. Delivering an intensity of 160,000 lux, the Nuvo Vu offers three color temperature settings (3800°K, 4800°K, and a Green Mode for minimally invasive procedures), a high Color Rendering Index (CRI) of 95, and an exceptional R9 value of 95 for accurate red tissue differentiation. With near-zero radiated heat, a depth of illumination of 32.3", five intensity settings, an ambient light mode, and controls conveniently located on the sterile handle, the Nuvo Vu provides surgeons with exceptional control and comfort. The LEDs boast a long life expectancy of 50,000 hours, ensuring reliable and energy-efficient performance.

[Read More](#)

SKU: AM-NVSD3D3



STERIS DUAL CENTER MOUNT VLED SURGICAL LIGHT

The **Steris Dual Center Mount Harmony vLED Surgical Lighting System** is a modular and configurable solution designed to provide high-quality illumination for a wide range of surgical and diagnostic applications. With a maximum central illuminance of 160,000 lux, a fixed or adjustable pattern size, and a color temperature of 4400K, this system delivers exceptional visibility and clarity for procedures ranging from simple diagnostics in outpatient centers to complex surgeries in acute care suites. Boasting a high Color Rendering Index (CRI) of 96 and a long lifespan of 30,000 hours, the Harmony vLED ensures accurate tissue discrimination and reliable performance in operating rooms, exam rooms, emergency departments, and intensive care units. Its efficient heat-to-light ratio of 3.4 mW/m²-lx contributes to a comfortable working environment for surgical teams.

[Read More](#)

SKU: AM-SVLED-DCM



STERIS HARMONY LA 585 LED DUAL HEAD OR LIGHTS

The **Steris Harmony LED 585** surgical lighting and visualization system is engineered to enhance surgical precision and patient outcomes. Its high-performance optics deliver exceptional surgical site illumination with brighter, whiter light, superior color rendition, and virtually shadow-free visibility thanks to 84 overlapping beams. Beyond illumination, the system's HD-capable suspension provides stunning clarity for video signals, accommodating various HD technologies and allowing for future upgrades and flexibility. By reducing energy consumption, eliminating bulb changes, and minimizing beam heat, the Harmony LED 585 creates a cooler, more comfortable operating room environment while offering years of reliable performance.

[Read More](#)

SKU: AM-LA585-DCM



AMICO ICE LED SURGICAL LIGHT - DUAL CEILING, 30M/30M LOW CEILING HEIGHT MODEL

The **Amico iCE LED Surgical Light** offers advanced illumination solutions for surgical environments, **specifically designed for spaces with low ceiling heights**. Featuring dual ceiling-mounted light-heads, the **iCE 30M LED** model boasts an impressive 30-inch diameter with 80 LEDs, delivering a powerful 160,000 lux of bright, cool white light. This system ensures optimal visibility with an adjustable light pattern ranging from 7.6 to 12 inches, catering to various surgical needs. With a minimum LED lifespan of 50,000 hours, the lights are both energy-efficient and long-lasting. The iCE LED system also includes smooth surfaces for easy cleaning and comes standard with sterile, disposable covers for the adjustment handle and in-light HD camera, enhancing hygiene and ease of use. The system's exceptional color rendering index of 95 ensures accurate color representation crucial for surgical precision. The innovative low ceiling height configuration accommodates spaces under 9'6", with additional features such as touch switch handles and backup boards, making it a versatile and essential fixture for modern operating rooms. The purchase also requires a universal anchor plate for installation, with an optional round plastic ceiling cover available upon request, and includes a five-year warranty on the bulbs.

[Read More](#)

SKU: AM-SL-D-30M-30M-LCH



MIDMARK RITTER 255 FLOOR STAND LIGHT

The **Midmark Ritter 255 Floor Stand Light** is designed to deliver superior, white, shadow-free illumination, ensuring safety and coolness to the touch. The innovative yoke design allows for smooth and effortless maneuverability, making it easy to position exactly where needed. With an exceptional color rendering index of greater than 95 and an impressive 7,500 footcandles, the light provides excellent visibility. Its color temperature of 4,400 K ensures clarity and true color representation. The large reflector, featuring 1,040 reflecting facets, guarantees even light distribution and minimal shadows, offering uninterrupted illumination. Additionally, the integrated controls on both sides of the light head allow for convenient adjustments, enhancing user experience and efficiency.

[Read More](#)

SKU: AM-R255-FS



MIDMARK RITTER 355 FLOOR STAND LIGHT

The **Midmark Ritter 355 Single Procedure Light**, a perfect blend of quality and affordability tailored to meet the needs of medical professionals. This light is meticulously counterbalanced, ensuring precise positioning without any drifting, which facilitates accurate procedures. With its 4,200 K color temperature and a large pattern size, the light provides an even distribution that enhances shadow control and delivers superior color rendition, crucial for improved diagnosis. The five-caster base design requires less than two pounds of force to move, allowing for effortless repositioning. Additionally, the handle is fully autoclavable, ensuring high standards of hygiene. The Ritter 355's true color temperature generates less heat, which is beneficial for both patient comfort and prolonged use. Its impressive specifications, including the ability to achieve up to 5,000 foot-candles at a 24-inch focal length, make it an indispensable tool in any medical setting.

[Read More](#)

SKU: AM-R355-FS



SKYTRON ST9 SINGLE FLOOR MOUNT

The Skytron ST9 Exam Light provides high quality, high performance illumination within an extremely durable light weight polymer light head, creating the perfect balance between performance, maneuverability, and durability. The ST9 features advanced high intensity, shadow free optics without glare. Applications include LDR, LDRP, OB/GYN, Exam, Emergency Room & Special Procedures.

[Download Brochure](#)

[Read More](#)

SKU: AM-ST9-SFM



AMICO MIRA LED 65, CEILING MOUNT, DOUBLE, STANDARD ARM

The Amico Mira LED 65, Ceiling Mount, Double, Standard Arm, exemplifies versatile and advanced lighting technology suitable for a variety of critical medical settings. As a true cross-platform design, this lighting system seamlessly integrates into diverse environments such as Emergency Departments (ED), Trauma centers, Intensive Care Units (ICU), Labor and Delivery Rooms (LDR), and spaces designated for Minor Surgery. The dual-arm configuration enhances its functionality, allowing for adjustable and focused illumination critical for precise medical tasks. The use of LED technology not only ensures energy efficiency but also provides consistent, bright light crucial for medical procedures, reducing eye strain and improving the overall visibility for healthcare professionals. This adaptability and high performance make the Mira LED 65 an essential tool in any modern medical facility.

[Read More](#)

SKU: AM-L-MLED65-CM-DC-ST

AMICO MIRA LED 65, CEILING MOUNT, SINGLE, STANDARD ARM



The **Amico Mira LED 65**, with its ceiling mount configuration, exemplifies versatility and efficiency in medical lighting technology. Designed as a true cross-platform solution, the Mira LED Series Lights are ideal for a wide range of medical settings, including Emergency Departments (ED), Trauma centers, Intensive Care Units (ICU), Labor and Delivery Rooms (LDR), and environments suited for minor surgical procedures. This adaptability ensures that healthcare facilities can maintain consistent lighting performance across various departments, enhancing both the functionality and cost-effectiveness of their operations. The Mira LED 65 also offers exceptional illumination quality and ergonomic design, making it a top choice for medical professionals seeking reliable, high-quality lighting solutions.

[Read More](#)

SKU: AM-L-MLED65-CM-SC-ST



AMICO VEGA 80, CEILING MOUNT EXAM PROCEDURE LIGHT

The **Amico Vega 80** offers a powerful and compact procedure examination light for use in various applications throughout the hospital environment. The Vega 80 provides a small environmental footprint and flexible mounting options for ultimate convenience. Its 18 precisely positioned LEDs provide 80,000 lux of high-quality illumination to support caregivers with advanced visual acuity. At the same time, a color temperature of 4500K and 95+ CRI allows for greater precision and color accuracy. Clinicians can easily customize light intensity with five dimming levels and an endo mode functionality located on the light head for easy access.

[Read More](#)

SKU: L-VLED80-CM-ST



STERIS HARMONY AIR M SERIES DUAL HEAD LED SURGICAL LIGHT

Exceptional visualization means better patient outcomes

The HarmonyAIR M-Series Surgical Lighting System provides exceptional surgical site illumination and streamlines suspension affordably without compromising quality.

The M-Series surgical light delivers up to 160,000 lux of pure, white LED light with overlapping beams that virtually eliminate shadows. Sterile surgical light handle controls adjust the focus and brightness intensity to consistently assess and interpret tissue and vessel appearance. The M-Series suspension systems are made to the highest quality standards. Their flexible nature integrates with your workflow and case load, providing you with the ideal solution for your unique challenges.

[Read More](#)

SKU: AM-HARMONYAIR-M-2L



SPECTRA G8 HOSPITAL KIT; G8 HEADLIGHT, 6 BATT, 4-BAY CHARGER, WALL PLUG/USB 45W, PADDING CONV KIT, CARRY CASE

Adjustable, advanced beam technology offering the light intensity to illuminate deep cavities. Versatile, lightweight design and patented AirPad™ comfort padding system in the headband improve the surgeon's experience compared to bulky, heavy headlamps. Hard and Velcro-top options that provide optimal comfort and customization. Dual Everlast™ batteries provide hours of cordless power and can be hot-swapped for infinite use. Integrated lightweight batteries in the head strap eliminate tethering cords and provide a comfortable fit

[Read More](#)

SKU: RXMG8HK



SPECTRA G8 STANDARD KIT; G8 HEADLIGHT, 2 BATT., 45W WALL PLUG/USB 45W, PADDING CONV KIT, CARRY CASE

Adjustable, advanced beam technology offering the light intensity to illuminate deep cavities. Versatile, lightweight design and patented AirPad™ comfort padding system in the headband improve the surgeon's experience compared to bulky, heavy headlamps. Hard and Velcro-top options that provide optimal comfort and customization. Dual Everlast™ batteries provide hours of cordless power and can be hot-swapped for infinite use. Integrated lightweight batteries in the head strap eliminate tethering cords and provide a comfortable fit

[Read More](#)

SKU: RXMG8SK



SPECTRA G6 HOSPITAL KIT; G6 HEADLIGHT, 6 BATT, 4-BAY CHARGER, 45W WALL PLUG/USB, PADDING CONV KIT, CARRY CASE

Adjustable, advanced beam technology offering the light intensity to illuminate deep cavities. Versatile, lightweight design and patented AirPad™ comfort padding system in the headband improve the surgeon's experience compared to bulky, heavy headlamps. Hard and Velcro-top options that provide optimal comfort and customization. Dual Everlast™ batteries provide hours of cordless power and can be hot-swapped for infinite use. Integrated lightweight batteries in the head strap eliminate tethering cords and provide a comfortable fit

[Read More](#)

SKU: RXMG6HK



SPECTRA G6 STANDARD KIT; G6 HEADLIGHT, 2 BATT., 45W WALL PLUG/USB 45W, PADDING CONV KIT, CARRY CASE

Adjustable, advanced beam technology offering the light intensity to illuminate deep cavities. Versatile, lightweight design and patented AirPad™ comfort padding system in the headband improve the surgeon's experience compared to bulky, heavy headlamps. Hard and Velcro-top options that provide optimal comfort and customization. Dual Everlast™ batteries provide hours of cordless power and can be hot-swapped for infinite use. Integrated lightweight batteries in the head strap eliminate tethering cords and provide a comfortable fit

[Read More](#)

SKU: RXMG6SK



SPECTRA G3 HOSPITAL KIT; G3 HEADLIGHT, 6 BATT, 4-BAY CHARGER, 45W WALL PLUG/USB, PADDING CONV KIT, CARRY CASE

Adjustable, advanced beam technology offering the light intensity to illuminate deep cavities. Versatile, lightweight design and patented AirPad™ comfort padding system in the headband improve the surgeon's experience compared to bulky, heavy headlamps. Hard and Velcro-top options that provide optimal comfort and customization. Dual Everlast™ batteries provide hours of cordless power and can be hot-swapped for infinite use. Integrated lightweight batteries in the head strap eliminate tethering cords and provide a comfortable fit

[Read More](#)

SKU: RXMG3HK



SPECTRA G3 STANDARD KIT; G3 HEADLIGHT, 2 BATT., 45W WALL PLUG/USB, PADDING CONVERSION KIT, CARRYING CASE

Adjustable, advanced beam technology offering the light intensity to illuminate deep cavities. Versatile, lightweight design and patented AirPad™ comfort padding system in the headband improve the surgeon's experience compared to bulky, heavy headlamps. Hard and Velcro-top options that provide optimal comfort and customization Dual Everlast™ batteries provide hours of cordless power and can be hot-swapped for infinite use. Integrated lightweight batteries in the head strap eliminate tethering cords and provide a comfortable fit.

[Read More](#)

SKU: RXMG3SK



INTEGRA LUXTEC MLX 300W LIGHT SOURCE

The Integra Luxtec MLX 300 is a xenon light source that delivers 300 watts of cool white infrared (IR) filtered light. The lamp is housed in a lamp module that can be easily serviced without special tools. The Luxtec MLX 300 is designed to supply high-intensity white light through a fiber optic cable for illumination of a surgical field during a surgical and/or medical procedure.

[Read More](#)

SKU: AM-LUXTECMLX



LUXTEC INTEGRA LED HEADLAMP PACKAGE

Comes with the LED Headlight, Adjustable headband, 2 Lithium Ion Batteries, Battery Holster, Battery Charger, Operating Manual. Super light, super bright with a 25,000 hour lamp life and 4 intensity settings. The beam is variable spot and there is an integrated negative flow cooling fan in the headset. Completely portable with the two included batteries

[Read More](#)

SKU: AM-LINTEGRA-LEDH



SKYTRON AURORA LED 2 HEAD SURGICAL LIGHT SYSTEM

Skytron's Aurora is a LED surgical lighting system equipped with high-performance hybrid LED optics that provides surgeons with direct control of intensity and focus.

[Read More](#)

SKU: AM-SKY-AURORA-2L



SKYTRON AURORA LED SINGLE SURGICAL LIGHT SYSTEM

Skytron's Aurora is a LED surgical lighting system equipped with high-performance hybrid LED optics that provides surgeons with direct control of intensity and focus.

[Read More](#)

SKU: AM-SKY-AURORA-1L



STERIS LED385 EXAM LIGHT ON FLOOR STAND

The **Steris LED385 Exam Light on a Floor Stand** is a highly efficient medical examination light designed to provide exceptional visibility and color accuracy for healthcare professionals. With a light intensity that reaches up to 50,000 Lux, it ensures that examination areas are well-illuminated. Its 6-inch fixed spot pattern size is ideal for focused assessments, allowing precise illumination without the need to frequently adjust the light's position. The LED385 also boasts a color temperature of 5000° K, simulating natural daylight to aid in accurate color discernment of tissues and surfaces. Its impressive Color Rendering Index (CRI) of 90 further enhances the ability to accurately recognize color nuances, which is crucial during medical examinations. Additionally, the light maintains a low heat-to-light ratio of 3.4 mW/m²-lx, minimizing patient discomfort from heat exposure. Durability is also a key feature, with an average LED bulb life of 30,000 hours, making it a reliable choice for busy clinical environments. This combination of features makes the Steris LED385 Exam Light on Floor Stand an optimal solution for medical settings requiring high-quality, focused lighting.

[Read More](#)

SKU: AM-LED385




AMICO VEGA 80, FLOOR STAND EXAM PROCEDURE LIGHT

The **Amico Vega 80** offers a powerful and compact procedure examination light for use in various applications throughout the hospital environment. The Vega 80 provides a small environmental footprint and flexible mounting options for ultimate convenience. Its 18 precisely positioned LEDs provide 80,000 lux of high-quality illumination to support caregivers with advanced visual acuity. At the same time, a color temperature of 4500K and 95+ CRI allows for greater precision and color accuracy. Clinicians can easily customize light intensity with five dimming levels and an endo mode functionality located on the light head for easy access.

[Read More](#)

SKU: L-VLED80-FS-XX-ST

Nova Motorized LED Exam Light



AMICO NOVA MOTORIZED EXAM PAIR 4000K, SINGLE KEYPAD

- 360 Degrees of rotation provides targeted lighting to the patient bed for increased functionality
- Compact, clutter free design provides a visual comfort to the patient
- Fully recessed exam light reduces potential hazards in low-ceiling environments and provides more space for critical healthcare equipment
- 55,000 Lux of light at task delivers high level illumination to support clinicians with advanced visual acuity
- Multiple ceiling configurations increase functionality and ensures maximum light output onto the patient bed
- Hospital grade material and finish
- Low voltage wall control keypad fits into a single gang box
- Fully recessed installation
- Long lasting LEDs with various color temperature options
- Low voltage wall control keypad fits into a single gang box
- One keypad to control up to 2 Nova exam lights

[Read More](#)

SKU: L-NVA-P-277-LED40-1-00



AMICO MIRA LED 90, FLOOR STAND, SINGLE

Being a true cross-platform design, the Amico Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

[Read More](#)

SKU: AM-L-MLED90-FS-XX-ST



AMICO ICE LED SURGICAL LIGHT -TRIPLE CEILING, 30M/30M/30M

iCE 30M LED SPECIFICATIONS:

Light-head Diameter - 30in. (76cm)

Light Intensity (Central Luminance) - 160,000 lux (14,870 fc)

Brightness Control -5 - 100%

Light Field Diameter - 7.6 - 12in. (19 - 31cm)

Depth of Illumination (L1 +L2) at 60% - 27.25 - 48in.

Depth of Illumination (L1 +L2) at 20% -42in. (107cm)

Color Rendering Index (CRI) -95

R9 >90

R13 - 97

[Read More](#)

SKU: AM-SL-T-30M-30M



**No Image
Available**



**No Image
Available**

AMICO SINGLE WALL CONTROL PANEL

The Wall Control Panel will have the exact same features as the Lighthouse Control Panel.

1. Turning the light ON and OFF: The LED will light up when the lamp is in the ON position.
2. Adjusting the light brightness: The brightness can be adjusted to 50% - 100% (total of 5 steps). Brightness level indicators show the current brightness level of the LED.
3. Switching to Endo mode: To deactivate Endo mode you can either press the Endo button or press any brightness buttons.
4. The color can be adjusted from warmer to cooler from 3500 K - 5000 K. Use the + and - keys to adjust color.

[Read More](#)

SKU: AM-SL-WC-S



AMICO ICE LED SURGICAL LIGHT - DUAL CEILING, 30M/SINGLE MONITOR HOLDER

Amico iCE LED Surgical Light - Dual Ceiling, 30M/Single Monitor Holder

[Read More](#)

SKU: AM-SL-D-30M-MH1



AMICO ICE LED SURGICAL LIGHT - SINGLE CEILING, 30M

The iCE LED surgical lighting system was designed to illuminate with a bright, white light and an adjustable pattern size that creates a cool and comfortable surgical field. The Amico ceiling mount plate allows for easy installation for new or upgraded lighting.

The LEDs in the iCE surgical lights produce a bright white light that stays cool and lasts a minimum of 50,000 hours. The surfaces of the iCE LED lights are smooth making them convenient and easy to clean.

Each light-head comes standard with sterile, disposable covers for the non-sterilizable pattern of light beam adjustment handle, sterile, disposable covers for the in-light HD camera and sterilizable handles, which allow the user to autoclave the handle before each use.

[Read More](#)

SKU: AM-SL-S-30M



AMICO MIRA LED 90, CEILING MOUNT, DOUBLE, STANDARD ARM

Being a true cross-platform design, the Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

550mm reach for standard arms

[Read More](#)

SKU: AM-L-MLED90-CM-DC-ST

AMICO MIRA LED 90, CEILING MOUNT, SINGLE, STANDARD ARM



Being a true cross-platform design, the Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

550mm reach for standard arm

[Read More](#)

SKU: AM-L-MLED90-CM-SC-ST

AMICO MIRA LED 90, CEILING MOUNT, SINGLE, EXTENDED ARM



Being a true cross-platform design, the Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

[Read More](#)

SKU: AM-L-MLED90-CM-SC-EX

AMICO GAMMA LED 60, FLOOR STAND, SINGLE



Powerful LED's provide 60,000 Lux in this new innovative exam light. The Gamma 60 LED is maintenance free and suitable for any type of exam room. The light comes with four dimming levels, two color shades, an antimicrobial finish and a 360° rotary lighthead. The Gamma 60 LED can be ceiling mounted, rail mounted or on a roll stand to meet all your room requirements.

[Read More](#)

SKU: AM-L-GLED60-AB-FX-X



AMICO GAMMA LED 50, FLOOR STAND, SINGLE

LEDs provide 50,000 lux of light for this new, innovative exam light. The Gamma 50 LED is maintenance-free and suitable for any type of exam room. There are five dimming levels with a simple switch on the head of the fixture. The articulating arm system offers simple positioning and movement of the light. Also, heat from the light is dissipated to the back of the head and away from the user.

[Read More](#)

SKU: AM-L-GA50-AR-FS-00-00



AMICO MIRA LED 50, FLOOR STAND, SINGLE

Being a true cross-platform design, the Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

550mm reach for standard arm

[Read More](#)

SKU: AM-L-MLED50-FS-XX-ST

AMICO MIRA LED 50, CEILING MOUNT, SINGLE, EXTENDED ARM



Being a true cross-platform design, the Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

Extended arm reaches 720 mm over 550mm for standard arm

[Read More](#)

SKU: AM-L-MLED50-CM-SC-EX

AMICO MIRA LED 65, CEILING MOUNT, SINGLE, EXTENDED ARM



Being a true cross-platform design, the Mira LED Series Lights can be used in your ED, Trauma, ICU, LDR and Minor Surgery environments.

Extended arm reaches 720 mm over 550mm for standard arm

[Read More](#)

SKU: AM-L-MLED65-CM-SC-EX

NUVO/MEDICAL ILLUMINATION SHUTTLE TRACK DUAL TROLLEY



69" of Travel
Can mount two Ceiling Lights to Track. Must be same models
Does not include Light Pictured -
Supports Nuvo EX, PR, PRO exam and procedure lights.

[Read More](#)

SKU: AM-NV064012

NUVO/MEDICAL ILLUMINATION SHUTTLE TRACK SINGLE TROLLEY



69" of Travel
Includes: Drop tube and Trolley included

[Read More](#)

SKU: AM-NV064011



STERIS HARMONY AIR G SERIES DUAL HEAD LED SURGICAL LIGHT

Exceptional visualization means better patient outcomes

The HarmonyAIR G-Series Surgical Lighting System provides exceptional surgical site illumination and streamlines suspension affordably without compromising quality.

The G-Series surgical light delivers up to 160,000 lux of pure, white LED light with overlapping beams that virtually eliminate shadows. Sterile surgical light handle controls adjust the focus and brightness intensity to consistently assess and interpret tissue and vessel appearance. The G-Series suspension systems are made to the highest quality standards. Their flexible nature integrates with your workflow and case load, providing you with the ideal solution for your unique challenges.

[Read More](#)

SKU: AM-HARMONYAIR-G-2L



AMICO ICE LED SURGICAL LIGHT - DUAL CEILING, 30M/30M

The iCE LED Surgical Lighting System comes in two light-head options - iCE 30M LED, with 30" (762mm) diameter, 80 LEDs, 160,000 lux, 7.6"-12" (193-305 mm) light pattern and laminar flow opening or iCE 25LED, with 25" (635mm) diameter, 52 LEDs, 160,000 lux or 130,000 lux and 7.6"-11" (193-279mm) light pattern. The iCE LED surgical lighting system was designed to illuminate with a bright, white light and an adjustable pattern size that creates a cool and comfortable surgical field. The Amico ceiling mount plate allows for easy installation for new or upgraded lighting.

[Read More](#)

SKU: AM-SL-D-30M-30M



AMICO ICE LED SURGICAL LIGHT

The iCE LED Surgical Lighting System comes in two light-head options - iCE 30M LED, with 30" (762mm) diameter, 80 LEDs, 160,000 lux, 7.6"-12" (193-305 mm) light pattern and laminar flow opening or iCE 25LED, with 25" (635mm) diameter, 52 LEDs, 160,000 lux or 130,000 lux and 7.6"-11" (193-279mm) light pattern. The iCE LED surgical lighting system was designed to illuminate with a bright, white light and an adjustable pattern size that creates a cool and comfortable surgical field.

Installation Plate required when purchasing.

[Read More](#)

SKU: AM-SL-T-30M-30M-MH1

Performance	Safety	Workflow	Patient Experience	Interoperability	Cybersecurity	Ma
Good	Good	Excellent	Not evaluated	Not evaluated	Not evaluated	Go
Fair	Good	Excellent	Not evaluated	Not evaluated	Not evaluated	Go
Excellent	Good	Good	Not evaluated	Not evaluated	Not evaluated	Go
Good	Good	Good	Not evaluated	Not evaluated	Not evaluated	Go
Fair	Good	Excellent	Not evaluated	Not evaluated	Not evaluated	Go
Good	Good	Good	Not evaluated	Not evaluated	Not evaluated	Go
Fair	Good	Excellent	Not evaluated	Not evaluated	Not evaluated	Go



NUVO VU TRIPLE HEAD SURGICAL W/ DUAL LIGHT, ONE CAMERA, AND ONE MONITOR ARM

The **Nuvo VU Triple Head Surgical system, equipped with Dual Light, One Camera, and One Monitor Arm**, represents a significant advancement in surgical lighting technology. Featuring innovative Cross Focus Technology (CFT), this system allows the LED lights to operate independently, enabling precise control over the illumination pattern size, which can be adjusted from 9.5 inches to 14 inches (229 mm to 356 mm). This flexibility ensures consistent light output across both large and small surgical fields. The system boasts an impressive intensity of 160,000 lux and a color temperature of 4300°K, providing natural-like light that enhances visibility. With a high Color Rendering Index (CRI) and R9 Value, both at 95, the lighting system offers exceptional color fidelity, crucial for accurate tissue differentiation during procedures. Additional features include ambient light mode, sterile handle controls for easy adjustments, and a long LED life expectancy of 50,000 hours. The inclusion of a monitor arm enhances the system's functionality, though the monitor itself is sold separately. This high-specification setup is designed to meet the demanding needs of modern surgical environments, ensuring optimal illumination and improved operational efficiency.

[Read More](#)

SKU: AM-NVSD3D3M1



NUVO VU LED SURGICAL LIGHT

When you're performing surgery, it's critical that you can see what you're doing down to the most minute detail. That's why the Nuvo Vu LED Surgical Light should be at the top of the list for your hospital operating rooms. The Nuvo Vu LED Surgical Light is brighter, whiter, cooler and more capable in the operating room than any previous technology. So you get the most accurate colors with unparalleled shadow control and virtually no heat emission. All this for a value that's hard to match.

[Read More](#)

SKU: AM-NVSD3

NUVO VU OR LIGHTING TOUCH SCREEN TABLET CONTROL (OPTIONAL)



Touchscreen wall control enhances efficiency in the OR

Nuvo gives you the option of a touchscreen wall control. The wall control panel communicates wirelessly, managing one, two, or three lights as well as the HD Video Camera. While you focus on operating, your nurse or assistant can control light intensity, color temperature, Green mode, plus camera functions — all remotely.

[Read More](#)

SKU: AM-NV1000775

Price: \$2,511.50



NUVO VISTOR PRO LED SINGLE HEAD CEILING PROCEDURE LIGHT

Adding unsurpassed brilliance in the OR
Ideal within the surgical environment, the VistOR Pro offers unsurpassed longevity and flexibility in a package that's incredibly affordable. The large diameter LED array provides excellent shadow control - even with just a single light head. Pressing the center button on the sterilizable handle allows 5-stage dimming plus ambient light mode (for minimally-invasive procedures) for easy management of brightness within the illumination field.

[Read More](#)

SKU: AM-NV061524

NUVO VISTOR EX LED EXAMINATION LIGHT WALL MOUNT



Significantly brighter, whiter and cooler than any previous technology—at the greatest value.

Large flood style pattern delivering a 6-inch pattern of light.

The light can be paired with the VistOR PR and VistOR PRO models to provide a combination of lighting.

[Read More](#)

SKU: AM-NV61322

NUVO VISTOR EX PORTABLE LED EXAMINATION LIGHT FLOOR MOUNT



Significantly brighter, whiter and cooler than any previous technology—at the greatest value.

Large flood style pattern delivering a 6-inch pattern of light.

The light can be paired with the Vistor PR and Vistor PRO models to provide a combination of lighting.

[Read More](#)

SKU: AM-NV61323

STERIS VLED SINGLE HEAD CEILING MOUNTED SURGICAL LIGHT



The **Steris Harmony vLED Surgical Lighting System** stands out as a highly versatile and technologically advanced solution for various medical environments. Designed to cater to a broad range of surgical and diagnostic needs, this modular and configurable system offers unparalleled illumination quality, essential for precision in procedures. With a maximum central illuminance of 160,000 lux and an adjustable pattern size, it ensures optimal visibility in both small outpatient surgery centers and large acute care suites. The system's color temperature of 4400K and an impressive Color Rendering Index (CRI) of 96 enable medical professionals to discern subtle variations in tissue color and texture, critical during intricate surgeries. Moreover, its low heat-to-light ratio minimizes the thermal discomfort often experienced in intensely lit settings, enhancing the comfort of both the surgical team and the patient. With a long lifespan of 30,000 hours, the Harmony vLED system is not only effective but also an economical lighting solution for operating rooms, emergency departments, and intensive care units.

[Read More](#)

SKU: AM-SVLED-SCM

STERIS HARMONY LC 500 OR LIGHTS



The cost-effective Harmony® LC Surgical Lighting System provides bright halogen lighting and a flexible, easy-to-maneuver suspension system. The Harmony® LC Surgical Lighting system is a contemporary product recently emerging on the aftermarket and is offered refurbished on a limited basis.

[Read More](#)

SKU: AM-LC500



RITTER MIDMARK 354 LIGHT

The Ritter 354 Procedure Spotlight has everything healthcare professionals need - but without the heat. It's an ideal task light for all medical environments where illumination is required for external examinations and procedures. And with its pure, crisp and defined light pattern, it can be adjusted to suit individual needs.

[Download Brochure](#)

[Read More](#)

SKU: AM-R354-SF



SKYTRON ST9 SINGLE CEILING MOUNT

The Skytron ST9 Exam Light provides high quality, high performance illumination within an extremely durable light weight polymer light head, creating the perfect balance between performance, maneuverability, and durability. The ST9 features advanced high intensity, shadow free optics without glare. Applications include LDR, LDRP, OB/GYN, Exam, Emergency Room & Special Procedures.

[Download Brochure](#)

[Read More](#)

SKU: AM-ST9-SCM



RITTER MIDMARK 355 PROCEDURE LIGHT

The Ritter 355 Procedure Light is the perfect balance between quality and affordability. It's counterbalanced to provide precise positioning without drifting. A 4,200 K color temperature, large pattern size with an even distribution of light, and unsurpassed shadow control provides the most accurate color rendition for improved diagnosis.

[Download Brochure](#)

[Read More](#)

SKU: AM-R355-SM



DUAL VISTOR PRO LED CEILING MOUNT

The VistOR PRO Procedural Light features a large 16-inch, highly-polished faceted reflector which keeps shadowing to a minimum and maximizes the light output. This pre-focused light also combines a great depth of field, smooth, drift-free positioning and an easy-to-change bulb that doesn't require the use of tools.

[Download Brochure](#)

[Read More](#)

SKU: AM-NV061525

NUVO VISTOR MSLED DUAL ARM LIGHTING SYSTEM



The NUVO Surgical Light provides superior performance for all aspects of surgery. The light head provides excellent shadow reduction, a high-intensity light pattern and optimal color rendition. Its unique, patented interface allows the user to fully control the light's intensity, pattern, position and orientation from the light head. Appropriate for the vast majority of general and specialized surgical procedures including Micro Surgery, Orthopedics, OB/GYN, Ophthalmology, and Cosmetic Surgery.

[Download Brochure](#)

[Read More](#)

SKU: AM-NVMSLED-DLM

STERIS SQ240



The Steris® SQ240™ Surgical Lighting Fixture offers a two-arm unit mounted on a central hub. The rotary suspension shaft of the surgical light head moves on tapered roller bearings at one end of the horizontal arm. The unit produces up to 10,000-foot candles and provides continuous lighting by automatically switching to an auxiliary bulb when needed. The SQ240™ is also available in Single Mount Lights.

[Technical Specs](#)

[Download Brochure](#)

[Read More](#)

SKU: AM-SQ240-DCM



STERIS HARMONY SERIES

Steris® offers a wide variety of high-quality surgical lighting products to meet the needs of any operating room - from large teaching hospitals to smaller Ambulatory Surgery Centers - because every patient deserves the best possible light. Harmony lights maintain the truest color for quick assessment of tissues and vessels to help surgeons provide the best possible patient outcomes. The optical coatings filter out 98% of heat-carrying infrared rays, which creates a cooler, more comfortable surgical environment.

[Download Steris Harmony 585 Brochure](#)

[Download Steris Harmony vLED Brochure](#)

[Download Steris Harmony LC Brochure](#)

[Download Steris Harmony LA500 Brochure](#)

[Download Steris Harmony LA300 Brochure](#)

[Read More](#)

SKU: N/A

SKYTRON STELLAR SERIES



Stellar® provides advanced lighting performance with high intensity, shadow-free and focusable optics, offering precise surgeon focus control for every procedure while retaining maximum intensity. Extended working distances, beyond one meter, are not uncommon for many surgical procedures such as OB-GYN, Orthopedics, Vascular, and Urology. Stellar's advanced VSRD Optics provide surgeons with sterile focus control. Just a 1/4 turn of the sterile center focus handle captures maximum intensity at any distance without having to physically reposition lightheads.

[Download Brochure](#)

[Read More](#)

SKU: N/A

SKYTRON AURORA SERIES



Skytron's Aurora® LED Surgical Lighting provides advanced bright, cool and color correct illumination to the surgical field and a selectable color temperature of 4,000K (soft white) or 4,500 (bright white) based upon surgeon preference for every procedure. Focusable LEDs extend illumination range while retaining maximum intensity at any distance. Typical applications include General, Orthopedic, Plastics, Vascular, Cardiovascular, OB/GYN, Trauma, Neurosurgery, and Urology.

[Download Brochure](#)

[Read More](#)

SKU: N/A

INDEX

A

- Amico Clarity 50 Exam Light, Roll Stand 3
- Amico Gamma LED 50, Floor Stand, Single 23
- Amico Gamma LED 60, Floor Stand, Single 22
- Amico iCE LED Surgical Light 27
- Amico iCE LED Surgical Light - Dual Ceiling, 30M/30M 26
- Amico iCE LED Surgical Light - Dual Ceiling, 30M/30M Low Ceiling Height Model 7
- Amico iCE LED Surgical Light - Dual Ceiling, 30M/Single Monitor Holder 20
- Amico iCE LED Surgical Light - Single Ceiling, 30M 20
- Amico iCE LED Surgical Light -Triple Ceiling, 30M/30M/30M 18
- Amico ICE Series Surgical Light Disposable Handle Cover for Non-Sterilizable Handle (1 box - 2pcs/pkg, 80pkg/box) 2
- Amico Leo Series MG 6in Clinical Recessed Downlight 5
- Amico Mira Disposable Handle Covers Box of 50 2
- Amico Mira LED 50, Ceiling Mount, Single, Extended Arm 24
- Amico Mira LED 50, Floor Stand, Single 23
- Amico Mira LED 65, Ceiling Mount, Double, Standard Arm 9
- Amico Mira LED 65, Ceiling Mount, Single, Extended Arm 24
- Amico Mira LED 65, Ceiling Mount, Single, Standard Arm 10
- Amico Mira LED 90, Ceiling Mount, Double, Standard Arm 21
- Amico Mira LED 90, Ceiling Mount, Single, Extended Arm 22
- Amico Mira LED 90, Ceiling Mount, Single, Standard Arm 21
- Amico Mira LED 90, Floor Stand, Single 18
- Amico Nova Motorized Exam Pair 4000K, Single Keypad 17
- Amico Nuvo Vu Dual Head LED Surgical Light 5
- Amico Single Wall Control Panel 19
- Amico Vega 80, Ceiling Mount Exam Procedure Light 10
- Amico Vega 80, Floor Stand Exam Procedure Light 17

D

- Dual Vistor PRO LED Ceiling Mount 33

I

- Integra Luxtec MLX 300W Light Source 14

L

- Luxtec 90600 Integra DUO LED Headlamp package 4
- Luxtec Integra LED Headlamp Package 15

M

- Midmark Ritter 255 Floor Stand Light 7
- Midmark Ritter 355 Floor Stand Light 8

N

- Nuvo Vistor EX LED Examination Light Wall Mount 29
- Nuvo Vistor EX Portable LED Examination Light Floor Mount 30
- NUVO Vistor MSLED Dual Arm Lighting System 33
- NUVO Vistor PRO LED Single Head Ceiling Procedure Light 29

NUVO Vu LED Surgical Light 28
Nuvo Vu OR Lighting Touch Screen Tablet Control (Optional) 28
Nuvo VU Triple Head Surgical w/ Dual Light, One Camera, and One Monitor Arm 27
Nuvo/DRE/Medical Illumination Surgical and Exam Light Disposable Sterile Handle Covers Case of 25 3
Nuvo/Medical Illumination Shuttle Track Dual Trolley 25
Nuvo/Medical Illumination Shuttle Track Single Trolley 25

R

Ritter Midmark 354 Light 31
Ritter Midmark 355 Procedure Light 32

S

Skytron Aurora LED 2 Head Surgical Light System 15
Skytron Aurora LED Single Surgical Light System 16
Skytron Aurora Series 35
Skytron ST9 Single Ceiling Mount 32
Skytron ST9 Single Floor Mount 9
Skytron Stellar Series 35
Spectra G3 Hospital Kit; G3 Headlight, 6 Batt., 4-Bay Charger, 45W Wall Plug/USB, Padding Conv Kit, Carry Case 13
Spectra G3 Standard Kit; G3 Headlight, 2 Batt., 45W Wall Plug/USB, Padding Conversion Kit, Carrying Case 14
Spectra G6 Hospital Kit; G6 Headlight, 6 Batt, 4-Bay Charger, 45W Wall Plug/USB, Padding Conv Kit, Carry Case 12
Spectra G6 Standard Kit; G6 Headlight, 2 Batt., 45W Wall Plug/USB 45W, Padding Conv Kit, Carry Case 13
Spectra G8 Hospital Kit; G8 Headlight, 6 Batt, 4-Bay Charger, Wall Plug/USB 45W, Padding Conv Kit, Carry Case 11
Spectra G8 Standard Kit; G8 Headlight, 2 Batt., 45W Wall Plug/USB 45W, Padding Conv Kit, Carry Case 12
Steris Dual Center Mount vLED Surgical Light 6
Steris Harmony Air G Series Dual Head LED Surgical Light 26
Steris Harmony Air M Series Dual Head LED Surgical Light 11
Steris Harmony LA 585 LED Dual Head OR Lights 6
Steris Harmony LA 585 LED Single Head OR Light 4
Steris Harmony LC 500 OR Lights 31
Steris Harmony Series 34
Steris LED385 Exam Light on Floor Stand 16
Steris SQ240 34
Steris VLED Single Head Ceiling Mounted Surgical Light 30