

# TABLE OF CONTENTS

**Surgical Power Equipment** ..... 2

    Stryker 940 Cast Cutter - Stainless Steel Blade (#STB043) ..... 3

    Stryker 940 Cast Cutter - Titanium Nitrided Blade (#STB045) ..... 3

    Stryker System 7 7215 Battery ..... 3

    Fully Refurbished Stryker Cordless Driver 3 and attachments ..... 4



## SURGICAL POWER EQUIPMENT

We know that precision, reliability, and familiarity are important when considering your surgical power equipment. All of our refurbished drivers, drills, and support equipment are rebuilt and certified through a rigorous quality assurance process to assure our customers that they are receiving the highest quality of their favorite brands.

If you do not see what you're looking for in this category, we may still carry it in our 25,000 SKU catalogue. Complete our contact form at the footer and a sales representative will reach out to you today!

**Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:**

**[BRAND NEW](#)**

**[RECONDITIONED](#)**

**[REFURBISHED](#)**



## STRYKER 940 CAST CUTTER - STAINLESS STEEL BLADE (#STB043)

The **Stryker 940 Cast Cutter Stainless Steel Blade (Replacement #STB043, OEM #940-25)** is designed for use with the Stryker 940 cast saw. This blade is specifically **optimized for plaster casts**, offering dependable performance for that material. It features a 2.36-inch oval blade and is made from stainless steel, providing solid durability and reliable cutting for plaster applications. Each replacement package includes 1 blade.

**QTY: 1 Blade**

[Read More](#)

**SKU:** 0940-025-000

**Price:** \$114.95



## STRYKER 940 CAST CUTTER - TITANIUM NITRIDED BLADE (#STB045)

The **Stryker 940 Cast Cutter Titanium Nitrided Blade (Replacement #STB045, OEM# 940-035)** is a precision blade designed for cutting synthetic casts. The STB045 blade fits the Stryker 940 cast saw and features a durable titanium-nitrided coating for enhanced longevity and smooth performance. It has a 2.36-inch oval blade and is **specifically optimized for synthetic cast materials**, delivering reliable cuts with improved strength and corrosion resistance.

**QTY: 1 Blade**

[Read More](#)

**SKU:** 0940-035-000

**Price:** \$114.95



**No Image  
Available**

## STRYKER SYSTEM 7 7215 BATTERY

Stryker System 7 7215 Battery  
Will not fit System 6 handpieces

[Read More](#)

**SKU:** AM-STRY-7215

**Price:** Request Quote



## FULLY REFURBISHED STRYKER CORDLESS DRIVER 3 AND ATTACHMENTS

Fully Refurbished Stryker Cordless Driver 3 Cordless Driver 3 extends the time-honored tradition of power and versatility in the Sports Medicine, Trauma and Orthopaedics arenas. It's compact and ergonomically designed with a smaller, lighter handpiece for optimal control and comfort.

[Read More](#)

**SKU:** AM-4300-DR3

**Price:** Request Quote

# INDEX

## F

Fully Refurbished Stryker Cordless Driver 3 and attachments 3

## S

Stryker 940 Cast Cutter - Stainless Steel Blade (#STB043) 2

Stryker 940 Cast Cutter - Titanium Nitrided Blade (#STB045) 2

Stryker System 7 7215 Battery 3