

TABLE OF CONTENTS

Surgical Tables	4
STANDARD SERIES - SURGICAL TABLE PAD SETS - FOR: 3080RC - 3080RL - 3080SP - 3085SP	5
Mindray V6 Surgical Table Extension	5
Skytron 6600B Surgical Table	5
Biodex 058-878 30in Wide Top Imaging Table	6
Morgan MEDesign FLTLT-FW (Float-Top, Free-Wheel Configuration)	7
STI EconoMAX 4	7
Steris 5085 Surgical Table	8
Patient Transfer Roller Board - 67" x 15"	8
Surgical Tables Inc. SpineXR Table	9
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide and Carbon Fiber Extension (Stainless Steel Base)	10
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide (Stainless Steel Base)	11
Mindray HyBase V6 Surgical Table (Stainless Steel Base)	11
Mindray HyBase V6 Surgical Table (Plastic Base)	12
Mayfield Headrest - Basic Headrest	13
Steris Narrow Headrest Lightweight, Radiolucent Head Support	14
Steris 9 inch Urology Extension	14
Mizuho OSI 6875 Hana Table	15
Steris SurgiGraphic 1027 C-arm Imaging Table	16
OSI Hana Equipment Cart	16
Skytron 3603 Ultra Slide Surgical Table	17
Steris Amsco Model 3085 Surgical Table	17
Biodex 840 Floating Top Imaging Table	17
Table Protective Drape CFUR Protective Drape- Clear PVC- attaches to T-Rail	18
Table Protective Drape CFLU Protective Drape- Clear PVC- attaches to T-Rail	18
CFPM300 Table 3 motion- Rectangular Top- 2in Pad, Hand Control	19
CFUR401 Urology Table 4 motion- Carbon Fiber Top- 2in Pad, Hand Control & Foot Control	19
CFUR301 Urology Table 3 motion- Carbon Fiber Top- 2in Pad, Hand Control & Foot Control	20

CFPMB301 Table 3 motion- 750 lbs rated, Rectangular Top- 2in Pad, Hand Control	20
CFPM401 Table 4 motion- Rectangular Top- 2in Pad, Hand Control	21
CFPM401 Table 4 motion- Integrated Headrest Top- 2in Pad, Hand Control	21
CFPM400 Table 4 motion- Integrated Headrest Top- 2in Pad, Hand Control	22
CFPM302 Table 3-motion- Rectangular Top- 2in Pad, Hand Control	22
CFPM302 Table 3-motion- Integrated Headrest Top- 2in Pad, Hand Control	23
CFPM301 Table 3 motion- Rectangular Top- 2in Pad Hand Control	23
CFPM301 Table 3 motion- Integrated Headrest Top- 2in Pad, Hand Control	24
CFPM300 Table 3 motion- Rectangular Top- 2in Pad, Hand Control	24
CFPM300 Table 3 motion- Integrated Headrest Top- 2in Pad, Hand Control	25
CFPM100 Table 1 motion- Rectangular Top- 2in Pad Foot Control	25
CFPM100 Table 1 motion- Integrated Headrest Top- 2in Pad, Foot Control	26
CFLU401 Lithotripsy/Urology Table 4 motion- Carbon Fiber Top- 2in Pad Hand Control & Foot Control	26
IDI 100RT Table	27
Skytron 6600 Surgical Table	27
Steris CMAX Surgical Table	28
Skytron 6702 Surgical Table	28
Berchtold OPERON OR Table D820	29
Berchtold OPERON OR Table D850	30
Maquet Alphamaxx 1133 Surgical Table	30
Nuvo Surgical V1000 Surgical Table Series	30
Skytron 3600B Surgical Table	31
Surgical Tables Inc. V-Max 500 lbs capacity (No Trend)	31
STI Uromax Series	32
Streamline 5 Pain Table	32
STI EconoMAX 3	33
STI Streamline 3 Pain Table	33
Streamline 2 Pain Table	34
Streamline 1 Pain Table	34
Streamline Pain Table, Fixed Height	35
Wilson Spinal Frame Radiolucent	35
NUVO Surgical V800	36

Steris 4085	36
Skytron 3501B	38
Skytron 3500 Elite	38
Skytron 6700	39
Steris Amsco 3080 Surgical Table	39
Skytron 6500	40
Skytron Hercules 6701 Surgical Table	40



SURGICAL TABLES


A superior surgical table is the foundation for better patient outcomes. Auxo Medical wide range of surgical tables, exam tables, and orthopedic tables provide a solution for every operating room and budget. Furthermore, on most typical tables, the head area is removable and may have an assortment of attachable head rests. In addition, our selection features operating tables ideal for use in hospital. Whether you're searching for flexibility, optimal imaging, or affordability, we carry the technology you need to improve safety.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



STANDARD SERIES - SURGICAL TABLE PAD SETS - FOR: 3080RC - 3080RL - 3080SP - 3085SP

The **Standard Series Surgical Table Pad Sets** are specially designed accessories compatible with **Steris Surgical Table models 3080RC, 3080RL, 3080SP, and 3085SP**. These pad sets provide enhanced comfort and support for patients during surgical procedures, ensuring optimal positioning and safety. Crafted from high-quality, durable materials, they are easy to install and maintain, offering reliable performance to meet the demands of various surgical environments. Whether for general surgery or specialized procedures, these pad sets are essential accessories that help improve patient care and procedural efficiency.

[Read More](#)

SKU: AM-SC-3320

Price: \$495.00



Mindray V6 Surgical Table Extension

MINDRAY V6 SURGICAL TABLE EXTENSION

The **Mindray V6 Surgical Table Extension** is a versatile accessory designed to expand the functionality and coverage of the Mindray V6 surgical table. Engineered to enhance surgical flexibility, this extension provides additional surface area, allowing for better patient positioning and support during complex procedures. Its robust construction ensures stability and durability, making it an essential component for surgical environments that require precise patient alignment and optimal access. Compatible with the Mindray V6 surgical table, this extension helps healthcare providers deliver more effective and efficient surgical care by accommodating a wider range of patient sizes and procedural needs.

[Read More](#)

SKU: 115-103847-00

Price: Request Quote




SKYTRON 6600B SURGICAL TABLE

The **Skytron 6600B Surgical Table** is a versatile, high-capacity operating table designed to support demanding surgical procedures with exceptional imaging access. Capable of lifting up to 1,000 lbs and articulating up to 600 lbs, it features removable Back and Leg Sections that accommodate Carbon Fiber components for enhanced imaging windows. Its 210° rotating top improves access for C-arm and microscope use, while the stable, rigid surface ensures patient positioning remains secure throughout surgery. The table can be upgraded with an x-ray top for unobstructed cassette passage, supporting films up to 14" x 17". Fully tested to four times its rated capacity for patient weights over 600 lbs, it includes a manual brake release, backup controls, and a 4-point auto-leveling braking system for added stability. With adjustable height, Trendelenburg, reverse Trendelenburg, and multiple tilt and flex functions, the Skytron 6600B offers comprehensive flexibility for various surgical needs.

[Read More](#)

SKU: AM-SKYTRON-6600B

Price: Request Quote



BIODEX 058-878 30IN WIDE TOP IMAGING TABLE

The **Biodex 058-878 30" W Top Imaging Table** is a C-Arm imaging table that is a high-stability, three-axis motorized platform designed for precise, quiet positioning and optimal image access. It features a low-attenuation carbon fiber 30" wide top, a 93" overall length, and a 500 lb weight capacity to accommodate a wide range of patients. The rectangular 90" L x 30" W top sits on high-speed actuators with 30"-40" height adjust, 0-20° Trendelenburg and reverse Trendelenburg, and 0-20° lateral roll for versatile patient positioning. The table is equipped with hand control for all movements, locking 5" swivel casters for mobility, and supports both 115 V and 230 V configurations. A compact base and optional battery operation enhance flexibility in the OR, facilitating rapid setup and stable, vibration-free imaging during procedures.

[Read More](#)

SKU: AM-BD-058-878

Price: Request Quote



MORGAN MEDESIGN FLTTLT-FW (FLOAT-TOP, FREE-WHEEL CONFIGURATION)

The **Morgan MEDesign FLTTLT-FW Float-Top, Free-Wheel Imaging Table** is a high-precision fluoroscopy-ready C-Arm compatible platform featuring a float-top, radiolucent surface and multidirectional motion. Its free-floating longitudinal and lateral movement enables micro-adjustments during procedures without repositioning the patient, while a powered height-adjust system provides ergonomic, effortless vertical changes. With a free-wheel base and integrated locking, the table moves easily yet stays securely in place when needed. The radiolucent tabletop, coupled with Trendelenburg and reverse Trendelenburg options on select configurations, supports a wide range of interventional and diagnostic applications. Designed to handle up to approximately 500 lbs (227 kg), it measures about 80" (203 cm) in length and 24" (61 cm) in width, with a height range of 28-40" (71-102 cm). The table weighs around 1,250 lb (567 kg), reflecting its robust construction and suitability for demanding clinical environments.

[Read More](#)

SKU: AM-MORGAN-FLTTLT-FW

Price: Request Quote



STI ECONOMAX 4

The **STI EconoMAX 4** is an economical, versatile mobile imaging table designed for high efficiency and patient comfort in medical facilities such as imaging centers, surgical suites, and pain management clinics. Equipped with motor-driven movements for height adjustment (28" to 38"), lateral tilt (+/-17°), and Trendelenburg tilt (+/-13°), it allows precise positioning tailored to various procedures. The table features a 24-inch wide top with a high-quality, 2-inch thick patient comfort pad, eliminating the need for arm boards in many cases. Its design includes locking swivel casters for easy mobility, a narrow base for straightforward patient transfer, and a facial cut-out suitable for prone procedures. With a standard 500-pound load capacity, the EconoMAX 4 offers excellent value, combining multiple adjustable movements with a low boarding height, making it an ideal choice for busy medical environments seeking affordability without sacrificing functionality.

[Read More](#)

SKU: AM-EC-4

Price: Request Quote



STERIS 5085 SURGICAL TABLE

The **Steris 5085 Surgical Table** is a versatile, adaptable solution designed to meet the needs of virtually every surgical specialty. Its higher weight capacity, exceptionally low height, and expansive imaging area eliminate the need for inflexible, dedicated specialty tables. With quick, one-touch patient positioning, the 5085 leverages an intuitive hand control featuring pre-programmed positions for rotation and beach chair setup, simplifying procedures and reducing setup time. An easy-to-read LCD display communicates precise table positioning data to caregivers, even in dimly lit operating rooms. Key features include a generous patient weight capacity of up to 1,200 lbs, a wide height range from 22.8" to 42.5", and a four-section tabletop that contours to virtually all patient postures. The table also provides unrestricted perineal access with a removable leg section, a powered radiolucent kidney elevator for lateral procedures, and compatibility with a broad range of existing accessories. Designed for safety and efficiency, it includes self-leveling floor locks, auto limit sensors to prevent collisions, and a manual override option. Whether supporting minimally invasive surgery, cardiothoracic, orthopedic, or neurosurgical procedures, the Steris 5085 ensures outstanding access, imaging, and patient support in a reliable, user-friendly package.

[Read More](#)

SKU: AM-5085-R

Price: Request Quote

PATIENT TRANSFER ROLLER BOARD - 67" X 15"

67" x 15"




Mid Central Medical's Patient Transfer Roller Board (67" x 15") is designed to ease staff workloads and reduce lifting during patient transfers. At 67" long, it accommodates larger patients or scenarios requiring a full-length transfer board, allowing smooth movement between beds, gurneys, operating tables, X-ray tables, exam tables, and cast tables. The board features smooth-rolling rollers, a protective cover, and is intended to facilitate safe, effortless transfers with minimal friction and effort.

[Read More](#)

SKU: AM-MCM130

Price: Request Quote



MIDCENTRAL
MEDICAL



SURGICAL TABLES INC. SPINEXR TABLE

Designed to meet the demanding needs of spine procedures, the **STI SpineXR Table** offers superior positioning, safety, and versatility—comparable to the OSI Jackson Table. Custom-built for precision, this table features motorized patient rotation programmable up to 360 degrees, making it ideal for procedures requiring height adjustment, Trendelenburg positioning, and lateral rotation. Its proprietary adjustable length accommodates patients from 4'8" to 6'6", ensuring optimal comfort and support across a wide range of body types. Equipped with carbon fiber support rails compatible with other brand accessories, four swivel, total-lock castors, and user-friendly controls, the SpineXR Table delivers both safety and convenience. Key features include an audible power-up message, emergency stop, illuminated power switch, and a removable power cord for ease of use. With a robust 500 lb weight capacity and dimensions of 88"-102"L x 29.5"W, this table is a reliable addition to any surgical suite.

[Download Brochure](#)

[Download the STI SpineXR Accessory List](#)

***** Accessories in photo sold separately, inquire for more information to assemble your table today. *****

[Read More](#)


SKU: AM-STI-SPINE-XR

Price: Request Quote



MINDRAY HYBASE V6 SURGICAL TABLE W/ LONGITUDINAL SLIDE AND CARBON FIBER EXTENSION (STAINLESS STEEL BASE)

The **Mindray HyBase V6 Surgical Table with Longitudinal Slide and Carbon Fiber Extension (Stainless Steel Base)** is designed to meet the complex demands of modern surgical procedures. With an impressive load capacity of up to 460 kg (1014 lbs), it accommodates a wide range of patients, ensuring stability and safety during operations. The table permits all positions up to an overall load of 250 kg (551.1 lbs), offering versatility for various surgical needs. Its advanced design includes a longitudinal shift of up to 350mm (13.8"), facilitating seamless access to C-arm imaging, complemented by head-end and foot-end shifts of 6.3" and 7.5" respectively. The table's height adjusts from a low 600mm (23.6") for neurosurgery to a high 1,050 mm (41.3") for orthopedic procedures. It features 12 electric movements, including table top adjustments and tilts, all controlled via a Bluetooth hand control for ease of use in the operating room. Ensuring maximum safety, the table




Carbon fiber tabletop



Radiolucent carbon fiber tabletop provides better imaging.

Huge Imaging Access



includes an independent manual override system, a triple-layer decompression pad for patient comfort, and IPX5 water resistance. Its bellows-free design, coupled with four double-swivel castors, enhances maneuverability. The package includes arm boards, an anesthesia frame, memory foam pads, a stainless steel base cover, a longitudinal slide, a kidney lift, and a carbon fiber extension that optimizes imaging areas at both the head (64.9") and foot (56.7") ends, making it a comprehensive solution for healthcare providers.

[Read More](#)

SKU: AM-121-002775-00

Price: Request Quote

MINDRAY HYBASE V6 SURGICAL TABLE W/ LONGITUDINAL SLIDE (STAINLESS STEEL BASE)

The **Mindray HyBase V6 Surgical Table with Longitudinal Slide (Stainless Steel Base)** is designed to meet the diverse needs of surgical teams. With an impressive load capacity of up to 460 kg (1014 lbs) and the ability to accommodate all positions up to an overall load of 250 kg (551.1 lbs), this table ensures stability and safety during complex procedures. It features a longitudinal shift of up to 350mm (13.8"), providing easy access to C-arm imaging, crucial for precision in surgical interventions. The table offers a head-end shift of 6.3" and a foot-end shift of 7.5", with a versatile height range from 600mm (23.6") for neurosurgical procedures to 1,050mm (41.3") for orthopedic operations. With 12 electric movements, including table top adjustments and tilt functions, the HyBase V6 adapts seamlessly to various surgical requirements. Its Bluetooth hand control offers unobstructed management in the operating room, while an independent manual override system enhances safety. The table is equipped with a triple-layer decompression pad for optimal patient comfort and weight distribution, and boasts IPX5 protection against water jets. Its bellows-free design and four double-swivel castors enhance mobility. Inclusions such as arm boards, an anesthesia frame, memory foam pads, a stainless steel base cover, a longitudinal slide, and a kidney lift add to its functionality, with additional options like a carbon fiber table top available for customization.

[Read More](#)

SKU: AM-121-002774-00

Price: Request Quote



Bidirectional longitudinal shift



Head-end shift: 160 mm (6.3")
Foot-end shift: 190 mm (7.5")

Built-in body elevator



Built-in body elevator for kidney surgery.



MINDRAY HYBASE V6 SURGICAL TABLE (STAINLESS STEEL BASE)

The **Mindray HyBase V6 Surgical Table**, featuring a robust **stainless steel base**, is a marvel of modern medical engineering designed to accommodate diverse surgical needs with precision and safety. With an impressive extreme load capacity of up to 460 kg (1014 lbs), it ensures stability and reliability during complex procedures. Its versatile design allows for all positions to be maintained with an overall load of up to 250 kg (551.1 lbs), catering to a wide array of surgical requirements. The table's height adjustability, ranging from a low 600 mm (23.6") ideal for neurosurgery to a high 1,050 mm (41.3") optimized for orthopedic procedures, ensures ergonomic access for surgeons. The table offers 12 electric movements, allowing seamless adjustments such as tabletop elevation, slope, and tilt, all of which can be controlled unobstructedly via a Bluetooth hand control. Safety is enhanced with an independent manual override system, complemented by a triple-layer decompression pad that redistributes the patient's weight for comfort. Its IPX5 protection rating guards against water jets from any direction, ensuring durability. The bellows-free design, stainless steel base cover, and four double-swivel castors facilitate easy maneuverability. Standard inclusions such as a pair of arm boards, an anesthesia frame, and memory foam pads enhance functionality, while optional features like a longitudinal slide, kidney lift, and carbon fiber tabletop offer additional customization.



[Read More](#)

SKU: AM-121-002773-00

Price: Request Quote


MINDRAY HYBASE V6 SURGICAL TABLE (PLASTIC BASE)

The **Mindray HyBase V6 Surgical Table with a plastic base** is a state-of-the-art solution designed to meet the diverse needs of modern surgical environments. With an impressive load capacity of up to 460 kg (1014 lbs), this surgical table accommodates a wide range of patients. It permits all positions up to an overall load of 250 kg (551.1 lbs), ensuring versatility and safety. The table features a remarkable height adjustment range, from a low of 600mm (23.6") suitable for neurosurgery to a high of 1,050 mm (41.3") ideal for orthopedic procedures. Its 12 electric movements, including tabletop up/down, slope, and tilt, provide adaptability for various surgical requirements. A Bluetooth hand control allows for unobstructed operation within the operating room, while an independent manual override system offers dual safety assurance. The triple-layer decomposition pad ensures optimal weight redistribution for patient comfort. Additionally, the table is equipped with IPX5 protection against water jets, ensuring durability and hygiene. The bellows-free design and plastic base cover, complemented by four double-swivel castors, enhance mobility and ease of use. Standard inclusions such as a pair of arm boards, an anesthesia frame, and memory foam pads offer comprehensive support, while optional features like a stainless steel base cover, longitudinal slide, kidney lift, and carbon fiber tabletop provide further customization.

[Read More](#)

SKU: AM-121-002772-00

Price: Request Quote





MAYFIELD HEADREST - BASIC HEADREST

The **Mayfield Headrest - Basic Headrest** is expertly designed to offer superior cranial stabilization in both prone and supine positions, ensuring optimal patient comfort and safety during medical procedures. It features both lateral and vertical mobility, allowing for precise positioning of the patient's head, which aids in better distribution of weight and forces to effectively combat pressure necrosis. The headrest comes with an A1011/1012 Adult Horseshoe headrest and a powder-coated frame for enhanced durability. To maintain hygiene and extend the product's life, replacement gel pads are available individually. Additionally, the headrest is equipped with a new vinyl foam pad, providing an extra layer of comfort. The product is backed by a six-month warranty covering any manufacturer's defects, underscoring its quality and reliability.

[Read More](#)

SKU: AM-MFHR-Basic

Price: Request Quote



STERIS NARROW HEADREST LIGHTWEIGHT, RADIOLUCENT HEAD SUPPORT

BF592- The radiolucent, Narrow Headrest with its elongated tapered edge provides easier temporal and over-the-brow access to head area for ophthalmic, ENT, neuro and plastic surgery. This fixed height headrest fits into the STERIS surgical table's back section or leg section by way of a required Eye-ENT-Neuro Adaptor. The Eye-ENT-Neuro Headrest Adaptor is sold separately. This accessory includes a 2 1/4" pressure management pad.

Special Note(s): The Eye-ENT-Neuro Headrest Adaptor is required (sold separately). Order BF185 for 3085 tables, BF173 for 5085, 4085, and Cmax tables, or BF309 for 5085SRT.

[Read More](#)

SKU: AM-BF592

Price: Request Quote



STERIS 9 INCH UROLOGY EXTENSION

This lightweight 9-inch long radiolucent extension provides additional fluoro imaging length for imaging of kidney, ureter and bladder (KUB) anatomy of tall patients. It easily attaches to the seat section of the table once the leg section is removed, and includes standard side rails for applying leg holders or Fluid Collector System. The perineal cutout will accommodate a metal drain pan.

For Steris 4085, 5085 Tables only

[Read More](#)

SKU: AM-BF578

Price: Request Quote



MIZUHO OSI 6875 HANA TABLE

The **Mizuho OSI 6875 Hana Table** is a versatile and specialized surgical table designed to enhance orthopedic procedures, particularly those involving the lower extremities. Equipped with a comprehensive range of accessories, it ensures precision, safety, and ease of use during complex surgeries. Key components include the Femoral Hook Support and Classic attachments for both left and right positioning, facilitating optimal femoral management. The Femur Lift Foot Pedal and its extensions provide controlled elevation of the femur, critical for accurate surgical access and alignment. Additional features like the Traction Hook Extension and Radiolucent Leg Spars support seamless traction and imaging compatibility. The table's ergonomic design includes Arm Boards with Pads, the Hana Knee Flexion System for total knee arthroplasty, and a Well Leg Support Adaptor, ensuring patient stability and comfort. Safety is prioritized with elements such as the Pelvic Pad, Perineal Post and Pad, and the Patent Safety Strap. With its Articulating Bracket and Well Leg Upright Support Assembly, the Hana Table embodies a comprehensive solution for modern orthopedic surgery, enhancing both surgical outcomes and operational efficiency.

[Read More](#)

SKU: AM-OSI-6875

Price: Request Quote



STERIS SURGIGRAPHIC 1027 C-ARM IMAGING TABLE

The **Steris SurgiGraphic 1027 C-Arm Imaging Table** is a highly advanced surgical platform designed for optimal patient imaging and versatility in the operating room. Constructed with a carbon fiber tabletop supported by dual columns and a mobile base assembly, this table allows for comprehensive 62" (1575 mm) patient body scans, enhancing image-guided procedures. Its unique 45° edge contour facilitates easy accessory attachment while enabling 360° patient imaging capability with a fluoroscopic C-Arm. The table measures 28" (711 mm) in width, 91" (2311 mm) in length, with an adjustable height ranging from 28" to 41" (711 mm to 1041 mm). It supports a maximum patient weight of 500 lbs (226 kg), ensuring stability and reliability. The standard configuration boasts features such as a Tri-Level Technology (TLT) tabletop pad, adjustable height, Trendelenburg and lateral tilt functions, electric power, pendant hand control, and a mobile base with locking casters. Additionally, it includes a 2-inch standard comfort pad and operates at 115 VAC, 60 Hz, 3.5 Amp. The table's motion capabilities include height adjustments from 28" to 41" (711 mm to 1041 mm), Trendelenburg/reverse Trendelenburg from 0-10°, and lateral tilt from 0-17°, making it a versatile and essential tool for modern surgical environments.

[Read More](#)

SKU: AM-SG-1027

Price: Request Quote



OSI HANA EQUIPMENT CART

Convenient Storage of OSI 6875 Accessories

[Read More](#)

SKU: 6875-500

Price: Request Quote



SKYTRON 3603 ULTRA SLIDE SURGICAL TABLE

The **Skytron Ultra Slide 3603** represents the visionary benchmark in Top Slide surgical table design, including even greater imaging versatility, lift and articulation power. Expanded features provide optimal patient positioning and C-arm access with 26" Top Slide, 800lb lift, articulation and full table functions.

[Read More](#)

SKU: AM-SKYTRON-3603

Price: Request Quote



STERIS AMSCO MODEL 3085 SURGICAL TABLE

The Amsco 3085 General Surgical Table is one of the most reliable and popular surgical tables. Recognized as the industry standard due to its sheer versatility and reliability, the Steris Amsco 3085 has been trusted for years by healthcare facilities across the world.

[Download Brochure](#)

[Technical Specs](#)

[Read More](#)

SKU: AM-3085SP

Price: Request Quote



BIODEX 840 FLOATING TOP IMAGING TABLE

With an advanced low attenuation float top for cardiovascular, endoscopic and pain management procedures, the Surgical C-Arm Table 840 is the ideal choice when stability, access, and precise, quiet, vibration-free positioning are essential. A cantilevered low attenuation carbon-fiber tabletop, with choice of rectangular or contoured design, accommodates portable or ceiling-suspended C-Arms.

[Read More](#)

SKU: AM-BD-058-840

Price: Request Quote




TABLE PROTECTIVE DRAPE CFUR PROTECTIVE DRAPE- CLEAR PVC- ATTACHES TO T-RAIL

[Read More](#)

SKU: 69390

Price: \$135.00




TABLE PROTECTIVE DRAPE CFLU PROTECTIVE DRAPE- CLEAR PVC- ATTACHES TO T-RAIL

[Read More](#)

SKU: 69389

Price: \$135.00



CFPM300 TABLE 3 MOTION- RECTANGULAR TOP- 2IN PAD, HAND CONTROL

- Integrated Head Rest
- Large Adjustable Height Range
- 56 in. Metal Free Imaging Space
- Cost Effective 3 Movement Table
- 500 lb. Patient Weight Capacity

[Read More](#)

SKU: 81806

Price: Request Quote



CFUR401 UROLOGY TABLE 4 MOTION- CARBON FIBER TOP- 2IN PAD, HAND CONTROL & FOOT CONTROL

36 month warrantee period

[Read More](#)

SKU: 83867

Price: Request Quote



CFUR301 UROLOGY TABLE 3 MOTION- CARBON FIBER TOP- 2IN PAD, HAND CONTROL & FOOT CONTROL

Carbon Fiber Top - 53 in. (135cm) Metal free imaging space with Patient Transfer Extension

Top Dimensions - 23 in. (58cm) wide x 47 in. (119cm) long

[Read More](#)

SKU: 83868

Price: Request Quote



CFPMB301 TABLE 3 MOTION- 750 LBS RATED, RECTANGULAR TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 81812

Price: Request Quote



CFPM401 TABLE 4 MOTION- RECTANGULAR TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures. Available with an integrated headrest or rectangular top.

[Read More](#)

SKU: 81809

Price: Request Quote



CFPM401 TABLE 4 MOTION- INTEGRATED HEADREST TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures. Available with an integrated headrest or rectangular top.

[Read More](#)

SKU: 81802

Price: Request Quote



CFPM400 TABLE 4 MOTION- INTEGRATED HEADREST TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures. Available with an integrated headrest or rectangular top.

[Read More](#)

SKU: 81772

Price: Request Quote



CFPM302 TABLE 3-MOTION- RECTANGULAR TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 84532

Price: Request Quote



CFPM302 TABLE 3-MOTION- INTEGRATED HEADREST TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 84530

Price: Request Quote



CFPM301 TABLE 3 MOTION- RECTANGULAR TOP- 2IN PAD HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 81811

Price: Request Quote



CFPM301 TABLE 3 MOTION- INTEGRATED HEADREST TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 81804

Price: Request Quote



CFPM300 TABLE 3 MOTION- RECTANGULAR TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 81810

Price: Request Quote



CFPM300 TABLE 3 MOTION- INTEGRATED HEADREST TOP- 2IN PAD, HAND CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.

[Read More](#)

SKU: 81803

Price: Request Quote



CFPM100 TABLE 1 MOTION- RECTANGULAR TOP- 2IN PAD FOOT CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures. Available with an integrated headrest or rectangle top.

[Read More](#)

SKU: 81813

Price: Request Quote




CFPM100 TABLE 1 MOTION- INTEGRATED HEADREST TOP- 2IN PAD, FOOT CONTROL

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures. Available with an integrated headrest or rectangle top.

[Read More](#)

SKU: 81805

Price: Request Quote



CFLU401 LITHOTRIPSY/UROLOGY TABLE 4 MOTION- CARBON FIBER TOP- 2IN PAD HAND CONTROL & FOOT CONTROL

Carbon Fiber Top - 53 in. (135cm) Metal free imaging space with Patient Transfer Extension
Top Dimensions - 23 in. (58cm) wide x 47 in. (119cm) long

[Read More](#)

SKU: 82014

Price: Request Quote

IDI 100RT TABLE

C-Arm table for GI and pain management procedures where motorized elevation, tilt, and lateral roll are required.

[Read More](#)

SKU: AM-IDI-100RT-3

Price: Request Quote





SKYTRON 6600 SURGICAL TABLE

The **Skytron 6600 general purpose surgical table** provides unparalleled imaging access for today's most demanding cases while lifting 1000 pounds and articulating 600 pounds. The removable Back and Leg sections facilitate the use of a carbon fiber Back and Leg sections to further extend the imaging window. The rotating top of 210° further enhances C-arm and microscope access along with more leg room for a seated surgeon and in many cases, better surgical access.

[Read More](#)

SKU: AM-SKYTRON-6600

Price: Request Quote



STERIS CMAX SURGICAL TABLE

Table Length 81" (2057 mm)

Table Width 20" (508 mm) (28" with extenders)

Table Height Range 26" to 45" (660 to 1143 mm)

Table Slide Range 9" (227 mm) to head, 9" to foot (18" total)

Patient Weight Capacity 1,100 lbs. (500 kg) patient support, including raise/lower (centered on the column)

[Read More](#)

SKU: AM-CMAX

Price: Request Quote

SKYTRON 6702 SURGICAL TABLE



Currently being sold new by Skytron, the 6702 model has been manufactured since 2013

Skytron 6702 Hercules allow 210 degree top rotation at the turn of a handle, making procedure and imaging setup easier than ever. Seamlessly rotate the table toward the C-arm or transition the patient between the surgeon and anesthesiologist during head-end procedures. Featuring 1,200 pound lift and 1,000 pound articulation capacity combined with our lowest table height, the 6702 Hercules can also support bariatric patients with comfortable surgical access. Its radiolucent top permits sharp image amplification while the sturdy side rails will accept a full range of complementary accessories. From its ease of setup to its impressive power, the 6702 Hercules helps any operating room manage a greater variety of procedures with only one surgical table!

The 6702 hercules also brings confidence to the operating room. Recognized for their reliable performance, Skytron surgical tables are consistently rated as having a lower cost of ownership.

[Read More](#)

SKU: AM-6702

Price: Request Quote

BERCHTOLD OPERON OR TABLE D820



The Stryker Berchtold Operon D820 is an advanced and feature-rich surgical table. The Operon D820 offers the best-in-class range of vertical travel, full weight capacity throughout its entire range of motion and superior imaging access. This leads to increased patient positioning options which translate into greater comfort, productivity, and safety for the surgical team.

[Read More](#)

SKU: AM-BERCH-D820

Price: Request Quote

BERCHTOLD OPERON OR TABLE D850



The Berchtold Operon 850 General Surgical Table is a very versatile, easy-to-use surgical table. With a weight limit of 1000 lbs in normal orientation and 500 lbs in reverse orientation, the Operon® D 850 can be used for almost any general surgical specialty, including bariatric applications. Features a small footprint for superior patient and equipment access. The Berchtold Operon D850 is supplied with a removable headrest unit, pendant hand control unit, and a set of electrically conductive table top cushions.

[Read More](#)

SKU: AM-BERCH-D850

Price: Request Quote

MAQUET ALPHAMAXX 1133 SURGICAL TABLE



With a few positioning adjustments and optional accessories, the Maquet Alphamaxx 1133 is ideal for operations in gynecology, urology, orthopedics, traumatology and even neurosurgery. The table features electro-hydraulic actuators that control longitudinal shift, lateral tilt, and various table angles. The Alphamaxx 1133 adjusts from 23.4 to 41.6 inches in height, allowing clinicians to work comfortably from a sitting or standing position.

An easy click mounting system allows for a wide range of accessories to be quickly added to the table, efficiently transforming the chair for each unique procedure. The durable table's maximum load capacity is rated at 990 pounds.

[Read More](#)

SKU: AM-AMAXX-1133

Price: Request Quote

NUVO SURGICAL V1000 SURGICAL TABLE SERIES



Designed to accommodate surgical procedures in nearly all fields, including urology, neurology, orthopedics, and cardiology, the Nuvo V1000 series surgical tables set new standards for versatility, innovative features and solid dependability.

With design and engineering thought through to a remarkable degree, they provide such benefits as unprecedented C-arm access, self-diagnosis with Instant Error Code display, and a control system that automatically prevents of sectional conflicts.

[Read More](#)

SKU: AM-B2AV1000K-1

Price: Request Quote



SKYTRON 3600B SURGICAL TABLE

The UltraSlide 3600B represents the visionary benchmark in Top Slide surgical table design, including even greater imaging versatility, lift and articulation power. Expanded features provide optimal patient positioning and C-arm access with 23" Top Slide, 1,000 lb. lift, 800 lb. articulation and full table functions. Additional features include removable back and leg sections, Auto Beach Chair position, 30° tilt, Trendelenburg and reverse Trendelenburg, and more! Optional Infrared (IR) Wireless Table Control is also available.

Skytron Part Number: 3600B Ultraslide

[Download Brochure](#)

[Read More](#)

SKU: AM-SKYTRON-3600B

Price: Request Quote



SURGICAL TABLES INC. V-MAX 500 LBS CAPACITY (NO TREND)

The **Surgical Tables Inc. V-MAX-T HEIGHT, 4 WAY FLOAT Lateral Travel Table** is a versatile, new surgical table designed for ease of use and optimal patient access. Featuring automatic leveling of the table top, a diving board style for unobstructed C-arm access, and a generous 500 lb weight capacity, it provides a reliable platform for various surgical procedures. The table includes a 2" high-quality comfort pad, facial cutout access, and a low height for easy patient boarding. Its 24" wide top (with an optional 27" width) eliminates the need for arm boards and allows for easy patient transfer. The manual 4-way float includes 36" of longitudinal travel (head to toe) and 14" of lateral movement, all controllable via a handle for smooth positioning. Compatible with any OR accessory, this lightweight, easy-to-move table is ideal for diverse surgical environments.

[Read More](#)

SKU: AM-STI-V-MAX

Price: Request Quote

STI UROMAX SERIES



Our URO-MAX C-arm imaging tables offer an excellent solution for real-time urological X-ray procedures. The cantilevered, carbon fiber, radiolucent top areas allow full fluoroscopic visualization and unobstructed C-arm positioning. The URO-Max is highly versatile with many movement options and accessories, and the carbon fiber leg extension accessory can quickly turn your Urology table into a procedure table for pain management and other C-arm guided procedures.

[Read More](#)

SKU: N/A

Price: Request Quote

STREAMLINE 5 PAIN TABLE



Designed for high patient out put for surgery centers, pain management, and other multidisciplinary medical facilities, the Streamline series X-ray imaging tables are the lowest priced tables in the imaging table industry.

[Read More](#)

SKU: AM-SL-5

Price: Request Quote



STI ECONOMAX 3

EconoMAX 3 moves (up/down, lateral tilt, and trendelenburg)

These tables come with a standard 500-pound load rating, a 2" thick, high quality patient comfort pad, a 24" wide top that can eliminate the need for arm boards, locking swivel castors, a narrow base for easy patient transfer, a facial cut-out for prone procedures, simultaneous motion movements, and extremely low boarding height. EconoMAX has the highest patient weight load rating at the best price. The EconoMax series offers the best value on the market. With many user-friendly accessories and options to choose from, the EconoMAX mobile economy imaging tables offer the perfect combination of features and price.

[Read More](#)

SKU: AM-EC-3

Price: Request Quote



STI STREAMLINE 3 PAIN TABLE

The **Streamline 3 Pain Table** from Surgical Tables Inc. is a cost-effective and versatile imaging table designed to meet the demands of high patient throughput in surgery centers, pain management clinics, and multidisciplinary medical facilities. Featuring a carbon fiber, diving board-style design, it allows for easy access for C-arm imaging, making it suitable for various diagnostic and interventional procedures. The table supports a patient weight capacity of up to 500 pounds and offers a spacious 84" x 24" tabletop with a 2-inch thick pad, including a facial cutout for patient comfort. Adjustable in height from 31" to 41", with lateral tilt options up to 17° and Trendelenburg tilt also up to 17°, the table provides flexible positioning to accommodate diverse surgical and pain management procedures. Its affordability combined with functionality makes it a practical choice for busy medical environments.

[Read More](#)

SKU: AM-SL-3

Price: Request Quote

STREAMLINE 2 PAIN TABLE



Designed for high patient out put for surgery centers, pain management, and other multidisciplinary medical facilities, the Streamline series X-ray imaging tables are the lowest priced tables in the imaging table industry.

[Read More](#)

SKU: AM-SL-2

Price: Request Quote

STREAMLINE 1 PAIN TABLE



Designed for high patient out put for surgery centers, pain management, and other multidisciplinary medical facilities, the Streamline series X-ray imaging tables are the lowest priced tables in the imaging table industry.

[Read More](#)

SKU: AM-SL-1

Price: Request Quote

STREAMLINE PAIN TABLE, FIXED HEIGHT



Designed for high patient throughput for surgery centers, pain management, and other multidisciplinary medical facilities, the Streamline series X-ray imaging tables are the lowest priced tables in the imaging table industry.

Features Include:

- Carbon fiber, diving board style for easy C-Arm access
- 500 lb. patient capacity
- 84" x 24" table top
- 2" thick pad with facial cutout

Fixed Height

[Read More](#)

SKU: AM-SL-FH

Price: Request Quote

WILSON SPINAL FRAME RADIOLUCENT



The Wilson Plus™ Radiolucent Frame delivers unrestricted radiolucency and C-arm integration to your OR Table with a convenient and stable method of maintaining patients in a flexed position for lumbar procedures.

[Download Brochure](#)

[Read More](#)

SKU: AM-WS-9015

Price: Request Quote

NUVO SURGICAL V800



Designed to accommodate surgical procedures in nearly all fields, including urology, neurology, orthopedics, and cardiology, the Nuvo V800 series surgical tables set new standards for versatility, innovative features and solid dependability.

With design and engineering thought through to a remarkable degree, these OR Tables provide such benefits as unprecedented C-arm access, self-diagnosis with Instant Error Code display, and a control system that automatically prevents sectional conflicts.

[Download Brochure](#)

[Read More](#)

SKU: AM-B2AV800SK-1

Price: Request Quote



STERIS 4085

The STERIS 4085 General Surgical Table's design gives a surgical team the flexibility, control and convenience they need for today's patients and procedures.

[Download Brochure](#)

To submit a quote request for this product please fill out the form below.
Product Quote Request

Name

*

Name

Name

Company Name

*

Email

*

Phone

*

Shipping Zip Code

*

Quantity

*

Quantity *

Submit Quote Request for Product

Δ

[Read More](#)**SKU:** AM-S4085**Price:** Request Quote

SKYTRON 3501B

The EZ Slide 3501B Specialty table provides an impressive range of full-body imaging and articulating functions without the need to reposition the patient or resort to multiple support accessories. This Surgical Table features a 45-inch maximum and 26-inch minimum table height with 21-inch top side travel.

[Download Brochure](#)[Read More](#)**SKU:** AM-SKYTRON-3501B**Price:** Request Quote

SKYTRON 3500 ELITE



The Skytron 3500 Elite features an incredible 21" top slide, 500 lb lift/articulation, return-to-level and return-to-center. An ultra modern operating table, able to satisfy all of your C-Arm and general surgery needs.

[Download Brochure](#)

[Read More](#)

SKU: AM-SKYTRON-3500

Price: Request Quote

SKYTRON 6700



The Skytron 6700 general surgical table represents the new standard in 21st Century Surgical Table design and performance, with unmatched imaging flexibility and even greater lift and articulation power. With it's 1,200 lb patient weight capacity and the ability to be widened up to 16", it's a very popular operating table for bariatric surgery.

[Download Brochure](#)

[Read More](#)

SKU: AM-SKYTRON-6700B

Price: Request Quote



STERIS AMSCO 3080 SURGICAL TABLE

The Amsco® 3080 Surgical Table features a radiolucent top and can be adjusted to the same variety of positions as the 2080 model. It comes standard with a remote control for easy and efficient patient positioning. This is one of the most durable, flexible, and reliable surgical tables available.

[Download Brochure](#)

[Read More](#)

SKU: AM-3080

Price: Request Quote



SKYTRON 6500

The Skytron® 6500 Surgical Table is a general surgery table filled with features. The 500 lb. (Elite) and 850 lb. (Hercules) weight capacity and the C-Arm compatibility allow this OR table to be put into any operating room. This table can also adjust to the trendelenburg, back, side, foot, and flex positions.

[Download Brochure](#)

[Read More](#)

SKU: AM-6500

Price: Request Quote



SKYTRON HERCULES 6701 SURGICAL TABLE

From ease of set-up to the reduced need for multiple support accessories, the Hercules 6701 is the most powerful, versatile, and image-friendly surgical table available. The Hercules 6701 surgical table represents the new standard in 21st Century Surgical Table design and performance, with unmatched imaging flexibility and even greater lift and articulation power. The Hercules 6701 is a operating table that can be used for general surgery or for bariatric patients.

[Download Quick Product Sheet](#)

[Download Manufacturer Brochure](#)

[Read More](#)

SKU: AM-SKYTRON-6701

Price: Request Quote

INDEX

B

Berchtold OPERON OR Table D820 26
Berchtold OPERON OR Table D850 27
Biodex 058-878 30in Wide Top Imaging Table 3
Biodex 840 Floating Top Imaging Table 15

C

CFLU401 Lithotripsy/Urology Table 4 motion- Carbon Fiber Top- 2in Pad Hand Control & Foot Control 24
CFPM100 Table 1 motion- Integrated Headrest Top- 2in Pad, Foot Control 23
CFPM100 Table 1 motion- Rectangular Top- 2in Pad Foot Control 23
CFPM300 Table 3 motion- Integrated Headrest Top- 2in Pad, Hand Control 22
CFPM300 Table 3 motion- Rectangular Top- 2in Pad, Hand Control 16, 22
CFPM301 Table 3 motion- Integrated Headrest Top- 2in Pad, Hand Control 21
CFPM301 Table 3 motion- Rectangular Top- 2in Pad Hand Control 21
CFPM302 Table 3-motion- Integrated Headrest Top- 2in Pad, Hand Control 20
CFPM302 Table 3-motion- Rectangular Top- 2in Pad, Hand Control 20
CFPM400 Table 4 motion- Integrated Headrest Top- 2in Pad, Hand Control 19
CFPM401 Table 4 motion- Integrated Headrest Top- 2in Pad, Hand Control 19
CFPM401 Table 4 motion- Rectangular Top- 2in Pad, Hand Control 18
CFPMB301 Table 3 motion- 750 lbs rated, Rectangular Top- 2in Pad, Hand Control 18
CFUR301 Urology Table 3 motion- Carbon Fiber Top- 2in Pad, Hand Control & Foot Control 17
CFUR401 Urology Table 4 motion- Carbon Fiber Top- 2in Pad, Hand Control & Foot Control 17

I

IDI 100RT Table 24

M

Maquet Alphamaxx 1133 Surgical Table 27
Mayfield Headrest - Basic Headrest 11
Mindray HyBase V6 Surgical Table (Plastic Base) 10
Mindray HyBase V6 Surgical Table (Stainless Steel Base) 9
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide (Stainless Steel Base) 8
Mindray HyBase V6 Surgical Table w/ Longitudinal Slide and Carbon Fiber Extension (Stainless Steel Base) 7
Mindray V6 Surgical Table Extension 2
Mizuho OSI 6875 Hana Table 12
Morgan MEDesign FLTLT-FW (Float-Top, Free-Wheel Configuration) 4

N

Nuvo Surgical V1000 Surgical Table Series 28
NUVO Surgical V800 33

O

OSI Hana Equipment Cart 13

P

Patient Transfer Roller Board - 67" x 15" 5

S

Skytron 3500 Elite 36
Skytron 3501B 35
Skytron 3600B Surgical Table 28
Skytron 3603 Ultra Slide Surgical Table 14
Skytron 6500 37
Skytron 6600 Surgical Table 25
Skytron 6600B Surgical Table 3
Skytron 6700 36
Skytron 6702 Surgical Table 26
Skytron Hercules 6701 Surgical Table 38
STANDARD SERIES - SURGICAL TABLE PAD SETS - FOR: 3080RC - 3080RL - 3080SP - 3085SP 2
Steris 4085 34
Steris 5085 Surgical Table 5
Steris 9 inch Urology Extension 12
Steris Amsco 3080 Surgical Table 37
Steris Amsco Model 3085 Surgical Table 14
Steris CMAX Surgical Table 25
Steris Narrow Headrest Lightweight, Radiolucent Head Support 11
Steris SurgiGraphic 1027 C-arm Imaging Table 13
STI EconoMAX 3 30
STI EconoMAX 4 4
STI Streamline 3 Pain Table 31
STI Uromax Series 29
Streamline 1 Pain Table 32
Streamline 2 Pain Table 31
Streamline 5 Pain Table 30
Streamline Pain Table, Fixed Height 32
Surgical Tables Inc. SpineXR Table 7
Surgical Tables Inc. V-Max 500 lbs capacity (No Trend) 29

T

Table Protective Drape CFLU Protective Drape- Clear PVC- attaches to T-Rail 16
Table Protective Drape CFUR Protective Drape- Clear PVC- attaches to T-Rail 15

W

Wilson Spinal Frame Radiolucent 33