

TABLE OF CONTENTS

Tourniquets	2
Stryker SmartPump Dual Channel Tourniquet	3
Zimmer ATS 4000 Tourniquet	3
Zimmer ATS Dual Blue Hose	4
ATS Dual Red Hose	4
Cuff Cylindrical 34 Inch ATS Zimmer	5
Cuff Large 42 Inch ATS Zimmer	5
Zimmer ATS1200 Single Channel Tourniquet System	6
Zimmer ATS 3000	6
Zimmer ATS 2000	7
Zimmer ATS 1500	8



TOURNIQUETS

Auxo Medical carries a wide range of tourniquets that help improve safety, precision, and speed during medical examinations and surgical procedures.

Offering a wide range of controls and features, get the compression you need with our industry leading tourniquets. It can assist in any emergent or non-emergency situation or procedure.

Our selection of factory level refurbished tourniquets offers numerous dual port, dual cuff systems with controls to measure pressure accurately. Also, several models offer precise microprocessors. This includes alarms to warn medical professionals of any hazardous conditions. In addition, it features like easy-to-read displays, touch pad controls, fast cuff inflation or deflation, and one-hour battery backup. Plus, with disposable and reusable cuffs to fit all patient sizes, our leading selection keeps you and your facility operating smoothly and efficiently in all situations.

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

RECONDITIONED

REFURBISHED



STRYKER SMARTPUMP DUAL CHANNEL TOURNIQUET

The **Stryker SmartPump Dual Channel Tourniquet System** (Reference #5920011000) offers advanced, continuous monitoring and precise control of tourniquet pressure during surgical procedures. Designed to support both single and dual-channel configurations, it seamlessly integrates with OR software systems, allowing for real-time data communication. This versatile system accommodates a variety of surgical needs, including single cuff, Bier Block, bi-lateral, and simultaneous upper and lower extremity procedures. Each cuff's pressure and inflation duration are independently adjustable, with real-time display and monitoring to enhance safety and efficiency. Post-operation, tourniquet data can be printed and incorporated into patient records, making the SmartPump system an essential tool for modern surgical environments.

[Read More](#)

SKU: AM-STRY-SP-DC

Price: Request Quote

ZIMMER ATS 4000 TOURNIQUET

The Zimmer ATS 4000 Tourniquet System doesn't just stop blood flow; we strive to personalize pressure for the patient by delivering the minimum tourniquet pressure at a specific time, with a specific cuff, for a specific patient limb. It has been found that when applying a lower tourniquet pressure is results in less post-operative pain. Using Personalized Pressure Technology the ATS 4000 can minimize tourniquet pressure and pressure gradients. This helps reduce the risk of nerve-related injuries.

[Read More](#)

SKU: AM-ATS4000

Price: Request Quote



ZIMMER ATS DUAL BLUE HOSE

Dual Channel Blue Hose

[Read More](#)

SKU: 60-4017-002QB

Price: Request Quote



ATS DUAL RED HOSE



Zimmer Dual Channel Red Hose

[Read More](#)

SKU: 60-4017-002QR

Price: Request Quote

CUFF CYLINDRICAL 34 INCH ATS ZIMMER



34 INCH DPSB ATS CYLINDRICAL CUFF AUTOMATIC TOURNIQUET SYSTEMS
(A.T.S.) REUSABLE CYLINDRICAL CUFFS (NON- STERILE)

[Read More](#)

SKU: 60-7500-006-00

Price: Request Quote

CUFF LARGE 42 INCH ATS ZIMMER



DPSB CONTOUR CUFF - XL (30 IN-39 IN) AUTOMATIC TOURNIQUET SYSTEMS (A.T.S.) XL CONTOUR CUFFS (REUSABLE, NON- STERILE)

[Read More](#)

SKU: 60-7500-008-00

Price: Request Quote

ZIMMER ATS1200 SINGLE CHANNEL TOURNIQUET SYSTEM



The Zimmer A.T.S. 1200 Tourniquet System provides accurate pressure monitoring for patient safety and confidence. Microprocessor controls and built in safety features with self-diagnostic tests make monitoring reliable in a single port, dual cuff surgical tourniquet system.

[Read More](#)

SKU: AM-ATS1200

Price: Request Quote

ZIMMER ATS 3000

The **Zimmer A.T.S.® 3000 Automatic Tourniquet System** is a state-of-the-art dual-port, dual-cuff system designed with microprocessor controls and dedicated ports that supply and measure pressure independently. This advanced system features the innovative Limb Occlusion Pressure (LOP) technology, which accurately detects, calculates, and reports the precise cuff pressure required to achieve complete blood occlusion in the operative limb. The A.T.S.® 3000 Tourniquet exemplifies Zimmer's long-standing commitment to safety, reliability, and convenience, boasting a four-hour battery backup to maintain pressure during power loss, along with visual and audible alarms for enhanced safety. Additional features include self-check calibration, a lightweight and portable design, color-coded cuffs ranging from 8 to 42 inches, and a convenient carrying handle. The system also ensures rapid inflation and deflation, independent pressure monitoring, and a cuff lockout feature that prompts users to confirm deflation commands, reducing the risk of complications from sudden cuff deflation. Built-in safety features and positive locking connectors for leak-free inflation make the A.T.S.® 3000 a versatile and dependable choice for use on a table or IV pole.

[Download Brochure](#)

[Read More](#)

SKU: AM-ATS3000

Price: Request Quote



ZIMMER ATS 2000

The built-in safety, reliability, and convenience of the A.T.S.® 2000 Tourniquet frees you to focus on surgery. Dual displays that monitor actual cuff pressure, precise microprocessor controls, and one-hour battery backup help assure patient safety. User-set alarms warn of unacceptable or hazardous conditions.

[Download Brochure](#)

[Read More](#)

SKU: AM-ATS2000

Price: Request Quote





ZIMMER ATS 1500

The Zimmer® A.T.S.® 1500 is an automatic tourniquet which has been found useful in producing a bloodless operating field in a surgical procedure involving extremities. It is indicated to temporarily occlude blood flow in the patient's extremities during a surgical procedure on those extremities.

[Download Brochure](#)

[Read More](#)

SKU: AM-ATS1500

Price: Request Quote

INDEX

A

ATS Dual Red Hose 4

C

Cuff Cylindrical 34 Inch ATS Zimmer 4

Cuff Large 42 Inch ATS Zimmer 5

S

Stryker SmartPump Dual Channel Tourniquet 2

Z

Zimmer ATS 1500 7

Zimmer ATS 2000 6

Zimmer ATS 3000 6

Zimmer ATS 4000 Tourniquet 3

Zimmer ATS Dual Blue Hose 3

Zimmer ATS1200 Single Channel Tourniquet System 5