

TABLE OF CONTENTS

Vaporizers 2

 Vetland D20 Veterinary Isoflurane Vaporizer 2



VAPORIZERS

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)

VETLAND D20 VETERINARY ISOFLURANE VAPORIZER

The **Vetland D20 Veterinary Isoflurane Vaporizer** is a high-precision, agent-specific device designed for continuous-flow veterinary anesthesia systems. Calibrated exclusively for isoflurane, it offers a concentration range of 0-5%, with automatic temperature and flow compensation to ensure consistent delivery regardless of ambient conditions or varying gas flow rates. Its universal Selectatec-style mounting allows for easy integration with a variety of anesthesia machines, with options for funnel-fill or keyed-fill systems—both featuring a large locking cap to prevent spills and incorrect fills. Built for durability, the rugged, corrosion-resistant construction ensures lifelong calibration, backed by a five-year warranty. Compact and lightweight at approximately 5.5 pounds with dimensions of 8" H x 4" W x 6" D, the D20 provides a reliable, precise, and versatile anesthesia vaporizing solution for veterinary use.

[Read More](#)

SKU: AM-520-11

Price: Request Quote



INDEX

V

Vetland D20 Veterinary Isoflurane Vaporizer 2