

TABLE OF CONTENTS

Vaporizers 2

 Mindray Desflurane Vaporizer 2

 Mindray Sevoflurane Vaporizer with quick fill adapter 3

 Datex Ohmeda Tec 6 Plus Vaporizer 3

 Datex Ohmeda Tec 7 Sevoflurane Vaporizer 4

 Datex Ohmeda Tec 7 Isoflurane Vaporizer 4

 Datex Ohmeda Tec 5 Sevoflourane Vaporizer 5

 Drager 19.1 Sevoflurane Vaporizer 5

 Drager 19.1 Isoflurane Vaporizer 6

 Drager Vapor 2000 Sevoflurane Vaporizer 6

 Drager Vapor 2000 Isoflurane Vaporizer 7



VAPORIZERS

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



MINDRAY DESFLURANE VAPORIZER

Desflurane Vaporizer for Mindray A4 Anesthesia Machine

[Read More](#)

SKU: AM-040-001997-00



MINDRAY SEVOFLURANE VAPORIZER WITH QUICK FILL ADAPTER

Includes Quick Fill Bottle Adapter 0004-00-0100
For Mindray A4 Anesthesia Machine

[Read More](#)

SKU: AM-0992-00-0148

DATEX OHMEDA TEC 6 PLUS VAPORIZER



The Tec 6 Plus is an electronic vaporizers that heats Desflurane to maintain constant temperature and vapour pressure for consistent output. Tec 6 Plus incorporates an LED display indicating vapouriser status: no output, low agent, warm-up, operational, and alarm battery low.

[Read More](#)

SKU: AM-TEC6P-DES-C

DATEX OHMEDA TEC 7 SEVOFLURANE VAPORIZER



The Tec 7 Vaporizer is compatible with all current generation Datex-Ohmeda Anesthesia machines and some of the Drager models. Features Easy-fil, which is designed to simplify the filling process while minimizing leaks. Sevoflurane only. Other agents require a separate vaporizer.

[Read More](#)

SKU: AM-TEC7-SEVO

DATEX OHMEDA TEC 7 ISOFLURANE VAPORIZER



The Tec 7 Vaporizer is compatible with all current generation Datex-Ohmeda Anesthesia machines and some of the Drager models. Features Easy-fil, which is designed to simplify the filling process while minimizing leaks. Isoflurane only. Other agents require a separate vaporizer.

[Read More](#)

SKU: AM-TEC4-ISO

DATEX OHMEDA TEC 5 SEVOFLOURANE VAPORIZER



The Tec 5 Vaporizer is compatible with all current generation Datex-Ohmeda Anesthesia machines and some of the Drager models.

Features Key Fill or Funnel Fill, which is designed to simplify the filling process while minimizing leaks. Sevoflurane only. Other agents require a separate vaporizer

[Read More](#)

SKU: AM-TEC5-SEVO



DRAGER 19.1 SEVOFLURANE VAPORIZER

The vaporizers are designed for out of circuit use in continuous flow techniques of inhalation anaesthesia. Temperature and flow compensated so that their output remains relatively constant despite cooling due to vaporization and variations in inlet flow. Each vaporizer is clearly labeled with the name of the anaesthetic agent for which it is designed.

[Read More](#)

SKU: AM-DR-19.1-SEV



DRAGER 19.1 ISOFLURANE VAPORIZER

The vaporizers are designed for out of circuit use in continuous flow techniques of inhalation anaesthesia. Temperature and flow compensated so that their output remains relatively constant despite cooling due to vaporization and variations in inlet flow. Each vaporizer is clearly labeled with the name of the anaesthetic agent for which it is designed.

[Read More](#)

SKU: AM-DR-19.1-ISO

DRAGER VAPOR 2000 SEVOFLURANE VAPORIZER



The Vapor 2000 is a vaporizer for volatile anesthetic agents with unsurpassed performance. Originally developed for the extreme requirements of mobile army surgical units, these demanding ambient specifications benefit the normal, daily clinical routine in anesthesia.

[Read More](#)

SKU: AM-DR2000-SEVO

DRAGER VAPOR 2000 ISOFLURANE VAPORIZER



The Vapor 2000 is a vaporizer for volatile anesthetic agents with unsurpassed performance. Originally developed for the extreme requirements of mobile army surgical units, these demanding ambient specifications benefit the normal, daily clinical routine in anesthesia.

[Read More](#)

SKU: AM-DR2000-ISO

INDEX

D

Datex Ohmeda Tec 5 Sevoflourane Vaporizer 4

Datex Ohmeda Tec 6 Plus Vaporizer 3

Datex Ohmeda Tec 7 Isoflurane Vaporizer 4

Datex Ohmeda Tec 7 Sevoflurane Vaporizer 3

Drager 19.1 Isoflurane Vaporizer 5

Drager 19.1 Sevoflurane Vaporizer 5

Drager Vapor 2000 Isoflurane Vaporizer 6

Drager Vapor 2000 Sevoflurane Vaporizer 6

M

Mindray Desflurane Vaporizer 2

Mindray Sevoflurane Vaporizer with quick fill adapter 2