

TABLE OF CONTENTS

Ventilators	2
Oxygen Enrichment Accessories for Ventway Sparrow Portable Ventilator	2
Pole Mount for Ventway Sparrow Portable Ventilator	3
Carrying Case for Ventway Sparrow Portable Ventilator	3
Power Supplies and Batteries for Ventway Sparrow Portable Ventilator	4
Filters for Ventway Sparrow Portable Ventilator	5
Patient Breathing Circuits for Ventway Sparrow Portable Ventilator (Adult and Pediatric Options)	5
VOCSN Pro Roll Stand	5
Ventec Life Systems VOCSN V+Pro Transport Ventilator	6
Ventway Sparrow - EMS and Transport Turbine Ventilator	7
Philips Respironics V60 Plus	8
Precision Medical EasyGoVac Aspirator (PM65)	9
Philips Trilogy 202 Ventilator	9
Dräger Evita XL Ventilator	10
Puritan Bennett 560 Ventilator	11
FLOWMETER - OXYGEN, 0-15 LPM, USA, Ohmeda Quick Connect	12
Medtronic Puritan Bennett 980 Respiratory Ventilator	12
Maquet Servo-i Universal ventilator	12
Covidien Puritan Bennett 840 Respiratory Ventilator	13
Autovent 3000 Automatic Transport Ventilator, 2-Second Version	13
Autovent 3000 Automatic Transport Ventilator, 1.5-Second Version	14
Autovent 2000 Automatic Transport Ventilator, 2-Second Version	14
Autovent 2000 Automatic Transport Ventilator, 1.5-Second Version	15
ALLIED Healthcare AHP300 Ventilator	15
Bipap Vision Ventilator System	16
Respironics PLV-102 Ventilator	16



VENTILATORS

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



OXYGEN ENRICHMENT ACCESSORIES FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- Oxygen Mixer (OEM: VWSP-580)
- Low Pressure O2 Enrichment Kit (OEM: VWSP-304)
- DISS-DISS Oxygen Tube (green) (OEM: VWSP-590)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$40.00



POLE MOUNT FOR VENTWAY SPARROW PORTABLE VENTILATOR

The **Pole Mount (OEM: VWMNT-100) for the Ventway Sparrow Portable Ventilator** provides a secure, space-saving solution for easily attaching your Ventway Sparrow Portable Ventilator to standard IV poles or medical stands. Engineered for optimal stability and quick accessibility, this durable mount ensures your critical equipment is safely positioned during routine bedside care and active patient transport.

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: AM-VWMNT-100

Price: \$200.00



CARRYING CASE FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- Ventway Sparrow Soft Carrying Bag (OEM: VWSP-910)
- Ventway Sparrow Rigid Case (OEM: VWSP-920)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$245.00



POWER SUPPLIES AND BATTERIES FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- Rechargeable Battery (OEM: VWSP-420)
- Non-Rechargeable Battery (OEM: VWSP-430)
- AC/DC Power Adaptor (OEM: VWSP-141)
- Desktop Charger (OEM: VWSP-860)
- DC/DC Car Power Adaptor (OEM: VWSP-142)
- DC Cable - Without Input Terminal (OEM: VWSP-929)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$125.00



FILTERS FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- HME Filter (OEM: VWSP-305)
- Bacterial/Viral Filter (OEM: VWSP-306)
- Filter Pack (5 pack) (EOM: VWSP-810)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$8.00



PATIENT BREATHING CIRCUITS FOR VENTWAY SPARROW PORTABLE VENTILATOR (ADULT AND PEDIATRIC OPTIONS)

Sold Per Box. Adult quantities are 10 per box. Pediatric quantities are 15 per box.

Available Options:

- Adult 1.8m Patient Breathing Circuit (OEM: VWSP-1450F) - QTY: 10/Box
- Adult 2.4m Patient Breathing Circuit (OEM: VWSP-1370) - QTY: 10/Box
- Pediatric 1.8m Patient Breathing Circuit (OEM: VWSP-1550F) - QTY: 15/Box
- Pediatric 2.4m Patient Breathing Circuit (OEM: VWSP-1390) - QTY: 15/Box

[Read More](#)

SKU: N/A

Price: \$229.00





VOCSN PRO ROLL STAND

The **VOCSN Pro Roll Stand** is a purpose-built mobility solution designed exclusively for the VOCSN respiratory system, enabling safe and efficient transport of this comprehensive ventilator across hospitals, long-term care facilities, and home care environments. Crafted with a durable, low-center-of-gravity five-spoke base and lockable chrome casters with 3-inch hard rubber wheels, the stand provides excellent stability and effortless maneuverability on most floor surfaces, while ensuring it stays securely in place during patient care. It features a secure, tool-free docking system with a precise T-slot for stable attachment of the VOCSN ventilator (sold separately), along with a built-in utility basket and integrated handle for organization and control. Versatile mounting points allow attachment of accessory arms, IV stands, or feeding pumps, creating a customizable care station, while dedicated shelves and brackets securely store the AC power supply, heater/humidifier units, and oxygen cylinders—up to two—making it the ideal solution for portable, reliable respiratory support in various care settings.

OEM Part Number: VENPRT-00569-001 | 24 Month Warranty

[Ventec Life Systems VOCSN V+Pro Transport Ventilator Sold Separately \(AM-VENPRT-01198-000\)](#)

[Read More](#)

SKU: AM-VENPRT-00569-001

Price: \$795.00



VENTEC LIFE SYSTEMS VOCSN V+PRO TRANSPORT VENTILATOR

The **Ventec Life Systems VOCSN V+Pro Transport Ventilator** is a versatile, all-in-one respiratory support device designed for hospital, transport, and home environments. Combining five essential therapies—ventilation (both invasive and non-invasive), oxygen concentration (internal, equivalent to 6 L/min), cough assist (Touch Button Cough™), suction (hospital-grade), and high-performance nebulization—VOCSN offers comprehensive respiratory management in one compact, lightweight unit. Weighing approximately 70% less and smaller than five separate devices, it provides up to 9 hours of onboard battery life and features an intuitive 7-inch color LCD touchscreen for quick setup and easy operation. Built for safety and efficiency, VOCSN delivers uninterrupted therapy with seamless switching between modes, reducing the risk of circuit misconnects, and supports scaling to individual patient needs—making it ideal for adult and pediatric patients over 5 kg. Its advanced VOCSN Multi-View system enables complete trending and monitoring across multiple therapies, ensuring precise patient management and safety during transport and in critical care settings.

[Ventec VOCSN V+Pro Brochure](#)

OEM Part Number: VENPRT-01198-000 | Warranty for 2 years or 10,000 hours

[Optional Roll Stand Available \(AM-VENPRT-00569-001 sold separately\)](#)

[Read More](#)

SKU: AM-VENPRT-01198-000

Price: \$12,395.00



VENTWAY SPARROW - EMS AND TRANSPORT TURBINE VENTILATOR

The **Ventway Sparrow** is a compact, impact-resistant emergency and transport ventilator designed for field and combat medicine, enabling effective ventilation from the point of injury to the ICU with a single, portable device. It supports multiple ventilation modes, including SIMV-VC-PS, SIMV-PC, CPAP, HFNC, and CPR mode, and is suitable for adult and pediatric patients over 5 kg. Key performance specs include a wide respiratory rate range (1-60 bpm), peak inspiratory pressure up to 60 cm H₂O, tidal volumes from 50-2000 mL, adjustable inspiratory time and I:E ratios, high peak flow for both mandatory and spontaneous breaths, FiO₂ from 21% to 95%, and PEEP from 0 to 20 cm H₂O. The device is operable across altitudes from -1200 ft to 25,000 ft and supports battery power (rechargeable 2850 mAh and non-rechargeable 3000 mAh), plus external AC/DC and DC-DC power adapters. Each package includes battery(s), a filter, and an AC/DC power adapter. This portable, user-friendly ventilator is designed to deliver reliable respiratory support in austere environments where space and weight are at a premium.

[Ventway Sparrow Ventilator & Accessory Brochure](#)

[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: AM-VWSP-100

Price: \$9,350.00



PHILIPS RESPIRONICS V60 PLUS

The **Philips Respironics V60 Plus** is a microprocessor-controlled ventilator providing bi-level, positive air ventilation for both non-invasive (NIV) and invasive spontaneous breathing support. It offers two ventilation therapy options: Non-Invasive Ventilation (NIV) and High Flow Therapy (HFT). The V60 Plus includes advanced features such as AVAPS (Average Volume Assured Pressure Support) to monitor tidal volume in a pressure-limited mode, PCV (Pressure Control Ventilation) for patients with reduced lung capacity, and CPAP (Continuous Positive Airway Pressure) to provide continuous pressure throughout the respiratory cycle. Suitable for both pediatric and adult patients, it incorporates Philips auto-tracking technology for therapy adjustments without machine changes. The ventilator features a user-friendly touchscreen display and a 6-hour internal battery, enabling intra-hospital transport. Additional capabilities include high-resolution graphics for waveform interpretation, standby mode for patient-clinician interaction, data communication for hospital connectivity, and Respi-Link for remote diagnostics and upgrades. CPAP

with C-Flex offers flow-based expiratory pressure relief for improved comfort, while PCV provides greater control over the patient's breathing pattern. The package includes a battery and stand, with light assembly potentially required.

[Read More](#)

SKU: AM-V60PLUS



PRECISION MEDICAL EASYGOVAC ASPIRATOR (PM65)

The **Precision Medical EasyGoVac Aspirator (PM65)** is a portable, battery-powered suction unit designed for use in hospitals, nursing homes, ambulances, and homecare settings. It comes equipped with a hospital-grade power cord, carry bag, DC power cord, charger, and canister. The EasyGoVac delivers continuous suction power, producing vacuum up to 483 mmHg, and its high-performance battery provides over 2 hours of runtime. User-friendly features include a low battery light and an easy-to-clean, closed design that guards against bacterial growth. Safety features, such as an overflow prevention system, are incorporated into the pump to prevent contamination. The unit weighs 9.6 lbs (4.4 kg) and includes a replaceable inline hydrophobic bacteria filter and an 800cc disposable suction canister with a float shut-off.

[Read More](#)

SKU: PM65DCA



PHILIPS TRILOGY 202 VENTILATOR

The **Philips Trilogy 202 Ventilator** is an advanced, dual-mode device offering both volume-control and pressure-control ventilation for a range of patient needs, whether invasive or noninvasive. Its versatile breath delivery system and flexible setup options eliminate the need for frequent equipment changes, ensuring seamless and continuous patient care. One of its standout features is its ability to compensate for leaks in both pressure and volume control modes, which allows for the use of simpler passive circuits, potentially saving significant time and costs. The Trilogy 202 is also designed for intra-hospital transport, equipped with a 3-hour internal battery and an optional 3-hour detachable battery, and weighs only 5.6 kg, making it highly portable. Its integrated carrying handle adds to the convenience, ensuring safe and easy transport across hospital settings. For patients transitioning from hospital to home care, the Trilogy 202 offers a consistent user interface and identical ventilation modes and breathing circuit configurations, ensuring continuity of care. The ventilator also provides leak-compensated volume modes, enabling consistent tidal volume delivery even with non-invasive ventilation or cuff leaks. Setup is straightforward, with a color screen and intuitive menu navigation, and the integrated blender allows for oxygen administration from 21% to 100%. Additionally, the device supports two different treatment prescriptions, accommodating patients who may require different ventilation settings during the day and night. Detailed data storage on a secure digital (SD) card enables comprehensive review of treatment trends, breath-by-breath waveforms, and usage patterns through DirectView software, making the Trilogy 202 a robust and versatile tool for both acute and chronic care environments.

[Read More](#)

SKU: AM-PH-T202



DRÄGER EVITA XL VENTILATOR

The **Dräger Evita XL Ventilator** is a sophisticated and versatile piece of medical equipment designed to provide advanced respiratory support for critically ill patients. Integrating state-of-the-art technology, it offers a wide range of ventilation modes to accommodate various patient needs, from neonatal to adult care. The intuitive user interface and customizable settings enhance ease of use for healthcare professionals, enabling precise control over ventilation parameters. Additionally, the Evita XL is equipped with advanced monitoring capabilities, ensuring real-time analysis and optimization of patient respiratory status. Its robust design, coupled with reliable performance, makes the Dräger Evita XL an indispensable asset in intensive care units and other critical care environments.

[Read More](#)

SKU: AM-EVITA-XL



PURITAN BENNETT 560 VENTILATOR

The **Puritan Bennett™ 560 portable ventilator** empowers patients to live with greater independence while offering clinicians the reassurance of a reliable and trusted system. This compact and lightweight device, weighing just 4.5 kg, operates quietly at less than 30 dBA at 1 meter, making it suitable for both hospital and home environments. Its internal lithium-ion battery provides up to 11 hours of operation, with a real-time battery use time indicator displaying the remaining hours and minutes, ensuring continuous support. The ventilator's USB memory device facilitates seamless data transfer from home to clinical settings, enhancing patient care. Additionally, ventilator parameters and alarm settings are linked to mitigate the risk of inappropriate configurations. With oxygen enrichment capabilities using a low flow source and an integrated measurement kit to monitor FiO₂ levels, the Puritan Bennett™ 560 is a versatile and essential tool for patient respiratory support. Its dimensions of 23.5 cm W x 31.5 cm D x 15.4 cm H, combined with a robust power supply and USB compatibility, make it an exceptional choice for both patients and healthcare providers.

[Read More](#)

SKU: AM-PB-560



FLOWMETER - OXYGEN, 0-15 LPM, USA, OHMEDA QUICK CONNECT

[Read More](#)

SKU: AM-FMO-15U-OM

Price: \$49.95



MEDTRONIC PURITAN BENNETT 980 RESPIRATORY VENTILATOR

The Medtronic Puritan Bennett 980 ventilator system helps patients breathe more naturally through innovative breath delivery technology. The Puritan Bennett 980 has an intuitive and highly configurable touch screen display that gives clinicians the options to see as little or as much information as needed including five graph and loop layout options. On the digital communication board, the PB 980 offers multiple USB ports for connections to external drives and an HDMI output port to project the touchscreen externally.

[Read More](#)

SKU: AM-PB-980



MAQUET SERVO-I UNIVERSAL VENTILATOR

The Maquet SERVO i ventilation system can be used on a variety of patients from neonatal, pediatric, to adults. It has the ability to handle the most acute phases of respiratory distress, through recovery, and to the waning phases of the ventilation process. The SERVO-i continuously delivers outstanding ventilator treatment as gently as possible. This is in part to performance, monitoring capabilities, treatment options, and tools that this ventilator offers. Since this ventilator can be used on a wide variety of patients, it is suitable for most environments from the ICU to the NICU department. The Maquet SERVO-i ventilator allows you to choose the treatment option based on your patients' needs without having to worry about the ventilator's capabilities.

[Read More](#)

SKU: AM-SERVO-I



COVIDIEN PURITAN BENNETT 840 RESPIRATORY VENTILATOR

The refurbished, Covidien Puritan Bennett® 840™ Ventilator System, for purchase or rental, is the flagship product in our line of critical care ventilators. It is responsive to patients and offers superior comfort, delivering sensitive, precise breaths to critically ill neonatal through adult patients.

[Read More](#)

SKU: AM-PB-840

AUTOVENT 3000 AUTOMATIC TRANSPORT VENTILATOR, 2-SECOND VERSION



The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. second-version 3000

[Read More](#)

SKU: AM-ALLIED-AV-3000-2.0

AUTOVENT 3000 AUTOMATIC TRANSPORT VENTILATOR, 1.5-SECOND VERSION



In some instances, product photos may vary from actual product models, sizes and/or colors listed on the product detail page. The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. 1.5E

[Read More](#)

SKU: AM-ALLIED-AV-3000-1.5

AUTOVENT 2000 AUTOMATIC TRANSPORT VENTILATOR, 2-SECOND VERSION



The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. second-version 2000

[Read More](#)

SKU: AM-ALLIED-AV-2000-2.0

AUTOVENT 2000 AUTOMATIC TRANSPORT VENTILATOR, 1.5-SECOND VERSION



In some instances, product photos may vary from actual product models, sizes and/or colors listed on the product detail page. The LSP AutoVent series integrates field experience and engineering excellence, producing a lightweight, transportable automatic ventilator. 1.5 Second-Version

[Read More](#)

SKU: AM-ALLIED-AV-2000-1.5

ALLIED HEALTHCARE AHP300 VENTILATOR



The AHP300 is a new broad-featured transport ventilator designed to meet the specific needs of modern EMS response. It features the broad set of ventilation modes of critical care ventilators and a single panel user interface that is easy to use, easy to clean rugged and visible even in bright light conditions. It has quick start programmable presets which let the institution pre-program up to five ventilation protocols into memory and recall them instantly with the press of separate rapid recall buttons. This feature allows the responder to typically begin ventilating a patient with their desired protocol in less than 30 seconds. It also features innovative ventilation characteristics like Dynamic Flow Optimization and Seamless Mode Transition.

[Read More](#)

SKU: AM-ALLIED-AHP300



BIPAP VISION VENTILATOR SYSTEM

Integrated Display Screen displays easy-to-read real-time graphics in waveform or bar scale format, providing the following measured and calculated numerical values to assist with ventilator management:

[Read More](#)

SKU: AM-RESP-BIPAP-VISION

RESPIRONICS PLV-102 VENTILATOR



Respironics LifeCare PLV-102 Ventilator (with internal oxygen blender) - The LIFECARE PLV-102 is a microprocessor controlled, electrically powered portable volume ventilator that may also be pressure limited. Specifically designed for long-term use in the home, it is equally suitable for transport, short-term and institutional use. The ventilator is small and lightweight and will operate using a 12 volt DC external battery as well as 120 volt AC power. Reliability, versatility, and ease of operation are all hallmarks of this machine, making it the ideal choice for this setting.

[Download Brochure](#)

[Read More](#)

SKU: AM-PLV-102

INDEX

A

ALLIED Healthcare AHP300 Ventilator 15
Autovent 2000 Automatic Transport Ventilator, 1.5-Second Version 14
Autovent 2000 Automatic Transport Ventilator, 2-Second Version 14
Autovent 3000 Automatic Transport Ventilator, 1.5-Second Version 13
Autovent 3000 Automatic Transport Ventilator, 2-Second Version 13

B

Bipap Vision Ventilator System 15

C

Carrying Case for Ventway Sparrow Portable Ventilator 3
Covidien Puritan Bennett 840 Respiratory Ventilator 12

D

Dräger Evita XL Ventilator 10

F

Filters for Ventway Sparrow Portable Ventilator 4
FLOWMETER - OXYGEN, 0-15 LPM, USA, Ohmeda Quick Connect 11

M

Maquet Servo-i Universal ventilator 12
Medtronic Puritan Bennett 980 Respiratory Ventilator 11

O

Oxygen Enrichment Accessories for Ventway Sparrow Portable Ventilator 2

P

Patient Breathing Circuits for Ventway Sparrow Portable Ventilator (Adult and Pediatric Options) 4
Philips Respironics V60 Plus 7
Philips Trilogy 202 Ventilator 9
Pole Mount for Ventway Sparrow Portable Ventilator 2
Power Supplies and Batteries for Ventway Sparrow Portable Ventilator 3
Precision Medical EasyGoVac Aspirator (PM65) 8
Puritan Bennett 560 Ventilator 10

R

Respironics PLV-102 Ventilator 16

V

Ventec Life Systems VOCSN V+Pro Transport Ventilator 6
Ventway Sparrow - EMS and Transport Turbine Ventilator 7
VOCSN Pro Roll Stand 5