

TABLE OF CONTENTS

- Ventilators** 2
 - Hallowell 2024 PRO Veterinary Anesthesia Ventilator 2
 - Oxygen Enrichment Accessories for Ventway Sparrow Portable Ventilator
..... 3
 - Pole Mount for Ventway Sparrow Portable Ventilator 3
 - Carrying Case for Ventway Sparrow Portable Ventilator 4
 - Power Supplies and Batteries for Ventway Sparrow Portable Ventilator
..... 4
 - Filters for Ventway Sparrow Portable Ventilator 5
 - Patient Breathing Circuits for Ventway Sparrow Portable Ventilator (Adult
and Pediatric Options) 5
 - VOCSN Pro Roll Stand 6
 - Ventec Life Systems VOCSN V+Pro Transport Ventilator 7
 - Ventway Sparrow - EMS and Transport Turbine Ventilator 8
 - Precision Medical EasyGoVac Aspirator (PM65) 8
 - FLOWMETER - OXYGEN, 0-15 LPM, USA, Ohmeda Quick Connect 9



VENTILATORS

Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)



HALLOWELL 2024 PRO VETERINARY ANESTHESIA VENTILATOR

The **Hallowell 2024 PRO Veterinary Anesthesia Ventilator** is a compact, reliable, and versatile pressure-limited, time-cycled ventilator specifically designed for small animal anesthesia. It delivers precise, controlled mechanical ventilation with adjustable tidal volume, breath rate, and inspiratory pressure relief, ensuring safe support across a wide range of small animal cases. Features like manual breath control, audible alarms for overpressure or circuit issues, and a visible breath indicator enhance clinician confidence and patient safety. Its compatibility with most anesthesia circuits, along with a durable, space-saving design, makes the 2024 PRO an ideal solution for general, specialty, emergency, and teaching veterinary practices seeking dependable ventilation tailored to small animal patients.

[Read More](#)

SKU: 000A6758

Price: Request Quote



OXYGEN ENRICHMENT ACCESSORIES FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- Oxygen Mixer (OEM: VWSP-580)
- Low Pressure O2 Enrichment Kit (OEM: VWSP-304)
- DISS-DISS Oxygen Tube (green) (OEM: VWSP-590)

[Ventway Sparrow Ventilator & Accessory Brochure](#)

[Ventway Sparrow Ventilator Sold Separately Here](#)

[View all Available Ventway Sparrow Accessories Here](#)

[Read More](#)

SKU: N/A

Price: \$40.00 - \$525.00



POLE MOUNT FOR VENTWAY SPARROW PORTABLE VENTILATOR

The **Pole Mount (OEM: VWMNT-100) for the Ventway Sparrow Portable Ventilator** provides a secure, space-saving solution for easily attaching your Ventway Sparrow Portable Ventilator to standard IV poles or medical stands. Engineered for optimal stability and quick accessibility, this durable mount ensures your critical equipment is safely positioned during routine bedside care and active patient transport.

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: AM-VWMNT-100

Price: \$200.00



CARRYING CASE FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- Ventway Sparrow Soft Carrying Bag (OEM: VWSP-910)
- Ventway Sparrow Rigid Case (OEM: VWSP-920)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$245.00 - \$450.00





POWER SUPPLIES AND BATTERIES FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- Rechargeable Battery (OEM: VWSP-420)
- Non-Rechargeable Battery (OEM: VWSP-430)
- AC/DC Power Adaptor (OEM: VWSP-141)
- Desktop Charger (OEM: VWSP-860)
- DC/DC Car Power Adaptor (OEM: VWSP-142)
- DC Cable - Without Input Terminal (OEM: VWSP-929)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$125.00 - \$480.00



FILTERS FOR VENTWAY SPARROW PORTABLE VENTILATOR

Available Options:

- HME Filter (OEM: VWSP-305)
- Bacterial/Viral Filter (OEM: VWSP-306)
- Filter Pack (5 pack) (EOM: VWSP-810)

[Ventway Sparrow Ventilator & Accessory Brochure](#)
[Ventway Sparrow Ventilator Sold Separately Here](#)
[View all Available Ventway Sparrow Accessories Here](#)
[Read More](#)

SKU: N/A

Price: \$8.00 - \$40.00



PATIENT BREATHING CIRCUITS FOR VENTWAY SPARROW PORTABLE VENTILATOR (ADULT AND PEDIATRIC OPTIONS)

Sold Per Box. Adult quantities are 10 per box. Pediatric quantities are 15 per box.

Available Options:

- Adult 1.8m Patient Breathing Circuit (OEM: VWSP-1450F) - QTY: 10/Box
- Adult 2.4m Patient Breathing Circuit (OEM: VWSP-1370) - QTY: 10/Box
- Pediatric 1.8m Patient Breathing Circuit (OEM: VWSP-1550F) - QTY: 15/Box
- Pediatric 2.4m Patient Breathing Circuit (OEM: VWSP-1390) - QTY: 15/Box

[Read More](#)

SKU: N/A

Price: \$229.00 - \$360.00





VOCSN PRO ROLL STAND

The **VOCSN Pro Roll Stand** is a purpose-built mobility solution designed exclusively for the VOCSN respiratory system, enabling safe and efficient transport of this comprehensive ventilator across hospitals, long-term care facilities, and home care environments. Crafted with a durable, low-center-of-gravity five-spoke base and lockable chrome casters with 3-inch hard rubber wheels, the stand provides excellent stability and effortless maneuverability on most floor surfaces, while ensuring it stays securely in place during patient care. It features a secure, tool-free docking system with a precise T-slot for stable attachment of the VOCSN ventilator (sold separately), along with a built-in utility basket and integrated handle for organization and control. Versatile mounting points allow attachment of accessory arms, IV stands, or feeding pumps, creating a customizable care station, while dedicated shelves and brackets securely store the AC power supply, heater/humidifier units, and oxygen cylinders—up to two—making it the ideal solution for portable, reliable respiratory support in various care settings.

OEM Part Number: VENPRT-00569-001 | 24 Month Warranty

[Ventec Life Systems VOCSN V+Pro Transport Ventilator Sold Separately \(AM-VENPRT-01198-000\)](#)

[Read More](#)

SKU: AM-VENPRT-00569-001

Price: \$795.00



VENTEC LIFE SYSTEMS VOCSN V+PRO TRANSPORT VENTILATOR

The **Ventec Life Systems VOCSN V+Pro Transport Ventilator** is a versatile, all-in-one respiratory support device designed for hospital, transport, and home environments. Combining five essential therapies—ventilation (both invasive and non-invasive), oxygen concentration (internal, equivalent to 6 L/min), cough assist (Touch Button Cough™), suction (hospital-grade), and high-performance nebulization—VOCSN offers comprehensive respiratory management in one compact, lightweight unit. Weighing approximately 70% less and smaller than five separate devices, it provides up to 9 hours of onboard battery life and features an intuitive 7-inch color LCD touchscreen for quick setup and easy operation. Built for safety and efficiency, VOCSN delivers uninterrupted therapy with seamless switching between modes, reducing the risk of circuit misconnects, and supports scaling to individual patient needs—making it ideal for adult and pediatric patients over 5 kg.

Its advanced VOCSN Multi-View system enables complete trending and monitoring across multiple therapies, ensuring precise patient management and safety during transport and in critical care settings.

[Ventec VOCN V+Pro Brochure](#)

OEM Part Number: VENPRT-01198-000 | Warranty for 2 years or 10,000 hours

[Optional Roll Stand Available \(AM-VENPRT-00569-001 sold separately\)](#)

[Read More](#)

SKU: AM-VENPRT-01198-000

Price: \$12,395.00



VENTWAY SPARROW - EMS AND TRANSPORT TURBINE VENTILATOR

The **Ventway Sparrow** is a compact, impact-resistant emergency and transport ventilator designed for field and combat medicine, enabling effective ventilation from the point of injury to the ICU with a single, portable device. It supports multiple ventilation modes, including SIMV-VC-PS, SIMV-PC, CPAP, HFNC, and CPR mode, and is suitable for adult and pediatric patients over 5 kg. Key performance specs include a wide respiratory rate range (1-60 bpm), peak inspiratory pressure up to 60 cm H₂O, tidal volumes from 50-2000 mL, adjustable inspiratory time and I:E ratios, high peak flow for both mandatory and spontaneous breaths, FiO₂ from 21% to 95%, and PEEP from 0 to 20 cm H₂O. The device is operable across altitudes from -1200 ft to 25,000 ft and supports battery power (rechargeable 2850 mAh and non-rechargeable 3000 mAh), plus external AC/DC and DC-DC power adapters. Each package includes battery(s), a filter, and an AC/DC power adapter. This portable, user-friendly ventilator is designed to deliver reliable respiratory support in austere environments where space and weight are at a premium.

[Ventway Sparrow Ventilator & Accessory Brochure](#)

[View all Available Ventway Sparrow Accessories Here](#)

[Read More](#)

SKU: AM-VWSP-100

Price: \$9,350.00



PRECISION MEDICAL EASYGOVAC ASPIRATOR (PM65)

The **Precision Medical EasyGoVac Aspirator (PM65)** is a portable, battery-powered suction unit designed for use in hospitals, nursing homes, ambulances, and homecare settings. It comes equipped with a hospital-grade power cord, carry bag, DC power cord, charger, and canister. The EasyGoVac delivers continuous suction power, producing vacuum up to 483 mmHg, and its high-performance battery provides over 2 hours of runtime. User-friendly features include a low battery light and an easy-to-clean, closed design that guards against bacterial growth. Safety features, such as an overflow prevention system, are incorporated into the pump to prevent contamination. The unit weighs 9.6 lbs (4.4 kg) and includes a replaceable inline hydrophobic bacteria filter and an 800cc disposable suction canister with a float shut-off.

[Read More](#)

SKU: PM65DCA

Price: Request Quote



FLOWMETER - OXYGEN, 0-15 LPM, USA, OHMEDA QUICK CONNECT

[Read More](#)

SKU: AM-FMO-15U-OM

Price: \$49.95

INDEX

C

Carrying Case for Ventway Sparrow Portable Ventilator 3

F

Filters for Ventway Sparrow Portable Ventilator 4

FLOWMETER - OXYGEN, 0-15 LPM, USA, Ohmeda Quick Connect 8

H

Hallowell 2024 PRO Veterinary Anesthesia Ventilator 2

O

Oxygen Enrichment Accessories for Ventway Sparrow Portable Ventilator 2

P

Patient Breathing Circuits for Ventway Sparrow Portable Ventilator (Adult and Pediatric Options) 5

Pole Mount for Ventway Sparrow Portable Ventilator 3

Power Supplies and Batteries for Ventway Sparrow Portable Ventilator 4

Precision Medical EasyGoVac Aspirator (PM65) 8

V

Ventec Life Systems VOCSN V+Pro Transport Ventilator 6

Ventway Sparrow - EMS and Transport Turbine Ventilator 7

VOCSN Pro Roll Stand 5