

# TABLE OF CONTENTS

<b>Warming Cabinets</b> .....	3
Enthermics EC230L Fluid Warmer .....	4
Steris Dual Compartment Blanket Warmer QDJ02 .....	4
Steris DJ02 Dual Compartment Blanket Warmer .....	5
MAC Medical SWC183024 Glass Door Warming Cabinet .....	6
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Left Hand Door Chrome Shelves .....	6
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Left Hand Door SS Shelves .....	7
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Right Hand Door Chrome Shelves .....	8
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Right Hand Door SS Shelves .....	9
Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Glass Door, 4" Base, Dual Chamber .....	10
Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Stainless Steel Door, 4" Base, Dual Chamber .....	11
Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Glass Steel Door, 4" Base, Single Chamber .....	12
Steris Dual Compartment Blanket Warmer DJ06 .....	13
Steris Tabletop Blanket Warmer Model DJ01 .....	13
Steris Tabletop Blanket Warmer - Model DJ03 .....	14
Single Chamber Blanket/Fluid Warmer Glass Door .....	15
MAC Medical Warming Cabinet (New) .....	15
Steris Dual Compartment Blanket Warmer QDJ05 .....	16
750 Bair Hugger without Cart .....	17
Arizant 775 Bair Hugger .....	17
MAC Medical Warmer, D-Series, 26.5inD X 30inW X 24.5inH .....	17
MAC Medical Warmer, D-Series, 26.5inD X 24inW X 24.5inH .....	18
Accucold 24" Wide Blanket Warmer .....	18
MAC Medical Warming Cabinet .....	19
Mac Medical Warmer Stainless Steel Doors .....	19
Steris DJ04 Dual Compartment Blanket Warmer .....	20
Single compartment heat table top warming cabinet .....	20
MAC Medical Warming Cabinet (Refurbished) .....	21

Blickmann Digital Tabletop Glass Front Warmer .....	21
GE Ohmeda 4400 Infant Warmer System .....	22
Blickmann 7929MG Tabletop Glass Door Blanket Warmer .....	22
Single Chamber Blanket/Fluid Warmer .....	23
Olympic Warmette 56920 Blanket Warmer .....	23
Blickmann Digital Full-size Warmer .....	24
Blickman 7927TG Warming Cabinet .....	24
Blickmann Digital Tabletop Warmer .....	25
Pedigo 2120 Fluid Warmer Refurbished with Caster Base .....	25
Steris table top blanket warmer model QDJ03 .....	26
Amsco QDJ01 Tabletop Warming cabinet .....	26
Steris QDJ04 Dual Compartment Blanket Warmer .....	27
CMP Blanket and Solution Warming Cabinet .....	28
Blickman Warming Cabinet models 7921TG and 7921T .....	28
Steris AMSCO QDJ06 Dual Compartment Warming Cabinets .....	28
Steris AMSCO QDJ05 Single Compartment Warming Cabinet .....	29
Skytron J2 Series Warming Cabinets .....	29
MAC Medical Warming Cabinets .....	30
Steris Analog Tabletop Warmer .....	31



## WARMING CABINETS

**Auxo Medical refurbishes and distributes warming cabinets to fit all medical practice or laboratory needs. With a focus on consistent, accurate, and efficient heating systems, Auxo Medical's inventory of medical warming cabinets is ready to serve your facility.**

Designed to store and heat sterile intravenous fluids, surgical irrigation fluids, linens, blankets, and more. Auxo Medical's selection of warming cabinets come from the most respected brands in the industry. Additionally, with features like open shelf designs, LED readouts, and over temperature alarms. Ultimately, our warming cabinets offer superior durability and quality. Some of our medical warming cabinet units also offer adjustable shelves for optimal storage of your rotating inventory. Also, we carry non-heated storage compartments for any other essential equipment.

**Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:**

[\*\*BRAND NEW\*\*](#)

[\*\*RECONDITIONED\*\*](#)

**REFURBISHED**



## ENTHERMICS EC230L FLUID WARMER

The **Fully Refurbished Entermics EC230L Fluid Warmer** is a reliable and efficient countertop device designed to keep IV fluids warmed, supporting patient recovery and reducing surgical site infection risks. Capable of holding up to 11 one-liter bottles or 16 one-liter bags, it provides consistent warming in a compact footprint measuring 20.7 inches high, 16.2 inches wide, and 28 inches deep. Constructed with heavy-duty, durable cabinets built for long-term use, the warmer features locked handles for tamper-proof security and an insulated glass door that allows staff to quickly assess fluid levels at a glance. Its stackable design enables multiple units to be easily combined, making it an ideal solution for busy medical environments seeking dependable, straightforward fluid warming.

[Read More](#)

**SKU:** AM-EC230L



## STERIS DUAL COMPARTMENT BLANKET WARMER QDJ02

The **Steris Model QDJ02 Dual Compartment Blanket Warmer** is a completely remanufactured and biomedically inspected unit designed to provide reliable, precise warming for blankets and IVs. It features two independent compartments with digital controls, allowing users to set temperatures within the range of 90 to 160°F for optimal patient comfort and safety. The upper compartment measures 23.5" W x 16" D x 13" H, and the lower compartment is 23.5" W x 16" D x 38" H, all housed within a modern stainless steel chassis. Enhanced with a new circulation fan, digital displays, and overheat indicators, this contemporary warmer offers durable, efficient performance suitable for medical facilities seeking a high-quality, dependable warming solution.

[Read More](#)

**SKU:** AM-QDJ02



## STERIS DJ02 DUAL COMPARTMENT BLANKET WARMER

The **Steris DJ02 Dual Compartment Blanket Warmer** is a fully remanufactured and biomedically inspected unit designed for reliable warming of blankets and IV fluids in clinical settings. It features two independent compartments, each with adjustable digital temperature controls—allowing precise regulation of warmers for both IV and blanket modes. The upper compartment measures 23.5" W x 18" D x 13" H, and the lower compartment is the same width and depth but taller at 38" high, providing ample space for patient warming needs. The overall unit dimensions are 30" W x 22" D x 75.5" H, housed in stainless steel for durability and easy cleaning. Key features include a new circulation fan for improved airflow, overheat indicators, and contemporary design. Temperature ranges are set to 90-110°F for IV mode and 90-160°F for blanket mode—regulating safety and effectiveness. Powered by 120V, 50/60Hz, 1-phase power, it draws 14 amps (1680 watts), making it suitable for standard clinical power supplies.

[Read More](#)

**SKU:** AM-DJ02



## MAC MEDICAL SWC183024 GLASS DOOR WARMING CABINET

The **MAC Medical SWC183024 Glass Door Warming Cabinet** combines clear visibility with dependable warming performance in a compact footprint. Measuring 30 inches wide, 20.5 inches deep, and 24.5 inches high, this cabinet is designed to fit constrained spaces while keeping contents easily accessible. The glass door offers quick visual confirmation of stored items, while a built-in data logger captures environmental data for compliant recordkeeping. The USB-style data logger is included, enabling straightforward data retrieval and offline analysis. Built for reliability in clinical settings, the cabinet provides consistent warming to help maintain instrument readiness and workflow efficiency.

[Read More](#)

**SKU:** AM-SWC183024-G-NB



## ACCUCOLD 24"W X 71"T FULL-SIZE WIDE BLANKET WARMER LEFT HAND DOOR CHROME SHELVES

The **Accucold 24"W x 71"T Full-Size Wide Blanket Warmer** is a specialized appliance designed to maintain blankets at optimal warmth, ranging from 90 to 140°F, ideal for enhancing patient comfort in hospitals, medical and physical therapy offices, and spas. This unit employs a passive heating system, utilizing fixed heaters at the bottom to minimize lint buildup and ensure efficient, safe operation. It features a digital thermostat with an electronic LED readout for easy and precise temperature management. The warmer operates quietly, and includes a double-pane tempered glass door with stainless steel trim, a self-closing door, and a stainless steel interior with six adjustable **chrome shelves**. It also features a **left-hand door swing**, a factory-installed lock, and a durable 304-grade stainless steel exterior.

[Read More](#)

**SKU:** PTHC155GLHD





## ACCUCOLD 24"W X 71"T FULL-SIZE WIDE BLANKET WARMER LEFT HAND DOOR SS SHELVES

The **Accucold 24"W x 71"T Full-Size Wide Blanket Warmer** is expertly designed to provide consistent warmth to blankets, with a temperature range of 90 to 140°F, making it suitable for hospitals, medical and physical therapy offices, and spas. Its passive heating system, which utilizes fixed heaters, minimizes lint buildup and ensures efficient operation. The unit features a digital thermostat with an electronic LED readout for easy temperature monitoring. Operating silently, it includes a double-pane tempered glass door with a stainless steel trim, a self-closing door, a stainless steel interior, and six adjustable **stainless steel shelves**. It is equipped with a **left-hand door swing**, a factory-installed lock, and a durable 304-grade stainless steel exterior.

[Read More](#)

**SKU:** PTHC155GCSSLHD





## ACCUCOLD 24"W X 71"T FULL-SIZE WIDE BLANKET WARMER RIGHT HAND DOOR CHROME SHELVES

The **Accucold 24"W x 71"T Full-Size Wide Blanket Warmer** is designed to keep blankets warm within a temperature range of 90 to 140°F, suitable for hospitals, medical and physical therapy offices, and spas. It utilizes a passive heating system with fixed heaters to minimize lint buildup and ensure efficient operation. The unit features a digital thermostat and an electronic LED readout for easy temperature monitoring. Its silent operation is achieved through the absence of moving parts. The warmer includes a double-pane tempered glass door with a stainless steel trim, a self-closing door, and a stainless steel handle. The interior is made of stainless steel and features six adjustable **chrome shelves**. Additional features include a factory-installed lock, a **right-hand door swing** (non-reversible), a flush back, a foamed-in-place cabinet, and a 12.17 cu.ft. capacity. It also has a hospital-grade cord with a 'green dot' plug and magnetic door gaskets. The exterior is made of 304-grade stainless steel. An optional dolly with locking casters is available. The unit's dimensions are approximately 24"W x 71"T x 24"D.

[Read More](#)

**SKU:** PTHC155GRHD





## ACCUCOLD 24"W X 71"T FULL-SIZE WIDE BLANKET WARMER RIGHT HAND DOOR SS SHELVES

The **Accucold Full-Size Wide Blanket Warmer** is designed to maintain blankets at a consistent temperature between 90 to 140°F for hospitals, medical offices, spas, and physical therapy facilities. It utilizes passive heating with fixed heaters at the bottom to minimize lint buildup. The unit features a digital thermostat with an electronic LED readout for easy temperature monitoring. It operates silently due to the absence of moving parts. The warmer has a double-pane tempered glass door with a stainless steel trim and a self-closing mechanism. It also includes a stainless steel interior, six adjustable **stainless steel shelves**, a factory-installed lock, and a **right-hand door swing**. The exterior is constructed from 304-grade stainless steel. It has a 12.17 cu.ft. capacity within a 24" wide footprint. Additional features include a hospital-grade cord with a 'green dot' plug, magnetic door gaskets, a flush back, and a foamed-in-place cabinet. An optional dolly with locking casters is available for mobility. The unit's dimensions are approximately 24"W x 71"T x 24"D.

[Read More](#)

**SKU:** PTHC155GCSSRHD





## **MAC MEDICAL WARMER, D-SERIES, 26.5"D X 30"W X 74.5"H, GLASS DOOR, 4" BASE, DUAL CHAMBER**

The **Mac Medical Warmer, D-Series**, is a **dual-chamber warming cabinet with a glass door**, designed for use in medical environments. Its overall dimensions are 26.5"D x 30"W x 74.5"H, providing a total of 17.64 cubic feet of storage space across 4 adjustable shelves. The upper chamber's internal dimensions are 23"D x 26"W x 25"H, offering 8.65 cubic feet of storage, while the lower chamber measures 23"D x 26"W x 26"H, providing 8.99 cubic feet of storage. Both chambers maintain a temperature range of 90°-160°F, with a tolerance of +/- 2°F. The unit is built with a 4" base and operates at 120V, 12.5 Amp, 50/60 Hz, consuming an average of 0.90 kWh and producing 3071 BTU/hr.

[Read More](#)

**SKU:** DWC243074E-G-4B



## MAC MEDICAL WARMER, D-SERIES, 26.5" D X 30" W X 74.5" H, STAINLESS STEEL DOOR, 4" BASE, DUAL CHAMBER

The **Mac Medical Warmer, D-Series**, is a dual-chamber warming cabinet designed for medical facilities. It has overall dimensions of 26.5"D x 30"W x 74.5"H and offers 17.64 cubic feet of storage space with 4 adjustable shelves. The upper chamber has internal dimensions of 23"D x 26"W x 25"H, providing 8.65 cubic feet of storage, while the lower chamber measures 23"D x 26"W x 26"H, offering 8.99 cubic feet of storage. Both chambers maintain a temperature range of 90°-160°F with a tolerance of +/- 2°F. The unit features a **stainless steel door** and a 4" base. Electrically, it operates at 120V, 12.5 Amp, 50/60 Hz, with an average consumption of 0.90 kWh and 3071 BTU/hr.

[Read More](#)

**SKU:** DWC243074E-4B



## MAC MEDICAL WARMER, D-SERIES, 26.5" D X 30" W X 74.5" H, GLASS STEEL DOOR, 4" BASE, SINGLE CHAMBER

The **Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Glass Steel Door, 4" Base, Single Chamber**, is a robust and efficient solution designed for healthcare settings where maintaining precise temperatures for medical supplies is crucial. With an overall size of 26.5"D x 30"W x 74.5"H, this warmer provides ample space through its two adjustable shelves and one fixed shelf. The interior chamber, measuring 23"D x 26"W x 61"H, offers 21.1 cubic feet of storage, ensuring sufficient room for a variety of medical supplies. The temperature range of 90°-160° F, with a tolerance of +/- 3° F, guarantees that blankets and solution bottles are kept at optimal conditions. Operating on 120V with an average energy consumption of .47 kWh and 1604 BTU/hr, this warmer is both powerful and energy-efficient. With an approximate capacity of 80 blankets or 120 one-liter solution bottles, the Mac Medical Warmer, D-Series, is an indispensable asset in maintaining readiness and efficiency in medical environments.

[Read More](#)

**SKU:** AM-SWC243074-4G

## STERIS DUAL COMPARTMENT BLANKET WARMER DJ06

The **Steris Dual Compartment Blanket Warmer DJ06** is a versatile and robust solution for healthcare facilities, offering a substantial 15.8 cubic feet capacity to accommodate a variety of warming needs. Its LED readout provides instant temperature confirmation, ensuring that the internal environment is consistently monitored for optimal performance. Constructed with superior durable materials, this blanket warmer is designed to withstand years of rigorous use. Adjustable shelving allows for flexible organization and storage, making it easy to adapt to varying inventory demands. The unit's reliable heating system ensures uniform and efficient heat transfer, crucial for maintaining the quality of warmed items. Options for visual access to inventory include dual-pane tempered glass or traditional solid steel doors, catering to different preferences and security needs. Designed to meet AORN Recommended Practices for a Safe Environment of Care, the DJ06 incorporates technology to warm IV fluids, aligning with RP VIII.b.3 standards. For added security, a temperature lock-out feature provides consistency and control, while the ability to select Celsius or Fahrenheit enhances convenience and accuracy. This versatile unit can be installed as a free-standing or recessed fixture, offering flexibility in placement. With a 24" (610 mm) depth and options for steel or dual-paned tempered glass doors, the Steris DJ06 is a comprehensive and reliable choice for any medical facility.

[Read More](#)

**SKU:** AM-DJ06





## STERIS TABLETOP BLANKET WARMER MODEL DJ01

The **Steris Tabletop Blanket Warmer Model DJ01** is an essential device for medical facilities seeking to enhance patient comfort and care. This model, fully refurbished and biomedically inspected, ensures reliability and safety in its operations. It features a user-friendly design with an easy-grip door handle and a straightforward power on/off rocker switch, making it accessible for all users. The temperature settings can be easily adjusted using a dial, and the current temperature is displayed on a clear LCD window. For added safety, there is an overtemperature warning light that alerts staff if the internal temperature exceeds safe limits. This unit includes a spacious 16 x 23-inch shelf, accommodating various blanket sizes within its 13-1/2 x 24 x 16-inch inner dimensions. Its compact outer dimensions of 29.25 inches wide and 21.5 inches deep allow it to fit efficiently in different spaces. Additionally, the 4.5-foot heavy-duty power cord offers flexibility in placement, making the Steris Tabletop Blanket Warmer Model DJ01 a versatile and crucial addition to any medical environment.

**Inner Dimensions: 13-1/2 x 24 x 16"**

**Outer Dimensions: 29.25"W x 21.5"Deep**

[Read More](#)

**SKU:** AM-DJ01



## STERIS TABLETOP BLANKET WARMER - MODEL DJ03

The **Steris Tabletop Blanket Warmer, Model DJ03**, is engineered for convenience and reliability, ensuring that medical blankets are warmed to an ideal temperature for patient comfort. It features an easy-grip door handle, allowing for effortless access. The user-friendly design includes a temperature set dial marked in Fahrenheit, complemented by an LCD temperature window for precise temperature monitoring. An overtemperature warning light ensures safety by alerting users if the unit exceeds the set temperature. This model is designed for durability and ease of use, with a 4.5-foot heavy-duty power cord and a 21 x 23-inch shelf included for ample storage. Operation is simplified with a power on/off rocker switch. Each unit comes fully refurbished and has been biomedically inspected to meet rigorous standards. The inner dimensions of 13-1/2 x 24 x 22-7/8 inches and outer dimensions of 29.25"W x 26.5"Deep make it a compact yet spacious solution for healthcare settings seeking to provide patients with the comfort of warm blankets.

**Inner Dimensions: 13-1/2 x 24 x 22-7/8"**

**Outer Dimensions:** 29.25"W x 26.5"Deep

[Read More](#)

**SKU:** AM-DJ03



## **SINGLE CHAMBER BLANKET/FLUID WARMER GLASS DOOR**

The **Mac Medical Blanket/Fluid Warmer** is a state-of-the-art single chamber unit designed to efficiently store and warm blankets and fluids, ensuring they are at the optimal temperature for patient care. With an overall size of 30 inches in width, 26.6 inches in depth, and 36 inches in height, it is compact enough for easy placement within medical facilities while offering a generous chamber size of 26 inches wide, 23 inches deep, and 25 inches high. This unit features two adjustable shelves, providing flexibility to accommodate various sizes of blankets and fluids, and boasts a substantial storage capacity of 8.63 cubic feet. The temperature range of this warmer can be set between 90 to 160 degrees Fahrenheit, allowing for precise temperature control. Additionally, its glass door enables quick inventory checks without the need to open the unit, minimizing energy loss and contamination risk. This combination of features makes the Mac Medical Blanket/Fluid Warmer an essential piece of equipment for hospitals and clinics aiming to enhance patient comfort and care.

**Overall size: 30"W x 26.6"D x 36"H**

[Read More](#)

**SKU:** SWC243036-G-2B

## MAC MEDICAL WARMING CABINET (NEW)



### New Mac Medical 18" x 24" Warming Cabinet

The MAC Medical Warming Cabinet is perfect for keeping IV, surgical irrigation fluids, linens, and blankets warm. The warming cabinet provides independent, digitally controlled heating chambers that offer actual temperature and setpoint displays. Data is recorded every 30 minutes and stored for up to two years

[Read More](#)

**SKU:** AM-SWC182424-G-NB



## STERIS DUAL COMPARTMENT BLANKET WARMER QDJ05

The **Steris Dual Compartment Blanket Warmer QDJ05** represents a pinnacle of innovation and efficiency in medical equipment, designed to enhance patient care and comfort. This state-of-the-art device features two separate compartments, allowing for the simultaneous warming of blankets at different temperatures to meet varied clinical needs. The precise temperature control ensures blankets are warmed to the ideal temperature, improving patient comfort and preventing hypothermia during surgical procedures. Its sleek, space-saving design makes it a perfect fit for any healthcare setting, from hospitals to private clinics, ensuring that warm blankets are always readily available. The Steris QDJ05 is not only a testament to the advancements in medical technology but also a crucial tool in enhancing patient care and recovery.

[Read More](#)

**SKU:** AM-QDJ05-FS



## 750 BAIR HUGGER WITHOUT CART

3M Bair Hugger Model 750 Warming Unit

[Read More](#)

**SKU:** AM-750B



## ARIZANT 775 BAIR HUGGER

Refurbished 775 Patient Warmer  
3M Bair Hugger Model 775 Warming Unit Features:

[Read More](#)

**SKU:** AM-775B



## MAC MEDICAL WARMER, D-SERIES, 26.5IND X 30INW X 24.5INH

Download Brochure

[Read More](#)

**SKU:** SWC243024-NB



## MAC MEDICAL WARMER, D-SERIES, 26.5IND X 24INW X 24.5INH

Download Brochure

[Read More](#)

**SKU:** SWC242424-NB

## ACCUCOLD 24" WIDE BLANKET WARMER



90 to 140°F temperature range  
Designed and set to keep blankets warm in hospitals, medical and  
physical therapy offices, spas, and other facilities

[Read More](#)

**SKU:** PTHC65GLHD

## MAC MEDICAL WARMING CABINET



MAC Medical Warming Cabinet

Overall Size  
30"W x 26.5"D x 36"T

[Read More](#)

**SKU:** AM-SWC 36-G-TL-R



## MAC MEDICAL WARMER STAINLESS STEEL DOORS

Mac Medical Warmer Stainless Steel Doors  
Overall Size  
26.5"D x 30"W x 74.5"H

[Read More](#)

**SKU:** AM-DWC243074T-4B-SS



## STERIS DJ04 DUAL COMPARTMENT BLANKET WARMER

Steris Model DJ04

[Read More](#)

**SKU:** AM-DJ04

## SINGLE COMPARTMENT HEAT TABLE TOP WARMING CABINET



CMP Table Top Single Compartment Heat Warming Cabinet

[Read More](#)

**SKU:** AM-DSW3AE-R



## MAC MEDICAL WARMING CABINET (REFURBISHED)



**Refurbished Mac Medical 18" x 24" Warming Cabinet**

The MAC Medical Warming Cabinet is perfect for keeping IV, surgical irrigation fluids, linens, and blankets warm. The warming cabinet provides independent, digitally controlled heating chambers that offer actual temperature and setpoint displays. Data is recorded every 30 minutes and stored for up to two years

[Read More](#)

**SKU:** AM-SWC1824-R



## BLICKMANN DIGITAL TABLETOP GLASS FRONT WARMER

The Blickman Warmer Cabinet Model 7922 is a Single-door warming cabinet and provide ideal flexibility for warming solutions and blankets in heated Shelves.

[Read More](#)

**SKU:** AM-B7922SGDR



## GE OHMEDA 4400 INFANT WARMER SYSTEM

The Ohmeda 4400 Infant Warmer uses a manual mode of operation that allows selection of the level of radiant heat output as indicated by the % power display on the control panel and the control circuit then maintains the selected level of radiant heat

[Read More](#)

**SKU:** AM-GE-4400



## **BLICKMANN 7929MG TABLETOP GLASS DOOR BLANKET WARMER**

Blickman line. Used for warming solutions or blankets where space is at a premium and/or where multiple locations have modest warming requirements—buy one 'Ready Warmer' for each location! It is suggested that this unit sits on an existing cabinet, or on a Blickman transport dolly or mobile stand. Comes complete with digital control from 90°F to 140°F with keyed lockout, fan for even heat distribution, glass door, 120 volt, UL-listed and CUL-certified. Solutions-only warming cabinets can be pre-programmed to heat between 90°F and 110°F.

[Read More](#)

**SKU:** AM-B7929MG

## **SINGLE CHAMBER BLANKET/FLUID WARMER**



Blanket/Fluid Warmer -Single Chamber -Overall size: 30W x 26.6D x 36H - Chamber size: 26W x 23D x 25H -2 adjustable shelves -8.63 cubic feet of storage -Temp range 90-160 F.

[Read More](#)

**SKU:** AM-SWC36-TL

## OLYMPIC WARMETTE 56920 BLANKET WARMER



The Olympic Warmette® Model 20 Warming Cabinet features a compact design which allows medical facilities to have warm blankets and solutions closer to where they are needed. The Warming Cabinet can be placed in or near rooms that are in need of warm blankets and solutions thereby eliminating the process of having to retrieve them from an inconvenient, centrally-located warmer. Convenient placement of the Blanket Warming Cabinet allows professionals to provide solutions and blankets to patients that are closer to the desired storage temperature.

[Read More](#)

**SKU:** AM-566920

## BLICKMANN DIGITAL FULL-SIZE WARMER



Model 7924D (solid door) double door warming cabinets (blanket warmers) are 74-1/2" tall and are the perfect answer for warming blankets and solutions in two separately heated compartments. Set temperature control via key and pass-code lockout. Self-closing, magnetic gasket, adjustable tension doors. No set temperature restriction between compartments. Each warming cabinet unit comes standard with integrated casters for ease of mobility.

[Read More](#)

**SKU:** AM-7924D

## BLICKMAN 7927TG WARMING CABINET



Models 7927TG (glass door) warming cabinets (blanket warmers) are a compact size for when space is at a premium. These single compartment countertop blanket warmers come standard with self-closing magnetic gasket, adjustable tension doors. An adjustable shelf is available as an option. It measures 20 - 5/8" deep, allowing for placement on most counter or tabletop surfaces. Unit is fully self contained with no heat vents, as are all Blickman blanket warmers. Temperature control from 90°F - 160°F. To set up a solutions only warming cabinet, the unit can be easily programmed by the user to heat only between 90F - 110F

[Read More](#)

**SKU:** AM-B7927TG

## BLICKMANN DIGITAL TABLETOP WARMER



The Blickman Warmer Cabinet Model 7922 is a Single-door warming cabinet and provide ideal flexibility for warming solutions and blankets in heated Shelves

[Read More](#)

**SKU:** AM-B7922SSDR



## **PEDIGO 2120 FLUID WARMER REFURBISHED WITH CASTER BASE**

Deluxe fluid warmers feature 4 digit LED display, window doors, control lock programming and auto time programming

[Read More](#)

**SKU:** AM-P2120-R-W



## **STERIS TABLE TOP BLANKET WARMER MODEL QDJ03**

The Steris / Amsco QDJ-03 series of tabletop blanket warmers is a perfect solution to save space in your department or surgery center. The QDJ-03 is a cost-saving solution for blanket warmer needs in your medical facility versus the cost of a comparable new unit. Comes with a rack and is completely refurbished to excellent operational and cosmetic condition.

Inner Dimensions: 13-1/2 x 24 x 22-7/8"

Outer Dimensions: 29.25"W x 26.5"Deep

[Read More](#)

**SKU:** AM-QDJ03

## AMSCO QDJ01 TABLETOP WARMING CABINET



STERIS's Amsco Warming Cabinets are designed to be more than just blanket warmers; they also store and warm sterile intravenous fluids (IV fluids) and surgical irrigation solutions, as well as other linens, to programmed digital temperature settings for patient applications.

18 Inches Deep x 30 Inches Wide / Countertop  
Single Compartment Warming Cabinet

[Read More](#)

**SKU:** AM-QDJ01



## STERIS QDJ04 DUAL COMPARTMENT BLANKET WARMER

The **Steris QDJ04 Dual Compartment Blanket Warmer** is a state-of-the-art medical appliance designed to enhance patient comfort by ensuring blankets are warm and ready for use. Featuring two separate compartments with independent temperature controls, this warmer allows for precise and customizable heating settings for both the top and bottom sections. Its contemporary design, coupled with robust stainless steel construction, ensures durability and a sleek appearance suitable for any medical setting. The unit is equipped with a new circulation fan and digital displays, providing efficient heat distribution and easy monitoring of temperature settings. Overheat indicators add an extra layer of safety, alerting users to any irregularities in temperature. Completely remanufactured and biomedically inspected, this blanket warmer meets high standards of quality and performance. The inner dimensions of the upper compartment are 23.5 inches in width, 22 inches in depth, and 13 inches in height, while the bottom compartment measures 23.5 inches in width, 22 inches in depth, and 38 inches in height. The overall dimensions of the unit are 30 inches in width, 28 inches in depth, and 75.5 inches in height, making it a spacious and efficient solution for medical facilities.

[Read More](#)

**SKU:** AM-QDJ04



## CMP BLANKET AND SOLUTION WARMING CABINET

CMP has been building high quality warming cabinets for over 65 years and has become the industry standard for hospital and emergency room blanket and solution storage. CMP stainless steel Warming Cabinets can be found in almost every hospital in the U.S. and can be configured for your unique requirements.

Models: DSW1AE-18, DSW1AE-18 (74H), DSW1AE-24, DSW1AE-24 (74 H), DSW2AE-1-18, DSW2AE-1-24, DSW2AE-2-18, DSW2AE-2-24, DSW2AE-2-24, DSW3AE-18 (18 W), DSW3AE-24

[Download Brochure](#)

[Read More](#)

**SKU:** AM-DSW2AE-2-R



## BLICKMAN WARMING CABINET MODELS 7921TG AND 7921T

Models 7921TG and 7921TS single door warming cabinets are 74-1/2" tall and ideal for warming blankets or solutions. An open and adjustable shelf design adds to flexibility and expanded use. Easy-to-use controls for setting internal temperature are conveniently located above the door. Each unit comes standard with integrated casters for ease of mobility.

[Download Brochure](#)

[Read More](#)

**SKU:** N/A



## STERIS AMSCO QDJ06 DUAL COMPARTMENT WARMING CABINETS

The Steris Amsco Warming Cabinet is the solution for all warming and storage needs within the operating room. This Steris Cabinet is designed as a blanket warmer as well as a warmer and storage container for sterile IV fluids and irrigation supplies. The newly designed system keeps linens and fluids warm, secure and compliant at the proper temperature settings.

[Download Brochure](#)

[Read More](#)

**SKU:** AM-QDJ06



## STERIS AMSCO QDJ05 SINGLE COMPARTMENT WARMING CABINET

The Steris Amsco Warming Cabinet is the solution for all warming and storage needs within the operating room. This Steris Cabinet is designed as a blanket warmer as well as a warmer and storage container for sterile IV fluids and irrigation supplies. The newly designed system keeps linens and fluids warm, secure and compliant at the proper temperature settings. The Steris QDJ-05 Single Compartment Warming Cabinet is equipped with digital temperature settings to ensure proper the proper warming temperature.

[Download Brochure](#)

[Read More](#)

**SKU:** AM-QDJ05

## SKYTRON J2 SERIES WARMING CABINETS

**Skytron's J2 Series Warming Cabinets** are designed to maintain clinically ideal temperatures for solution bottles, solution bags and blankets. Six models are available to provide optimum clinical storage requirements.

**Counter Top Models:** SS2201-J2, SS2201-J2G

**Free Standing Models:** SS2205-J2, SS2205-J2G, SS2207-J2, SS2207J2G

[Download Brochure](#)

[Read More](#)

**SKU:** N/A



## MAC MEDICAL WARMING CABINETS

These warmers are equipped with independent, digitally controlled heating chambers that offer actual temperature and set point displays. Simple plug and play data download requires no additional software.

[Download Brochure](#)

[Read More](#)

**SKU:** N/A



## STERIS ANALOG TABLETOP WARMER



The Steris Amsco M70 Series line of blanket warmers was one of the most popular lines ever manufactured. This series features a tabletop version with a rotary knob that controls temperature. The Amsco® M70® also features over-temperature alarms (visual and audible) and a heating indicator (visible). On full size models, the upper compartment of the Amsco® M70® is the same as the tabletop version, but with a separate lower cabinet that is larger in size than the upper. The lower cabinet of the full size model comes as either a non-heated storage compartment or a heated version, which would maintain a constant 120° F temperature for warming of blankets only.

[Download Brochure](#)

[Read More](#)

**SKU:** AM-M70WC-E



# INDEX

- 7**  
750 Bair Hugger without Cart 15
- A**  
Accucold 24" Wide Blanket Warmer 17  
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Left Hand Door Chrome Shelves 4  
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Left Hand Door SS Shelves 6  
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Right Hand Door Chrome Shelves 7  
Accucold 24"W x 71"T Full-Size Wide Blanket Warmer Right Hand Door SS Shelves 8  
Amsco QDJ01 Tabletop Warming cabinet 25  
Arizant 775 Bair Hugger 15
- B**  
Blickman 7927TG Warming Cabinet 23  
Blickman Warming Cabinet models 7921TG and 7921T 26  
Blickmann 7929MG Tabletop Glass Door Blanket Warmer 21  
Blickmann Digital Full-size Warmer 22  
Blickmann Digital Tabletop Glass Front Warmer 20  
Blickmann Digital Tabletop Warmer 23
- C**  
CMP Blanket and Solution Warming Cabinet 26
- E**  
Enthermics EC230L Fluid Warmer 2
- G**  
GE Ohmeda 4400 Infant Warmer System 20
- M**  
MAC Medical SWC183024 Glass Door Warming Cabinet 4  
Mac Medical Warmer Stainless Steel Doors 18  
Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Glass Door, 4" Base, Dual Chamber 9  
Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Glass Steel Door, 4" Base, Single Chamber 10  
Mac Medical Warmer, D-Series, 26.5"D X 30"W X 74.5"H, Stainless Steel Door, 4" Base, Dual Chamber 9  
MAC Medical Warmer, D-Series, 26.5inD X 24inW X 24.5inH 16  
MAC Medical Warmer, D-Series, 26.5inD X 30inW X 24.5inH 16  
MAC Medical Warming Cabinet 17  
MAC Medical Warming Cabinet (New) 14  
MAC Medical Warming Cabinet (Refurbished) 19  
MAC Medical Warming Cabinets 28
- O**  
Olympic Warmette 56920 Blanket Warmer 22
- P**  
Pedigo 2120 Fluid Warmer Refurbished with Caster Base 24

## **S**

Single Chamber Blanket/Fluid Warmer 21  
Single Chamber Blanket/Fluid Warmer Glass Door 13  
Single compartment heat table top warming cabinet 19  
Skytron J2 Series Warming Cabinets 28  
Steris AMSCO QDJ05 Single Compartment Warming Cabinet 27  
Steris AMSCO QDJ06 Dual Compartment Warming Cabinets 27  
Steris Analog Tabletop Warmer 29  
Steris DJ02 Dual Compartment Blanket Warmer 3  
Steris DJ04 Dual Compartment Blanket Warmer 18  
Steris Dual Compartment Blanket Warmer DJ06 11  
Steris Dual Compartment Blanket Warmer QDJ02 2  
Steris Dual Compartment Blanket Warmer QDJ05 14  
Steris QDJ04 Dual Compartment Blanket Warmer 25  
Steris table top blanket warmer model QDJ03 24  
Steris Tabletop Blanket Warmer - Model DJ03 12  
Steris Tabletop Blanket Warmer Model DJ01 12