

# TABLE OF CONTENTS

**Washer Disinfectors** ..... 2

    Steris Reliance 400XLS Washer ..... 3

    Miele PLW 6111 Laboratory washer with 3 wash levels in SlimLine design  
    ..... 3

    Steris Reliance 400/500XLS Washer Racks with Spray Arms ..... 4

    Steris Reliance 500XLS Washer ..... 5

    Tuttnauer TIVA 10 Lab Washer, 265 Liter, With 2 Automatic Vertical Door, 3  
    Phase ..... 6

    Tuttnauer TIVA 10 Lab Washer, 265 Liter, With 1 Automatic Vertical Door, 3  
    Phase ..... 6

    Tuttnauer TIVA 10 Lab Washer, 265 Liter, With 1 Manual Door, 3 Phase  
    ..... 7

    Tuttnauer Tiva 8 Laboratory Undercounter Washer Disinfector ..... 8



## WASHER DISINFECTORS

**Maximize your facility's productivity and safety with one of Auxo Medical's cost-effective and durable refurbished thermal washer disinfectors.**

We focus on the highest quality machines and dependability that you can count on. As a result, Auxo Medical distributes washer disinfector equipment from some of the most respected brands in the medical industry. Additionally, this includes models from [Steris](#), [Steelco](#), Scientek, and Miele. Furthermore, to help your facility equip itself with the best washer disinfector for your needs, Auxo Medical offers a variety of sizes and load capacities. Whether you're looking for lower utility costs, faster processing times, or less environmental impact. Ultimately, our washer disinfectors are ready for the task.

**Choose from applications for:**

[HEALTHCARE](#)

[LIFE SCIENCES](#)

**Our catalog consists of brand new, reconditioned and refurbished equipment. Would you like to filter the equipment by:**

[BRAND NEW](#)

[RECONDITIONED](#)

[REFURBISHED](#)

## STERIS RELIANCE 400XLS WASHER

The **Steris Reliance 400XLS Laboratory Glassware Washer** is a robust, cabinet-style washing system designed for efficient cleaning of laboratory glassware and instruments. Equipped with advanced Allen-Bradley® control systems, it offers five preprogrammed adjustable cycles with five additional customizable options, along with specialized descaling and priming cycles for routine maintenance. Built to meet seismic building codes, it is available in single- or double-door configurations, suitable for freestanding or recessed installation through a barrier wall. The washer features a spacious chamber measuring approximately 26" W x 25-1/8" H x 26" D, with overall dimensions of 42" W x 80" H x 33-7/8" D, and is powered by either steam or electric condensing systems. It requires compressed air (utilities supplied separately) and weighs around 1,500 lbs. Please note, loading equipment is not included and must be selected separately, making this washer ideal for demanding lab environments requiring reliable, customizable cleaning performance.

[Reliance 400XLS / 500XLS Laboratory Glassware Washers](#)

[Technical Data Sheet](#)

[Reliance 400XLS /500XLS Laboratory Glassware Washer](#)

[Accessories](#)

[Read More](#)

**SKU:** AM-400XLS

**Price:** Request Quote





## MIELE PLW 6111 LABORATORY WASHER WITH 3 WASH LEVELS IN SLIMLINE DESIGN

The **Miele PLW 6111 Laboratory Washer**, designed in a SlimLine format, is a high-performance equipment tailored for meticulous reprocessing of laboratory glassware. It boasts three wash levels and a suite of advanced features, including a steam condenser, boiler, printer, conductivity meter, and chamber lighting. This washer can handle a substantial load per cycle, efficiently cleaning 126 narrow-neck glassware or 121 pipettes, with a usable capacity of 7.416 cubic feet. Its ergonomic design, enhanced by telescopic runners, ensures ease during loading and unloading. The washer offers flexible operation with four level positions and a full-glass door for clear visibility. It is capable of reprocessing containers up to 13 gallons and comes equipped with two integrated hose dispenser pumps for liquid media. Additional features include a boiler for preheating deionized water, volume-flow control, conductivity monitoring, and a multi-component filter system. The unit supports comprehensive process documentation via an integrated printer and RS 232 serial interface. The wash cabinet, made from high-quality stainless steel (1.4404/316L), is illuminated for better visibility. With a control system of MP.CON 3, it offers 18 available programs, including preprogrammed cycles for thermal disinfection and reprocessing of laboratory glassware and plastics. The Miele PLW 6111 is a robust, versatile, and efficient solution for laboratory cleaning needs, ensuring high capacity and reliability.

[Read More](#)

**SKU:** AM-PLW6111

**Price:** Request Quote



## STERIS RELIANCE 400/500XLS WASHER RACKS WITH SPRAY ARMS

The **Steris Reliance 400/500XLS Washer Racks**, equipped with efficient spray arms, represent a significant advancement in the sterilization and cleaning process of medical instruments. These racks are designed to fit seamlessly into the Steris 400/500XLS washer models, ensuring optimal compatibility and performance. Each rack, capable of holding various surgical tools, features strategically placed spray arms that ensure thorough cleaning by distributing water and cleaning agents evenly across all instruments. The image displaying six racks highlights their capacity and modular design, which facilitates customization based on the specific needs of a medical facility. This system not only enhances the efficiency of the cleaning cycle but also contributes to the overall reliability and effectiveness of the sterilization process, crucial for maintaining high standards of hygiene and patient safety in any healthcare setting.

**IMPORTANT: Product is a Single Rack with Spray Arms (Photo shows Qty of 6)**

[Read More](#)

**SKU:** AM-REL-XLS-R

**Price:** Request Quote



## STERIS RELIANCE 500XLS WASHER

The **Steris Reliance 500XLS Laboratory Glassware Washers** represent a pinnacle of efficiency and customization in lab cleaning technology. These cabinet-type washers are designed with robust Allen-Bradley® control systems and offer a versatile array of ten programmable cycles, including five that are pre-set and five that can be customized to suit specific laboratory needs. An additional programmed descaling and a priming cycle ensure the equipment is maintained in optimal condition, enhancing longevity and performance. The washers conform to seismic building codes, providing additional safety and compliance assurance, with a seismic report included in the standard documentation. Available in both single- and double-door models, the 500XLS can be installed as a freestanding unit or recessed through a barrier wall, offering flexibility in laboratory design. The unit features two adjustable peristaltic pumps and a non-vented condensing system, ensuring efficient operation without the need for external venting. With a substantial chamber load capacity and a powerful 480 Volt, 60 Hz, 3-Phase electrical requirement, this all-electric washer does not require steam, making it an eco-friendly and efficient choice for modern laboratories. The operational weight of 1500lbs and generous dimensions further underscore its robust build and capacity for

handling substantial glassware loads, making the Reliance 500XLS an essential tool for any high-volume laboratory environment.

[Reliance 400XLS / 500XLS Laboratory Glassware Washers Technical Data Sheet](#)

[Reliance 400XLS / 500XLS Laboratory Glassware Washer Accessories](#)

[Read More](#)

**SKU:** AM-500XLS

**Price:** Request Quote



## TUTTNAUER TIVA 10 LAB WASHER, 265 LITER, WITH 2 AUTOMATIC VERTICAL DOOR, 3 PHASE

The **Tuttnauer TIVA 10 Lab Washer**, boasting a generous 265-liter capacity and equipped with 2 automatic vertical sliding doors, represents an advanced solution for the comprehensive cleaning, disinfection, and drying of a wide range of medical instruments. Measuring 26.77in x 25.2in x 72.83in with an interior dimension of 22.05in x 23.03in x 31.5in, this model is designed to fit snugly into limited spaces while maximizing its load capacity. Its high-quality AISI 316L washing chamber and capacity to accommodate up to 12 DIN baskets across 6 washing levels make it particularly suitable for surgical and sterilization departments. The TIVA 10-V also features a user-friendly touch panel with a color graphic display offering 40 preset and programmable washing cycles, ensuring adaptability to various load types. Moreover, its advanced features, such as air drying with a prefilter and HEPA H14 filter, a built-in steam condenser, three heating elements, and two dosing pumps for precise chemical application, underscore its efficiency and effectiveness. Additionally, the washer includes hands-free door operation, enhancing user convenience and safety. The package comes complete with 1 bottle each of T-Clean Detergent and T-Clean Neutralizing Agent, alongside necessary connections for hot and cold water, deionized water, drainage, and electrical requirements (208V-240V) for 3-phase power supply, making the TIVA 10-V an all-encompassing solution for the modern sterilization department.

[Read More](#)

**SKU:** AM-TIVA10-2V-MD

**Price:** Request Quote



## **TUTTNAUER TIVA 10 LAB WASHER, 265 LITER, WITH 1 AUTOMATIC VERTICAL DOOR, 3 PHASE**

The **Tuttnauer TIVA 10 Lab Washer**, with its 265 Liter capacity and an automatic vertical door, is a state-of-the-art solution designed to meet the rigorous demands of surgical and sterilization departments. Its dimensions of 26.77in (W) x 25.2in (D) x 72.83in (H), and an interior space of 22.05in x 23.03in x 31.5in, allow for a substantial loading capacity, accommodating up to 12 DIN baskets across 6 washing levels. This flexibility is further enhanced by the availability of trolleys for various load types. The unit's advanced touch control provides users with a color graphic display and access to 40 different programs, both preset and programmable, ensuring a broad range of applications. Its high-quality AISI 316L washing chamber, automatic sliding glass doors for space efficiency, air drying with prefilter and HEPA H14 filter, built-in steam condenser, and two dosing pumps, are among the advanced features that make the TIVA 10-V an ideal choice for those seeking a comprehensive decontamination cycle. Additionally, the machine is equipped with a hands-free door operation sensor, includes bottles of T-Clean Detergent and Neutralizing Agent, and supports multiple connections for water and drainage, including a 3/4" HW connection, 3/4" CW connection, 1/2" deionized water, and a 1" drain, along with a 208V-240V electrical 3-phase system, rounding off its extensive list of high-end features designed for efficiency, safety, and user convenience.

[Read More](#)

**SKU:** AM-TIVA10-1V-MD

**Price:** Request Quote



## TUTTNAUER TIVA 10 LAB WASHER, 265 LITER, WITH 1 MANUAL DOOR, 3 PHASE

The **Tuttnauer TIVA 10 Lab Washer**, boasting a generous 265-liter capacity, is an essential addition to any laboratory prioritizing efficiency and reliability in their cleaning processes. This model, equipped with a single manual door, offers a user-friendly interface that ensures ease of operation while maintaining the high standards of cleanliness required in a lab environment. Designed for heavy-duty use, the TIVA 10 operates on a 3-phase power supply, providing the robust performance needed for thorough cleansing of lab instruments and glassware. Its ample size and advanced features make it an ideal choice for laboratories looking to streamline their cleaning protocols without compromising on quality or performance. Whether for a research facility, a pharmaceutical lab, or a university science department, the Tuttnauer TIVA 10 Lab Washer stands out for its efficiency, capacity, and reliability, making it a valuable asset in maintaining a sterile and productive lab environment.

[Read More](#)

**SKU:** AM-TIVA10-1M-L

**Price:** Request Quote



## TUTTNAUER TIVA 8 LABORATORY UNDERCOUNTER WASHER DISINFECTOR

UP TO 4 washing levels of total 8 DIN baskets capacity, 175 Liter chamber

[Read More](#)

**SKU:** AM-TIVA8-L

**Price:** Request Quote

# INDEX

## M

Miele PLW 6111 Laboratory washer with 3 wash levels in SlimLine design 3

## S

Steris Reliance 400/500XLS Washer Racks with Spray Arms 4

Steris Reliance 400XLS Washer 2

Steris Reliance 500XLS Washer 4

## T

Tuttnauer TIVA 10 Lab Washer, 265 Liter, With 1 Automatic Vertical Door, 3 Phase 6

Tuttnauer TIVA 10 Lab Washer, 265 Liter, With 1 Manual Door, 3 Phase 7

Tuttnauer TIVA 10 Lab Washer, 265 Liter, With 2 Automatic Vertical Door, 3 Phase 5

Tuttnauer Tiva 8 Laboratory Undercounter Washer Disinfectors 7