



## Baxter Flo-Gard 6201 Infusion Pump - Fully Refurbished

The Flo-Gard 6201 Volumetric Infusion Pump from Baxter can deliver a wide variety of fluids over a broad range of infusion rate. Enhanced safety, greater convenience and versatility, and Baxter standard IV administration sets that eliminate the need for costly disposable sets and reduce the potential for waste make the Flo-Gard 6201 the perfect infusion pump.

### Features:

- Flow Check Occlusion Alarm-an in-line resistance display of incremental back pressure
- Flow Rate Calculation is automatic after volume and time are selected
- Slide Clamp Option offers an additional step to protect against inadvertent gravity free flow
- Front Panel Lock-out protects against tampering
- Programmed Delivery Profile allows up to 10 steps, individualizing control of infusion ramping and tapering
- Incremental Flow Rate can deliver 1-99.9 mL /hr in increments of 0.1 mL /hr
- Automatic Restart once occlusion clears
- Automatic Piggybacking of secondary medications
- Channels: Single
- Range of Flow: Rate 1 - 99.9 mL/hr in 0.1 mL/hr increments and 1 - 1,999 mL/hr in 1 mL/hr increments
- Maximum VTBI: 1 - 9,999 mL
- Alarm Intervals: 1-7 seconds
- Battery/Time: 12 volt, 2.0 ampere, sealed lead acid
- Height: 11.4"
- Width: 8"
- Weight: 11.6 lbs



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.