



## Midmark 630 Power Procedure Table - Fully Refurbished

---

No matter what type of procedures you perform, proper positioning is crucial to providing accessibility and flexibility at the point of care. Power height, back, foot and tilt gives the Midmark 630 virtually endless movement capabilities, allowing you to position patients correctly for the most effective procedure. And, adjustable height allows you to work in a position that is comfortable for you.

### Features:

- Patient Weight Capacity: 450 lb. (204 kg)
- Height: Minimum 19" (48.3 cm) Maximum 40" (101.6 cm)
- Seat Tilt Angle Range: 0° to 30°
- Back Support Angle Range: 0° to 85°
- Foot Support Angle Range: 0° to 90°
- Upholstery Width Dimensions:
- Headrest: 12" (30.5 cm)
- Back support: 24"-26" (61-66 cm) Seat: 28" (71.1 cm)
- Foot support: 18" (45.7 cm)
- Patient Support Surface Length:
- Flat position: 72" (182.9 cm)
- Flat position with headrest extended: 83" (210.8 cm)
- Paper Roll Holder: Holds 18" x 3.5" paper roll (45.7 x 8.9 cm)
- Stirrups are standard



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.