



Midmark 75L Exam Chair - Fully Refurbished

The Midmark Ritter 75L offers full eight-way power, essential low height, unprecedented patient comfort and a multitude of efficiency features. With a low height of 19 inches and an exclusive Premium Comfort System, the 75L provides easy access for patients who are pregnant, elderly, or have disabilities. This enables patients to access the table themselves, while helping you and your staff position patients efficiently.

Features:

- Table Up/Down
- Tilt
- Foot Up/Down
- Back Up/Down
- Adjustable Headrest
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder
- Table Paper Straps
- Foot Switch or Hand Control
- Accessory Side Rails
- Width: Seat 28" (68.6 cm), Back 24" - 26" (58.4 - 66.0 cm)
- Length with foot section extended: 72" (182.9 cm)
- Length with foot section and headrest extended: 83" (210.8 cm)
- Height: Minimum 19" (48.26 cm) Maximum 40" (101.5 cm)
- Patient Load Rating: 450 lb. (204.1 kg)
- Paper Roll Holder: Holds 18" x 3.5" paper roll (45.7 x 8.9 cm)



Why an Auxo Medical Refurbished Chair or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each chair and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each chair and stretcher is reassembled, tested and is warranted against defects.

