



## Midmark Ritter 222/223 Exam Table— Fully Refurbished

With an increasingly aging population, your patients may need assistance in accessing the exam table. With an incredible 18-inch low height, the Ritter 222/223 Barrier-Free Power Exam Table makes it easy for your practice to accommodate the elderly, expectant mothers or patients with disabilities easily and efficiently.

### Features:

- Length with footrest extended: 76.7" (194.7 cm)
- Height: Minimum 18" (45.7 cm) Maximum 37" (94.0 cm)
- Patient weight capacity: 400 lbs (181.4 kg)
- Paper roll holder (hidden): Holds 21" x 3.5" paper roll (53.3 x 8.9 cm)
- Electrical receptacle: (115 VAC 5 Amps Max.) Duplex, hospital grade, electrical receptacle on left side
- Pelvic tilt (5 degrees)
- Drawer heater (35 watt)
- Stirrups are standard
- Unprecedented low and high positioning heights
- Counterbalanced
- Gas spring assisted back section assists in placing the back in an infinite number of positions
- Easy access backrest mechanism
- Built-in pillow
- Hideaway stirrups that pull out to lock in one of four lateral positions
- Foot extension shelf to expand the table to its fullest length
- Polystyrene treatment pan
- Exam Assistant Drawer System
- Leveling screws



### Why an Auxo Medical Refurbished Table or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.