



## Physio-Control Lifepak 20 - Fully Refurbished

The world leader in defibrillation technology brings you the improved LifePak 20 defibrillator/monitor, created especially for hospitals and clinics. Designed with extensive input from clinicians around the world, the 20 is simple yet sophisticated, providing a flexible and effective tool to meet diverse defibrillation and monitoring needs.

### Features:

- AHA/ERC Guidelines 2005 consistent cprMAX™ technology
- Simple yet sophisticated for responders of many skill levels
- Highly intuitive for quick, effective AED use
- Easily converts to manual mode
- ADAPTIV™ biphasic technology adjusts shock duration and voltage
- Broad dosage capability up to 360J, when needed
- Compact and lightweight for easy carrying and versatile use
- Superb color graphics allow easy viewing from many angles

### Specifications:

- Weight:
  - Fully featured defibrillator/monitor (pacing and SpO2 ) 5.58 kg (12.3 lbs)
  - QUIK-COMBO cable: .20 kg (.43 lbs)
  - Standard (hard) paddles: .88 kg (1.95 lbs)
- Height: 21.3 cm (8.4 in)
- Width: 26.2 cm (10.3 in)
- Depth: 26.2 cm (10.3 in)

**Additional Accessories:** AED Case, Hard Paddles



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.