



# auxoMEDICAL

## Zoll AED Plus - Fully Refurbished

When sudden cardiac arrest occurs, the fact is that only half of the victims will need a shock. You need an automated external defibrillator (AED) that helps you all the time. The ZOLL AED Plus is the only AED that can actually see when you are doing CPR and help you do it well. The AED Plus features Real CPR Help, a unique CPR feedback tool that is able to monitor CPR quality and provide real-time feedback for depth and rate of chest compressions. Audio and visual prompts help you to rescue with a confidence and clarity unmatched by any other AED.

### Features:

- Roadmap for resuscitation combining pictures, text displays, and voice prompts
- One-Piece reconnected pair of electrodes
- Intelligent pediatric capability
- Interactive feedback for rate and depth
- Consumer batteries for convenience and savings
- Rectilinear biphasic waveform
- Supports full chain of survival best
- Improves CPR performance the most
- Includes 1 set of electrodes, 1 set pf lithium ion batteries, soft carrying case, and manual



### Specifications:

- Dimensions: 5.25" x 9.50" x 11.50"
- Weight: 6.7 lbs.

### Additional Options:

- Extra CPR-D Padz Electrodes
- Pedi Padz II Electrodes

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.