

## Steris AMSCO SQ240 Light - Fully Refurbished

The Steris AMSCO SQ240 combines the flexibility and positioning ease of a small-diameter lighthead with the superb optical performance of a large-diameter lighthead. It's a high-powered, easy-to-manuever lighting system designed for every size and type of surgical procedure.

### Features:

- Bright consistent illumination and constant pattern definition.
- Large depth of field eliminates the need to continually refocus during a procedure.
- Excellent shadow control and maximum cavity illumination.
- Optimal brightness is achieved without glare.
- Skin and operative tissue appear in a natural state without color distortion.
- Easy to grip sterile handle makes for easy positioning and adjustment of pattern size.
- Exceptional ease of maneuvering during procedures.
- Drift-free suspension and counterbalance arms.
- Easy access non sterile handles for easy positioning by the circulating nurse
- Good range of light intensity adjustment.
- Automatic-Bulb changeover mechanism ensures continuation of illumination even when the primary bulb fails.
- Easy change bulb that can be changed between cases in a matter of minutes.
- The high quality, high gloss finish facilitates ease of cleaning at the end of a long day of cases.
- Light intensity: 129,170 lux
- Color temperature: 4,400 K
- CRI: 92
- Depth of field: 26"
- Bulb life: 1000 hours



### Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.