



Steris Harmony 585 Lights - Fully Refurbished

Harmony LED Surgical Lighting and Visualization System optics provide exceptional surgical site illumination while its HD-capable suspension renders virtually every image signal in stunning clarity. Harmony LED lights also reduce energy consumption by 3/4, have no bulb change requirements and give you a cooler, more comfortable room with 34% less beam heat.

Features:

- Brighter light - up to 41% more illuminance than other LED lights
- Pure white light - 4400K color temperature \pm 300k
- A large, bright spot - 160,000 lux max
- 84 overlapping beams of light virtually eliminate shadows from the surgeons head and hands
- Better tissue discrimination with a CRI of up to 96
- Deep, saturated reds with an R9 value of up to 98
- Accommodates virtually all HD signals - regardless of manufacturer or signal type (DVI, HD-SDI, RGBHV, etc)
- Grows with you - add lightheads, ambient light, monitors and cameras
- Add Harmony ConnectPoint for video integration without costly or involved renovations
- Convert general rooms to specialty suites quickly and easily
- Count on years of reliability due to the unique cable management system
- Color temperature: 4400K
- Max Central Illuminance: 160,000 lux (14,870 fc)
- Life: 30,000 hrs
- Average pattern sizes: 7-11" (18-28 cm)



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.