



Steris Harmony LC Lights - Fully Refurbished

The Harmony LC Surgical Lighting System offers extraordinary value by providing advanced optics with patented features in a cost-effective system that can be configured to meet a variety of surgical lighting and video applications.

Features:

- Intensity (LUX): 160,000 lux
- Intensity (foot candles): 14,900 f.c.
- Pattern size range: 6"-11" (15-28 cm)
- Depth of field (small pattern): 42" (107 cm)
- Light head diameter (O.D.): 23" (58 cm)
- Color temperature: 4,400K
- CRI (Color Rendering Index): 94
- Performance and illumination of a large light head in a smaller size that is easier to handle
- Filter heat away from the surgical site through faceted cold-mirror reflectors
- Controlling Options:
 - Can be controlled by circulating nurse with a wall-mounted control panel
 - Can be controlled by the surgeon with the integrated controls at the light handle
- Wide range of configurations:
 - Up to 3 light heads on the same ceiling mount
 - Or, up to 4 flat screen monitors
 - One of the light heads can be camera-ready
 - Choose from one of four sizes of suspension system



Why purchase Auxo Medical Refurbished Lights?

Auxo Lights, halogen and LED, are refurbished to the highest standards in the industry. We completely disassemble each light and arm assembly, and completely rebuild or replace all worn parts and return each assembly to OEM standards. Each Light is reassembled, tested and is warranted against defects.