

Brewer Access Exam Table - Fully Refurbished

The Brewer Access Exam Table includes great features such as the largest and most secure four-leg patient step, pass-through side drawers, pneumatic backrest, and bariatric weight capacity of 500 pounds. The backrest is electrically adjustable. The Brewer Access Exam Table is one of the safest exam tables with the greatest storage capacity available.

Features:

- Height: 32" (81.3 cm)
- Width: 28" (71.1 cm)
- Length: 57.3" extends to 71" (145.5 cm – 180.3 cm)
- Leg extension length: 15" (38.1 cm)
- Step Height: 8" (20.3 cm)
- Step Length: 17.7" (45 cm)
- Step Width: 20.4" (51.8 cm)
- Hospital-grade duplex electrical outlets
- 115VAC, 60Hz, 5.8 amps max
- Outlets can be put on left or right side
- Easy to clean seamless upholstery
- Simple pneumatic backrest for easy patient positioning
- Multi-density, bolstered support upholstery
- Leg extension/plastic debris tray
- Adjustable ergonomic stirrups
- Simple pneumatic backrest for easy patient positioning
- Pass-through side drawers accessible from right or left



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.