



Abbott Triple Channel Plum A+ IV - Fully Refurbished

The Plum A+3 Infusion System is a single unit that allows clinicians to manage an industry-leading six-infusion line medication management system throughout the continuum of care. The wireless-enabled system assists healthcare professionals to better manage the complex medication dosing regimens typically seen in emergency rooms (ERs), intensive care units (ICUs) and oncology centers. Patients in these settings are often critically ill and may require multiple drug infusions delivered at precise doses and on varying schedules. The Plum A+3 enables clinicians to deliver these infusions either simultaneously or in specific tiered sequences, depending on each patient's needs.

Features:

- Minimizes the time required to assemble the appropriate technology at the patient's bedside
- Enhances operational efficiency
- Uses a distributed approach to provide safeguards against multiple channels being simultaneously compromised
- Optimizes asset management and productivity
- Six infusion lines within one unit
- Rechargeable battery
- Single channel (2 concurrent medications) or triple channel (6 concurrent medications) models available
- Deliver secondary infusions via IV container or syringe
- Reliable performance, steady, consistent and proven
- Automated piggyback and concurrent delivery
- Programmable standby settings
- Multistep delivery
- Loading dose automation
- Programmable delayed starts



Specs:

Dimensions: 13.75(H) x 12.2(W) x 7.5(D)

Weight: 20 lbs with batteries

Power cord: Hospital grade AC cord, 10ft

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.