



Alaris Gemini PC-1 IV Pump - Fully Refurbished

Gemini PC-1 is a single channel volumetric infusion pump and controller which provides accurate and automatic infusion of intravenous drugs and fluids. The PC-1 may be operated as either a pump or a controller, or it may be configured as only one or the other using a lockout feature. Secondary (or piggyback) fluid may also be automatically infused at delivery rates and volumes independent of the primary infusion parameters, with automatic changeover to the primary infusion parameters when the secondary infusion is complete and when using an ADD-I-MED administration set. Computer operation is provided when the PC-1 is interfaced with a host computer. The technical data necessary to interface the PC-1 with a host computer is provided in the C-2 Programming Guide.

Features:

- Dual rate programming
- Panel lock prevents inadvertent changes to infusion program.
- Micro/Macro Infusion Capability
- Rapid Rate On-Line Titration
- Programmable parameters include:
 - Pump/Controller Mode
 - Micro, Macro and Universal Rate Modes
 - VersaTaper® and AutoTaper® Modes
 - "10 PSI" Mode
 - Panel Lock
 - Computer Link
 - Language Selection
 - Biotech Mode

Specs:

- Delivers precise infusions at rates from 0.1 to 999 mL/hr
- Height: 10.87"
- Width: 8.4"
- Depth: 6.5"
- Weight 11.1 lbs



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.