



auxoMEDICAL

Alaris MedSystem 2865 - Fully Refurbished

The MedSystem III® (MSIII) multi-channel infusion pump brings clinical versatility to your drug infusion technology. In one compact infusion pump, the MSIII incorporates a customized drug library, three independent fluid delivery channels, and a drug dose calculator enabling users to safely and quickly deliver infusion therapy. The Medsystem III's small size allows for more flexibility when transporting patients and allows for a more organized bedside space. The lightweight MedSystem III® pump is ideal for ambulatory, EMT and air transport. It features flexible and system parameters that can be set to match the clinical needs of the patient, as well as a unique, easy-to-use pole clamp, allowing the pump to be attached virtually anywhere.

Features:

- Three independent infusion channels
- Small, compact size and lightweight design
- Weight: 5.1 lbs
- Set-based, free-flow protection with needle-free access
- Delivers intravenous, intra-arterial, epidural and subcutaneous applications
- Infuses a wide range of fluids out of various containers, including syringes, utilizing a syringe adapter set
- 6-watts AC power
- Dose rate calculation, single-page programming, a selectable drug library and titration of drugs by the user in either dose or volumetric rate
- Rechargeable NiCad battery pack memory
- Non-rechargeable lithium backup battery
- Delivery rang: 0.1-999 mm/hr
- KVO rate range: 0.1-20.0 mm/hr
- Charger and pole clamp included
- Dimensions: 8 X 6 X 2 inches



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.