



auxoMEDICAL

Alaris Medley 8015- Fully Refurbished

The Alaris® System is a modular system intended for adult, pediatric and neonatal care in today's growing professional healthcare environment. It consists of the PC Unit, the Guardrails® Suite MX, and up to 4 detachable infusion and/or monitoring modules (channels). Alaris® System wireless communication makes it easier than ever to increase the safety of iv medication and continuously improve clinical best practices, regardless of existing wireless infrastructure.

Features:

- Intended for adults, pediatric, neonatal care in today's growing professional healthcare environment
- The Point Of Care provides a common user interface for programming infusion monitoring
- The system provides the right protection for every infusion by double checking dose and duration as well as delivery rate
- Wireless connectivity
- Powerful analytics
- Exceptional versatility
- Increased protection for IV safety
- Guardrails Suite MX software is the most comprehensive safety software available

Specs:

- Dimensions: 6.9" W x 8.8" H x 9" D including pole clamp
- Weight: 7.2 lbs
- Battery Operation: Battery run time is a function of the number of modules attached and module activity
- Communication Data Port: RS-232 with an RJ45 connector
- Electronic Memory: System configuration parameters stored in volatile memory are retained for at least 6 months by an internal backup lithium battery
- Module specific parameters are stored for 8 hours when system is turned off
- Power Requirements: 100-240V, 50/60 Hz, 150VA Max



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.