



auxoMEDICAL

Alaris Medley 8100 - Fully Refurbished

The Vista Basic was designed to ensure simplicity in all aspects of its use: from set installation, to programming, to special mode use and alarm correction. Based on a design that has sold over 100,000 units worldwide, the Vista Basic provides high reliability packaged in one of the easiest pumps to use. As a workhorse in the outpatient market, it offers a rugged and reliable design, all relevant therapies such as continuous, piggyback and program mode, and a selection of cost-effective, straight line gravity sets.

Features:

- Auto-Restart
- Bolus Dose Channel Labels allows a bolus infusion to be programmed using either the Guardrails® Drug Library or drug calculation feature
- Concentration Limits
- Delay Options allow system to be programmed to delay start of an infusion a) for up to 120 minutes or b) for a specific time up to 23 hours 59 minutes.
- Drug Calculation allows entry of drug dose and entry of flow rate
- Dynamic Pressure Display graphically displays current patient-side occlusion pressure set point and current patient-side operating pressure for that module
- Free Flow Protection utilizes a unique clamping device, the Flo-Stop® Device, to prevent inadvertent free-flow when the administration set is removed from the instrument
- Downstream Occlusion Detection system minimizes nuisance patient-side occlusion alarms caused by momentary kinking of tubing, IV pushes, etc.
- Guardrails® Drug Library
- Guardrails® Limit
- KVO Rate Adjust
- Multidose Mode
- Occlusion Pressure
- Dimensions: 3.3"W x 8.9"H x 5.5"D
- Weight: 2.5 lbs



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.