



auxoMEDICAL

Atria 6100 - Fully Refurbished

The Atria 6100 is a practical, portable ECG / EKG machine for hospitals and cardiology clinics. The Atria 6100 portable ECG helps streamline workflow with a friendly user interface and clear, readable color display for easy reading of ECG waveform data.

Features:

- Large color monitor to view 3-, 6-, and 12-lead ECG printouts in real-time
- Portable 6100 ECG Stores up to 300 electrocardiograms
- Manual or automatic 3-, 4-, 6-, and 12-lead ECG printouts
- Standard with wired Ethernet, wireless 802.11, and Bluetooth; complies with FAA requirements
- Bi-directional communication using HeartCentrix connectivity software
- AccuPrint feature restricts ECG printing if lead wires are not properly attached
- The standard keyboard, dedicated function keys, large color monitor, and intuitive menu make this portable ECG easy to use.
- The ECG comes with Ethernet connectivity as a standard feature
- ECG It's compatible with wireless 802.11 and Bluetooth, complies with FAA requirements, and has bi-directional communication with our HeartCentrix connectivity solution.
- Pediatric applications
- The AccuPrint feature restricts ECG printing if lead wires are not properly attached, ensuring a clear printout and eliminating retakes.
- This durable ECG is covered by an industry-leading warranty and service program for worry-free usage.
- Dimensions: 5.5" (H) x 13.125" (D) x 15.0" (W) (14.0 cm x 33.4 cm x 38.1 cm)
- Weight: 11 lb (5 kg) (including external power supply)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.