



## B Braun Perfusor fm Pump IV - Fully Refurbished

B. Braun Space is a modular infusion system that is smaller, lighter, and easier to handle than any other existing configuration. Leading the way into the future, this system meets all demands of an outstanding safety concept. B. Braun Space can also be integrated into the data communications network of every advanced hospital operation. Performance data Automatic rate calculation by calculating volume over time

### Features:

- Dose rate mode:
  - Automatic calculation of delivery rate based on the entry of drug concentration in conjunction with the desired dose rate (e.g.: mg/kg/min)
- Drug Library (Must obtain From Braun)(not included):
  - Up to 1500 drug names including therapy data and information can be stored in 30 categories
  - Up to 50 care units and 15 patient profiles possible
- Drug specific soft and hard limits as well as default values can be specified
- Operator interface
- Standardized and intuitive user interface for all B. Braun Space infusion pumps due to menu operation and arrow keys the user has to look at the display when doing adjustments
- Safety concept
- Integrated piston brake prevents free-flow during syringe change
- Automatic drive technology ensures startup and convenience when changing syringes
- Automatic bolus reduction triggered by occlusion alarm
- Data lock on 3 safety levels (parameters and disposables lockable)
- LED alarm indicator with clear alarm message in display
- Dimensions: 36 x 20 x 11 cm.
- Rate: 1 mL / hr to 1800 ml / h
- Volume: 0.1 ml to 50.0 ml,



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.