



B. Braun Infusomat Space - Fully Refurbished

Designed for acute care adult and pediatric facilities, the compact 9" x 6" Infusomat Space large volume pump weighs just 3 pounds. The lightweight large volume pump allows for flexibility and easy transport within a number of care settings. Along with industry-leading safety features, the lightweight mobile design allows you to quickly and easily stack up to 24 pumps in a single bedside tower. Equipped with an extensive drug library, easy navigation and intuitive interface, this small infusion pump packs state-of-the-art versatility that promotes efficiencies across numerous care settings.

Features:

- Intuitive user interface
- Alarms are audible, visible and differentiated
- Strong lightweight pole clamp
- 1,500-entry drug library available
- Independent Modularity
- Easily Stackable
- KeyGuard™
- Wireless Integration: Allows for complete interoperability
- Single channel design



Specs:

- Dimensions: 249 x 68 x 152 mm (WxHxD)
- Weight: Approx. 1.4 kg
- Display: Backlit graphic display, ~40° read angle from all sides
- Flow Rates: 0.01-1800 ml/h
- Voltage: 11 - 16 V DC supplied by external Space Power Supply or by SpaceStation
- Battery operating time: Minimal power consumption of the devices and new battery technology ensure long operating times, e.g. 16.8 hrs @ 1 ml/h, 16 hrs @ 5 ml/h, 14 hrs @ 20 ml/h

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.