



BIODEX Echo Pro 058-700 Ultrasound Table - Fully Refurbished

The Echo Pro™ Echocardiography Table is engineered with the sonographer in mind. The ergonomic design provides a barrier-free imaging experience, resulting in less sonographer injury and better quality images. Features make it easier for the sonographer to access the patient from all sides. The cardiac scanning cushion now drops down and releases from either side of the table with nothing protruding. The sonographer has open access to the left thorax area for an unobstructed apical approach. New side rails are flush to the table and fold completely beneath the table out of the sonographer's way, eliminating potential obstruction.

Features:

- Ergonomic design helps prevent musculoskeletal injuries
- Extra wide top with 500-lb patient capacity
- Cardiac drop-down cushion
- Industry's largest cardiac scanning cutout – 9.5" w x 15.75" l
- Fowler positioning, 0-80 degrees, redesigned to accommodate natural body extension and sitting position without slide or shift
- Height range down to 23" accommodates wheelchair transfers and those with ambulatory difficulties
- Motorized height, Fowler positioning, Trendelenburg/Reverse Trendelenburg and auto-level adjustment
- Naugahyde® antimicrobial mattress cover with Advanced BeautyGard® provides protection against bacteria
- Central floor-locking system maintains table stability
- Durable hand controller
- Allows seated or standing work positions
- Overall Dimensions: 73"(L) x 30"(W); 35"(W) with optional side rails



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.