

Baxter 6301 - Fully Refurbished

Treat the varied conditions of patients with this advanced infusion system that offers one pump for all of your diverse infusion needs. The Baxter 6301 allows you to administer fluids precisely and continuously, freeing up time to work with other patients. Its anti-tampering system means you can have the confidence that your patients are receiving the exact treatment they need and nothing else. The Baxter 6301 enables you to implement patient-specific treatments and protocols with flawless precision. This is a dual-channel volumetric infusion pump. Once volume and time are selected, flow rate calculation is automatic in this model.

Features:

- Remote operation and monitoring
- Recharge battery to 80% within 8 hours
- Delivers a variety of fluids over a broad range of infusion rates
- Enhanced safety through flow check occlusion alarm
- Provides flow rate calculation for errorless administration of fluids
- Front panel lock outs to protect against tampering and leave you in control
- Automatic Restart once occlusion clears
- Automatic Piggybacking of secondary medications
- Two separate pump channels
- Configurable settings
- Five-hour memory retains infusion data after power-off.



Specs:

- Weight: 8.1 kg (17.9 lbs)
- Dimensions: 33 cm W × 13 cm D × 29 cm H (13" W × 5.1" D × 11.4" H)
- Range of Flow Rate:
 - 1-99.9 mL/hr in 0.1 mL/hr increments
 - 1-9,999 mL/hr in 1 mL/hr increments
- 12V battery: lasts 6 hours at 1,400 mL/hr

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.