

Baxter AS50 Syringe Infusion Pump - Reconditioned

The Auto Syringe AS50 Infusion Pump is designed to meet the fluid and drug delivery requirements of today's changing clinical environment. It provides accurate, continuous or intermittent infusion of intravenous solutions, drug solutions, whole blood, and packed red blood cells. The pump is indicated for infusion via intravenous (IV), intra-arterial (IA), epidural, or subcutaneous routes of administration.

Features:

- Digital syringe pump
- Standard disposable syringes from 1mL to 60mL
- Features flexible drug libraries
- Provides accurate infusion of intravenous solutions, drug solutions, whole blood, and packed blood cells
- Indicated for the following routes of administration: intravenous, intra-arterial, epidural or subcutaneous
- Dimensions: Approx. 3.4"W x 2.6"D x 10"H
- Weight: Approx. 3.2 lbs.
- Flow Rate Accuracy: $\pm 3\%$ of full scale plunger travel (not including syringe tolerance)
- Volume Accuracy: $\pm 3\%$, or 0.007" of travel, whichever is greater (not including syringe tolerance)
- Syringes:
 - B-D: 1, 3, 5, 10, 20, 30, 60 mL plastic
 - Monoject: 1, 3, 6, 12, 20, 35, 60 mL
 - Terumo: 1, 3, 5, 10, 20, 30, 60 mL
- Flow Rate Range: 0.01 mL/hr to 438 mL/hr, depending on the syringe

Additional Options:

- IV Pole, 5 caster 2 or 4 hook



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.