



Baxter Infus O.R.- Fully Refurbished

The Infus O.R.™ syringe infusion device aids in administering many intravenous agents given during anesthetic procedures. It conveniently delivers narcotics, muscle relaxants and vasoactive drugs via SMART LABELS – a system that displays appropriate infusion parameters and rates for specific medication, helping decrease the likelihood of calculation errors in fast-paced surgical environments.

Features:

- Alarms: Infusion complete, Battery, Malfunction.
- Smart Labels are magnetically coded /Drug Specific
- Adjust settings easily at any time during infusion
- Modify dose without multiple reprogramming steps
- Drug specific Smart Labels that help minimize the potential for medication errors in dose programming
- Weight 2 LB. with batteries.
- Operating Voltage 5-7 Volts
- Uses 4 "C" Alkaline Batteries.
- Size -9.2 x 4.5 x 2.0 in.
- Flow Rates -0-600 mL/HR.
- Designed to aid in many intravenous agents given during anesthetic procedures.
- Provides for the convenient delivery of narcotics, muscle relaxants, vasoactives and other drugs routinely given during anesthesia.
- Easy to use and set up.
- Drug delivery choices are easily selected and changed using front panel rotary switches.
- Change dial settings with a single hand for an immediate response.
- Save calculation time.
- Minimize risk of calculation error.
- Change Infusion Parameters/Rates with SMART LABELS
- SMART LABELS are magnetically coded/drug specific.
- Select from the available SMART Label library.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.