



## Baxter PCA II - Fully Refurbished

The Baxter PCA II pump is a syringe infusion device capable of both continuous and intermittent infusion of parenteral fluids in a hospital setting. The Baxter PCA II Syringe Pump offers the convenience of safe and precise Patient Controlled Analgesia in a flexible system designed with advanced technology. This is progressive, innovative, and effective pain management in a sophisticated, yet simple-to-operate, instrument.

### Features:

- Internal safety checks help minimize the risk of tampering and errors
- Standard programming routine helps reduce procedural errors and minimizes training
- Built-in history logs successful injections, attempts, alarms, and other operating details
- Dosage history available in hard copy when connected to printer
- Large, back-lit display easily readable even at night
- Easily programmed infusion modes:
  - PCA only
  - Basal plus PCA
  - Continuous Basal.
- Program by milliliters or milligrams.
- Select the syringe option to best suit your needs:
  - Pre-filled syringes for convenience
  - Standard syringes for cost savings.

### Specs:

- Flow Accuracy: 3 % average linear displacement
- Flow Rate Range: 0.1 to 150 ml/h
- Infusion Pump Type: Patient-Controlled Analgesic (PCA)
- Height: 13 in
- Length: 2.8 in
- Weight: 4.2 lbs with batteries
- Width: 6.3 in



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.