



# auxoMEDICAL

## Biolight M8000 - Fully Refurbished

Improve patient outcomes with a comprehensive patient monitoring system covering the needs of multiple departments within hospitals and outpatient facilities. Standard configuration monitors ECG, HR, RESP, NIBP, SpO2, PR, TEMP, and Battery. Optional Configurations include 2-TEMP, 2-IBP, Recorder, EtCO2(side stream, main stream), Nellcor SpO2, and IC.

### Features:

- 10.4" TFT LCD Display, High Resolution
- Portable design, stylish and convenient
- Store and review 1000 pieces of alarm events
- Auto identify ECG Waveforms of 3-lead or 5-lead, ideal for OR & CCU
- NIBP dual Over Pressure Protection
- Isolated float, anti defibrillation and anti high -frequency interference
- Networking capacity
- ST arrhythmia analysis
- Tabular and graphic trends data over 72 hours
- Low-end solution, short trend
- 7 lead ECG Waveform display in the same screen
- Large number display for long distance vision
- Digital SPO2 technology offers anti-motion and anti-low perfusion functions
- Up to 66 beds network capacity to connect with HIS 120 hours tabular and graphic trends
- Dimension: 318mm (W) x 264mm (H) x 152mm(D)
- Weight: 4.5kg
- Power Supply : 100~240V; AC: 50~ 60 Hz
- Alarm Indication: Visual and audible alarm
- 1000 groups of alarm events review
- Application: Adult, Pediatric and neonate
- Display Resolution 800x600 or higher
- Sealed Lead Acid cell battery
- Charging time of <10hrs



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.