



## Biolight M9000- Fully Refurbished

Biolight BLT M9000 Capnograph Patient Monitor with CO2 side stream. Suitable for adult, pediatric and neonate. It includes 7 parameters and has a large 12.4" screen LCD. It is recommended for use in hospitals and clinics for: ECG, Respiration-breathing also NIBP-non invasive blood pressure, Pulse rate, SPO2, Temperature, (optional IBP and Multi-gas/O2) parameters.

### Features:

- 12.1" LCD display- high resolution
- Digital SPO2 technology offers anti-motion and anti-low perfusion functions
- 168 hours storage
- Bed to bed view
- Product Specification/Models
- Safety: IEC60601 (GB9706)
- Dimension: 318mm (W) x 264mm (H) x 152mm(D)
- Power Supply : 100~240V; AC: 50~ 60 Hz
- Alarm Indication: Visual and audible alarm
- 1800 groups of alarm events review
- Application: Adult, Pediatric and neonate
- Max 7 waveform display simultaneously
- Big font convenient for long distance observation
- Supporting for ST, OTC and arrhythmia analysis
- 10 minutes ECG waveform reviewing
- 120h trend data/trend maps storage and reviewing
- 500 groups NIBP measurement
- 1000 groups Alarm events reviewing
- Supporting for nurse call, pitch tone, built-in recorder wireless or wired networking



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.