



Biolight M9500- Fully Refurbished

M9500 Human Multi-Parameter Patient Monitor adequately meets clinical needs with 7/12 leads in the same screen. Digital SpO2 technology with anti-motion, anti-low perfusion, and anti-electrosurgical interference. Can be connected to a central monitoring system that is wired or wireless.

Features:

- 15" Clear and Colorful TFT LCD Screen
- Store and Review options
- 30 minutes of ECG waveform review
- Automatic recognition of applicable ECG cables 3/5 leads
- 72 hour data/trend maps storage and review
- Built in detachable rechargeable lithium battery
- Dual alarm light, Physiological alarm light, Equipment technical alarm light
- Back light buttons
- Option for thermal recording with real-time and alarm trigger printing function
- 1000 groups NIBP measurement
- 1800 groups alarm events recording



Specs:

- Size 366mmX330mmx 162mm
- Weight 5 kg
- Battery
 - Type : Rechargeable Lithium Battery
 - Charge time: ~ 6 hours
 - Operating time under the normal use and full charge: 120 minutes

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.