



auxoMEDICAL

BladderScan BVI 3000 - Fully Refurbished

The BladderScan BVI 3000 is a portable, 3D ultrasound instrument that quickly, accurately, and non-invasively measures urinary bladder volume and post-void residual (PVR). The BladderScan BVI 3000 calculates bladder volume using patented Vmode ultrasound technology. A BVI 3000 exam is quick, easy, and comfortable for the patient. When you press the scan button, within seconds, the instrument measures ultrasonic reactions on multiple planes inside the body producing a 3D image. Bladder volume measurements made with Vmode ultrasound are more accurate than those from conventional 2-dimensional ultrasound, as they are based on this more complex, 3D image.

Features:

- LCD screen displays bladder position and volume; aiming icon assists in proper probe placement
- 3D display shows the bladder in 2 cross-sectional images to verify an accurate scan
- Easy-to-print exam results for patient records or reimbursement procedures
- Helps document cost savings and catheterizations-prevented based on user-determined statistical input
- Compact console with easy-to-read LCD screen
- Two NiMH rechargeable batteries
- Battery charger
- Onboard printer for patient records or reimbursement
- Ultrasound gel
- User's Manual, Quick Reference Guide and Quick Reference Cards



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.