



Braun Vista Basic Infusion Pump - Fully Refurbished

The Vista Basic was designed to ensure simplicity in all aspects of its use: from set installation, to programming, to special mode use and alarm correction. Based on a design that has sold over 100,000 units worldwide, the Vista Basic provides high reliability packaged in one of the easiest pumps to use. As a workhorse in the outpatient market, it offers a rugged and reliable design, all relevant therapies such as continuous, piggyback and program mode, and a selection of cost-effective, straight line gravity sets.

Features:

- Robust design of mechanical components ensures maximum durability
- Battery life of approximately 3 years before replacement is required
- Helpful prompts walk the user through the programming process
- Straight-line tubing loads quickly into the pump. Simply enter the rate and volume and press Start. Helpful prompts walk the user through the programming process.
- Easy to use, teach and learn. Perfect for home care settings.
- Pump based free-flow protection helps prevent accidental boluses when door is opened.
- Optional set-based free-flow protection helps prevent free-flow when set is removed from the pump.
- DEHP free sets protect patients.
- Needle-free sets protect both the patients and the clinician.
- Adjustable occlusion settings help protect patients from vein injuries. Handy bag hook saves space and provides ease of mobility with no need for an IV pole.
- Oversized display and ergonomic handle makes transport easy under difficult conditions.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.