



## Buffalo Filter Porta PlumeSafe 604- Fully Refurbished

Porta PlumeSafe® 604 is an ideal surgical smoke evacuator for veterinarian and non-acute care settings where moderate flow is needed. The unit is portable, easy to use, and features patented technology to track filter life. EZLink®, an automatic activation device, can be used to remotely turn the evacuator on and off, saving filter life. ViroSafe® Fluid Trap is designed to capture fluids to help maintain the integrity of the filtration system and assist with maintaining filter life.

### Features:

- Lightweight & Portable – Allows for easy transport or tabletop use.
- LED Display – Indicates filter life and suction power.
- Adjustable Flow Control – Suction setting can easily be adjusted in 10% increments.
- Filter Tracking Technology – Automatically tracks and records filter life through a patented technology.
- Quiet Operation – Addresses noise concerns.
- Economical – Provides 0.1 to 0.2 micron filtration @ 99.999% efficiency
- Simple filter changes
- Ceiling-mounted or Wall-mounted
- Germicidal
- High-efficiency
- Measures room pressure
- Portable
- Zone-specific minimum air changes per hour

### Specs:

- Electrical: 100/120V, 10 Amps, 50/60 Hz
- Maximum Flow Setting (cfm-U.S.): Standard Hose 1 3/8 in. = 79
- Noise Level, dBA: Max 49, Min. 69.5
- Weight: 19.5 lbs
- Dimensions: 13.5 H x 8 W x 17 D inches



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.