



Cardiac Science Powerheart AED G3- Fully Refurbished

In the chaos that follows sudden cardiac arrest (SCA), concerned but untrained bystanders are scared. Do they know what to do? The Powerheart™ AED G3 Plus with RescueCoach™ voice prompts to talk rescuers through the steps. When the rescuer applies the defibrillation pads, the G3 Plus Automatic analyzes the heart rhythm and “knows” when to deliver (or not deliver) the shock. The shock is delivered automatically, with no button to push, and no human intervention (semi-automatic version also available). After the shock, the unit prompts for CPR with a built-in metronome that sets the pace for proper chest compressions.

Features:

- RescueCoach
- Daily, weekly, and monthly Rescue Ready self testing
- Pads tested for presence and functionality
- Pediatric pad availability
- 7-year device warranty
- 4-year, full battery replacement guarantee
- Download patient event data
- Non-committed shock
- Pacemaker pulse detection
- Synchronized shock
- Fully automatic mode
- Semi-automatic mode
- Color ECG display, manual override, and 3-lead ECG monitoring capabilities.
- Rescue Ready technology to self-test all critical components
- A built-in metronome (which can be switched off) to set the pace for CPR compressions.
- Dimensions(H x D x W): 8 cm x 31 cm x 27 cm (3.3 in x 12.4 in x 10.6 in)
- Weight: 3.1 kg (6.6 lb)
- RHYTHMx™ analysis algorithm lets first responders proactively monitor for life-threatening arrhythmias



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.