



## CardioCap 5 with 5 Agent Gas with Printer- Fully Refurbished

GE™ CardioCap 5 monitor is a compact, all-in-one monitor designed for use in the operating room, ambulatory surgery unit, induction room and PACU. The Cardiocap/5 offers full vital signs monitoring capabilities from the monitoring of oxygenation and circulation to advanced airway gas analysis and ventilation measurement with Patient Spirometry." Compatible with Datex-Ohmeda S/5 Network and Central for monitor-to-monitor communication and centralized viewing. This CardioCap 5 monitor includes a printer.

### Features:

- Small footprint bedside monitor with large 10.4" active matrix color display with adjustable brightness
- Proven easy-to-use S/5 user interface
- Shows up to 6 waveforms
- Color-coded, prioritized alarms
- 4 user-configurable modes for care-specific preconfigurations
- Waveform snapshot to store and recall waveform samples for later analysis
- 5 or 30-minute minitrends continuously on screen for quick historical view
- 24-hour graphical trends of all parameters
- Optional built-in, three-channel thermal array recorder
- Wall mount, roll stand available for best workplace ergonomics
- Hemodynamics
- Impedance respiration
- Pulse oximetry (SpO<sub>2</sub>)
- NIBP
- Airway Gases: CO<sub>2</sub> (N<sub>2</sub>O), Patient O<sub>2</sub>, N<sub>2</sub>



### Specs:

- Dimensions: 13.0" (W) × 8.7" (D) × 11.8" (H)/330 × 220 x300 mm
- Weight: <11.2 kg, 24.8 lbs (F-MXG); <10.2 kg, 22.6 lbs (F-MX)
- Power: 100-240 Vac ±10%, 60/50 Hz
- Back-Up Battery: 15 min guaranteed, charging time typically 5 h

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.